# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

## User's Manual

# VMA ELISA

*Enzyme Immunoassay for the quantitative determination of  
Vanillylmandelic Acid (VMA) in human urine*

**DE1620**

**96 wells**

**Table of Contents**

1	NAME AND INTENDED USE .....	3
2	PRINCIPLE OF THE ASSAY.....	3
3	MATERIALS PROVIDED IN THE KIT .....	4
4	MATERIALS REQUIRED BUT NOT PROVIDED .....	4
5	PREPARATION FOR THE ASSAY .....	4
6	WARNING AND PRECAUTION .....	4
7	STORAGE AND STABILITY.....	4
8	SPECIMEN COLLECTION AND HANDLING .....	5
9	ASSAY PROCEDURE.....	5
10	PROCEDURAL NOTE.....	8
11	LIMITATIONS .....	8
12	QUALITY CONTROL.....	8
13	EXPECTED VALUES .....	9
14	PERFORMANCE CHARACTERISTICS .....	10

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## 1 NAME AND INTENDED USE

The VMA ELISA is designed for in vitro quantitative measurement of vanillylmandelic acid (VMA) concentration in patients' urine.

**For in vitro diagnostic use only.**

### 1.1 Background

Among spectrophotometric procedures the most widely used procedure was developed by Pisano et al (1) measuring vanillin directly at 360 nm. The colorimetric methods in general, however, suffer from interfering compounds and intensive labor requiring a number of manual steps. The use of paper chromatographic methods (2,3) employing diazotized p-nitroaniline and a thin layer chromatographic (TLC) procedure on silica gel (4) were reported for the analysis of extracted VMA but are considered to be of only historical interest. More recently, high performance liquid chromatography (HPLC) procedures are reported using a variety of detectors including ultraviolet light (5-7), amperometric detection (8-11), and postcolumn reaction (12, 13). These procedures are considered to be a considerable improvement for overall specificity and sensitivity of analysis.

### 1.2 SUMMARY AND EXPLANATION

Catecholamines include dopamine (found mostly in the central nervous system), norepinephrine (mainly in the sympathetic nervous system) and epinephrine (mainly in the adrenal medulla). They are stored as inactive complex. Released catecholamines, having a short half-life, are taken up by sympathetic nerve endings, or metabolized by the liver and kidney and excreted.

Vanillylmandelic acid (VMA) and 4-hydroxy-3-methoxy-mandelic acid (HMMA) are the end product of both epinephrine and norepinephrine catabolism. Quantitation of the acidic metabolites has long proven to be a reliable diagnostic as well as commonly used follow-up procedure for pheochromocytoma and other catecholamine-, producing tumors (14-16). The prevalence of pheochromocytoma is 0.1% to 0.2% of hypertensive patients (17, 18). Pheochromocytomas have a long record of misdiagnosis due to the metabolic, cardiac and gastrointestinal symptoms that can mimic many other diseases. Undetected, or mistreated it can be fatal. However, since this is a surgically curable disease, early diagnosis by demonstration of excess VMA excretion is critically important.

## 2 PRINCIPLE OF THE ASSAY

The VMA ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the competition between VMA coated on a microtiter well and that in urine for the monoclonal antibody. Outlined steps are:

1. *Sampling and reaction:* The samples are incubated in the wells with horseradish peroxidase conjugated anti-VMA monoclonal antibody.
2. *Washing:* Unbound VMA and the antibody bound to urinary VMA are removed by washing with 0.9% NaCl solution.
3. *Enzyme Reaction (Color Development):* The amount of bound peroxidase is inversely proportional to the concentration of the VMA present in the urine sample. Upon addition of the substrate (TMB), a blue color is developed, then it is changed to yellow by adding Stopping Solution. The intensity of this is inversely proportional to the concentration of VMA in the Calibrator or urine sample.
4. *Absorbance Detection:* After addition of Stopping Solution, absorbance is measured at 450 nm. And the readings are converted into the concentrations from the Calibration curve.

### 3 MATERIALS PROVIDED IN THE KIT

1. VMA Coated **Microwell Plate**:  
1 VMA coated 96-microwell plate.
2. **Anti-VMA-Enzyme Conjugate**:  
Horseradish peroxidase conjugated to an anti-VMA monoclonal antibody, 10 ml.
3. **VMA Color Developing Reagent**:  
Tetramethylbenzidine (TMB) solution, 20 ml.
4. **VMA Stopping Solution**:  
Mixture of acidic solution and Hydrochloric Acids (20 ml)
5. **VMA Calibrators**:  
of 0, 0.0625, 0.125, 0.25, 1, 2, 4 and 8 µg/ml  
in phosphate buffered saline, 0.01M, pH 7.4 (1.0 ml each).

### 4 MATERIALS REQUIRED BUT NOT PROVIDED

1. Plate reader with 450 nm filter
2. pH meter or pH paper with the range of 5.0-10.0
3. Pipettor with tips for 10, 50, and 100 µl
4. Pipettor with tips for 50 and 100 µl
5. Volumetric cylinders, 10 and 100 ml
6. Volumetric and serological pipettes, 10 ml
7. Disposable test tubes or vials
8. 5N NaOH solution
9. 5N HCl solution
10. Plate washer (optional)
11. Plate shaker (optional)
12. Sodium Chloride or Saline Solution
13. 0.01M phosphate buffered saline, pH 7.4

### 5 PREPARATION FOR THE ASSAY

1. Prepare 0.01 M phosphate buffered saline, pH 7.4.  
This solution is used to dilute all unknown urine samples prior to analysis.
2. Before beginning the test, bring all urine samples and reagents to room temperature (15-30 °C) and mix well.
3. Set up all reagents and urine samples before running the assay. The entire test procedure must be performed without any interruption in order to get the most reliable and consistent results.

### 6 WARNING AND PRECAUTION

- The VMA KIT is designed for in vitro use only.
- The components in this kit are intended for use as an integral unit.
- The components of different lots should not be mixed.
- Do not use the VMA Calibrators in this kit for other purposes (e.g. HPLC).
- Use a new pipette tip for each Calibrator or urine sample to avoid cross-contamination.

### 7 STORAGE AND STABILITY

- Store the kits at 2-8°C in refrigerator.
- Keep microwell plates in dry bag with desiccants. Open the bag only when needed.
- Expiration dates of the reagents are stated on their labels.
- Color Developing Reagent should be colorless.
- Protect the reagents and reaction mixture from exposure to direct sunlight.

## 8 SPECIMEN COLLECTION AND HANDLING

A 24-hour urine specimen should be collected with 10 ml of 6 M HCl as a preservative.

Overnight or randomly collected urine should be acidified to a pH between 2 and 3 immediately after collection. Record the total volume and save 1-5 ml for the analysis of VMA and total creatinine.

All samples should be refrigerated until tested. Centrifuge turbid urine samples containing crystals or sediment.

## 9 ASSAY PROCEDURE

### 9.1 Preparation of Reagents.

#### Washing Solution:

Dissolve 9 grams of NaCl in 1 liter of deionized or distilled water. Commercially available normal saline can also be used.

### 9.2 Preparation of Samples

1. Take 1.0 ml of acidified urine and transfer to a disposable tube in which the pH of urine sample can be readjusted.
2. Bring pH of all samples within the range of 6 and 9 by stepwise addition of small amounts of 5N NaOH (e.g. 5 µl) while checking pH either with a pH meter or using pH paper.
3. Dilute pH re-adjusted samples at a 1:10 ratio with phosphate buffered saline. The pH for diluted samples should be between 7.0 and 8.0.

### 9.3 Standard Procedure for the Assay

1. Make work sheet with Calibrators and sample identification.
2. *Sampling:*
  - a. Dispense 50 µl of VMA Calibrators into appropriately designated wells.
  - b. Dispense 50 µl diluted (1:10) controls or samples to respective wells.
3. *Addition of the Anti-VMA-Enzyme Conjugate:*

Dispense 50 µl of Anti-VMA-Enzyme Conjugate to each well, using a pipettor.
4. *Antigen-Antibody Reaction:*

Mix the plate by moving it back and forth slow horizontal movements for a minute. A plate shaker can be used for this purpose also. Allow the plate to stand at 15-30°C, room temperature for 1 hour.
5. *Washing:*


**Wash only once.** Removing incubation mixture by decanting the plate into a sink and blotting the plate on absorbent paper.  
Dispense 300 µl of normal saline into each well. Remove saline by decanting the plate and blotting it on absorbent paper. This also can be done by a plate washer.
6. *Enzyme Reaction:*

Dispense 100 µl of Color Developing Reagent to the well and allow it to stand at 15-30°C, room temperature for 25 minutes.
7. *Stopping Color Development:*

Dispense 100 µl of Stop Solution to the wells.
8. *Absorbance Measurement:*

Any microwell reader capable of detecting absorbance at 450nm may be used.

## PROCEDURE SUMMARY FLOW DIAGRAM

**9.4 Calculation of Results**

- Using semi-log linear paper (Fig. 3) or log-logit paper (fig. 4) the Calibration curve is generated by plotting VMA concentrations on the abscissa and the absorbance on the ordinate. VMA concentration for each unknown sample is obtained from the Calibration curve.

For example:

Description	Absorbance (450 nm)		Average of B/Bo(%)	VMA (µg/ml)
VMA Calibrator (µg/ml)				
0	1.487	1.466	100.0	
0.0625	1.255	1.258	85.1	
0.125	1.037	1.052	70.7	
0.25	0.818	0.839	56.1	
1	0.359	0.378	25.0	
2	0.245	0.261	17.1	
4	0.145	0.145	9.8	
8	0.078	0.084	5.5	
Sample A	1.299	1.218	82.9	0.072 x 10
Sample B	0.769	0.746	51.3	0.285 x 10
Sample C	0.522	0.522	35.4	0.600 x 10

The results obtained above indicate VMA concentration in µg/ml. When the total VMA in 24-hour urine sample is required;

$$\text{VMA } (\mu\text{g/ml}) \times \frac{\text{urine volume (ml)}}{1000} = \text{VMA mg / 24 hours}$$

or, when VMA/Creatinine value is required;

$$\text{VMA } (\mu\text{g/ml}) \div \frac{\text{Creatinine (mg/dl)}}{100} = \text{VMA } \mu\text{g/ml Creatinine or}$$

$$= \text{VMA mg/g Creatinine}$$

- Calculation can be made with a computer set so as to draw calibration curves based on 4 coefficient log-logit.


Fig.3 Calibration curve of VMA using semi-log linear paper


Fig.4 Calibration curve of VMA using log-logit linear paper

**NOTE:** This is a representative curve and should not be used to calculate results for unknown samples.

## 10 PROCEDURAL NOTE

1. It is very important to wash the microwells thoroughly, yet uniformly and remove any residual liquid from the wells to achieve optimal results.
2. Pipette Calibrators or urine samples into the bottom of each well. Vortex-mixing or shaking of wells after each pipetting is not required.
3. Absorbance is a function of the incubation time and temperature. It is, therefore, recommended to ensure the equally elapsed time for each pipetting without interruption.

## 11 LIMITATIONS

1. This VMA ELISA kit is designed for the quantitative determination of VMA in urine only.
2. All samples with VMA concentrations greater than 8 µl/ml should be repeated on much larger dilution(s), e.g. 1:20 or more.
3. Interference by Sodium Azide: As Sodium Azide inhibits the enzyme reaction, urine or any buffer used to dilute urine samples containing sodium azide as an antiseptic can not be used.

## 12 QUALITY CONTROL

In order to monitor precision of the analytical performance, it is recommended that commercially available urine control samples be included in every run.


### 13 EXPECTED VALUES


Each laboratory should determine a normal range to conform to the characteristics of the population being tested. The range given here was determined from 24-hour urine collections on 280 subjects.

Urinary creatinine was measured on Astra to assess the completeness of each collection, and mg VMA per gram creatinine was calculated.

	VMA mg/day	VMA mg/g Creatinine
<b>Number of samples (n)</b>	280	280
<b>Mean Value (x)</b>	7.92	6.87
<b>± 2 S.S. range</b>	0.40 – 15.44	0.00 – 14.37
<b>Reference Range</b>	up to 15.0	up to 14.0
	(Fig. 5)	(Fig. 6)

**Fig. 5** Distribution of total VMA excretion per day in 280 cases of normal controls.

**Fig. 6** Distribution of total VMA concentration per day in 280 cases of normal controls.

## 14 PERFORMANCE CHARACTERISTICS

### 14.1 Recovery

Recovery studies were performed in urine samples from four healthy adults, to which VMA was added at various concentrations. Recovery was determined by dividing observed value by expected value for each prepared sample.

Sample	Baseline (µg/ml)	Added (µg/ml)	Expected Value (µg/ml)	Observed Value (µg/ml)	Recovery (%)
A	1.76	0.00			
	0.88	9.60	10.48	10.10	96
	1.32	4.80	6.12	5.92	97
	1.54	2.40	3.94	3.90	99
	1.65	1.20	2.85	2.78	98
B	6.16	0.00			
	3.08	9.60	12.68	12.85	101
	4.62	4.80	9.42	9.55	101
	5.39	2.40	7.79	8.00	102
	5.77	1.20	6.97	7.29	104
C	11.71	0.00			
	5.85	9.60	15.45	15.85	102
	8.78	4.80	13.58	18.77	101
	10.24	2.40	12.64	12.61	99
	10.98	1.20	12.18	13.30	109
D	8.43	0.00			
	4.21	9.60	13.81	15.14	109
	6.32	4.80	11.12	11.34	102
	7.37	2.40	9.77	10.34	105
	7.90	1.20	9.10	9.27	101

### 14.2 Linearity and Parallelism

Four urine samples and two commercial urine controls Lyphochek 1 and 2 (Bio-Rad) were serially diluted with phosphate buffered saline.

The ratio (B/Bo) of absorbance for each dilution (B) to the absorbance of 0 µg/ml Calibrator (Bo) was calculated and plotted on log-logit paper.

### 14.3 Precision

#### 14.3.1 Intra-assay

Intra-assay coefficient of variation was evaluated in three urine samples at different VMA concentrations.

Intra-Assay	Sample A	Sample B	Sample C
<b>N</b>	26	24	24
<b>Mean (µg/ml)</b>	0.96	8.25	11.24
<b>S.D. (µg/ml)</b>	0.06	0.52	0.71
<b>C.V. (%)</b>	6.0	6.3	6.0

#### 14.3.2 Inter-assay

Inter-assay coefficient of variation was evaluated at three different concentrations, by analyzing the samples in 16 to 19 separate occasions.

Inter-Assay	VMA 1	VMA 2	VMA 3
<b>N</b>	16	19	19
<b>Mean (µg/ml)</b>	1.193	2.821	9.698
<b>S.D (µg/ml)</b>	0.086	0.232	0.686
<b>C.V. (%)</b>	5.7	8.2	7.1

#### 14.4 Specificity

The following substances were tested for cross-reactivity of the assay. Cross-reactivity is expressed in terms of percentage of the concentration of each substance that produced 50% displacement.

Substance	Cross-reactivity (%)
	VMA
Vanillylmandelic Acid	100
Homovanillic Acid	<0.01
DL-3,4-Dihydroxymandelic Acid	4
3,4-Dihydroxyphenylacetic Acid	<0.01
Metanephrine	<0.01
Vanillylpyruvic Acid	4
Vanillic acid	<0.01
Dopamine	<0.01
5-Hydroxy-3-indolacetic Acid	<0.01
Vanillyllactic Acid	<0.01
3-methoxy-4-hydroxyphenyl Glycol	<0.01

#### 14.5 Sensitivity

The sensitivity of this test is higher than 0.0625 µg/ml. The minimal detectable concentration of VMA is estimated to be 0.035 µg/ml. The minimal detectable concentration is defined as the concentration of VMA which corresponds to the absorbance that is two standard deviations from the mean absorbance of 20 determinations of zero dose VMA.

#### 14.6 Sample Stability

Sample stability was studied in two different urine samples at 4 °C and -20 °C for up to 5 days and 50 days, respectively. The results confirm that VMA is stable at 4 °C storage at least for 5 days tested, and for up to 50 days tested when stored at -20 °C.

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**15 REFERENCES /LITERATUR / REFERENCIAS**

1. Pisano, J.J., Crout, J.R. and Abraham, D.: Determination of 3-methoxy-4-hydroxymandelic acid in urine. *Clin.Chim. Acta* 7: 285-291, 1962.
2. Vahidi, H.R., Roberts, J.S., San Filippo, J., Jr. and Sankar, D.V.: Paper chromatographic quantitation of 4-hydroxy-3-methoxymandelic acid (VMA) in urine. *Clin. Chem.* 17: 903, 1971.
3. Gitlow, S.E., Mendlowitz, M. and Bertain, L.: The biochemical techniques for detecting and establishing the presence of pheochromocytoma. *Am. J. Cardiol.* 26: 270, 1970.
4. Badelly, M., Routh, M.W., Gump, B.H. and Gigliotti, H.J.: Thin-layer chromatographic method for urinary 4-hydroxy-3-methoxymandelic acid (vanilmandelic acid). In Gooper, G.R., editor: *Selected methods in clinical chemistry*, Washington, D.C., 1977. American Association for Clinical Chemistry, pp. 139-145.
5. Yoshida, A., Yoshioka, M., Tanimura, T. and Tamura, Z.: Determination of vanilmandelic acid and homovanillic acid in urine by high speed liquid chromatography. *J. Chromatogr.* 116: 240, 1976.
6. Felice, L.J. and Kissinger, P.R.: A modification of the Pisano method for vanilmandelic acid using high pressure liquid chromatography. *Clin. Chim. Acta* 76: 317, 1977.
7. Bertani-Dziedzic, L., Krstulovic, A.M., Ciriello, S. and Gitlow, S.E.: Routine reverse-phase high performance liquid chromatographic measurement of urinary vanillylmandelic acid in patients with neural crest tumors. *J. Chromatogr.* 164:345, 1979.
8. Soldin, S.J. and Hill, J.G.: Simultaneous liquid chromatographic analysis for 4-hydroxy-3-methoxymandelic acid and 4-hydroxy-3-methoxyphenylacetic acid in urine. *Clin. Chem.* 26: 291, 1980.
9. Soldin, S.J. and Hill, J.G.: Liquid chromatographic analysis for urinary 4-hydroxy-3-methoxyphenylacetic acid and its use in investigation of neural crest tumors. *Clin. Chem.* 27: 502, 1981.
10. Moleman, P. and Borstok, J.J. M.: Determination of urinary vanillylmandelic acid by liquid chromatography with electro-chemical detection. *Clin.Chem.* 29: 878, 1983.
11. Fujita, K., Maruta, K., Ito S. and Nagatsu, T.: Urinary 4-hydroxy-3-methoxymandelic (vanillylmandelic) acid, 4-hydroxy-3-methoxyphenylacetic (homovanillic) acid, and 5-hydroxy-3-indoleacetic acid determined by liquid chromatography with electrochemical detection. *Clin. Chem.* 29: 876, 1983.
12. Rosano, T.G. and Brown, H.H.: Liquid chromatographic assay for urinary 3-methoxy-4-hydroxymandelic acid as measured by liquid chromatography, with on-line
13. Flood, J.G., Granger, M. and McComb, R.B.: Urinary 3-methoxy-4-hydroxymandelic acid as measured by liquid chromatography, with on-line post-column reaction. *Clin. Chem.* 25: 1234, 1979.
14. Manger, W.M. and Gifford, R.W., Jr.: *PHEOCHROMOCYTOMA* p.202, Springer-Verlag, N.Y., 1977.
15. Forman, D.T., Grayson, S.H. and Kirkpatrick, K.: Rapid determination of Vanilmandelic acid in urine by ion-exchange chromatography. *Clin. Chem.* 19:646, 1973.
16. Forman, D.T.: Measurement of Vanillylmandelic acid (VMA) by Column chromatography. *Clinical Pathology of Cancer of the Endocrine Glands and Target Organs*. Sunder, F.W., ed. Philadelphia, Institute for Clinical Sciences, Inc., p. 167, 1978.
17. Caraway, W.T. and Kammeyer, C.W.: Chemical interference by drugs and other substances with clinical laboratory test procedures. *Clin. Chim. Acta* 41:395, 1972.

## SYMBOLS USED WITH DEMEDITEC ASSAY'S

Symbol	English	Deutsch	Français	Español	Italiano
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	European Conformity	CE-Konformitätskennzeichnung	Conformité aux normes européennes	Conformidad europea	Conformità europea
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Numéro de catalogue	Número de catálogo	Numero di Catalogo
	Lot. No. / Batch code	Chargen-Nr.	Numéro de lot	Número de lote	Numero di lotto
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore
Content	Content	Inhalt	Conditionnement	Contenido	Contenuto
Volume/No.	Volume / No.	Volumen/Anzahl	Volume/Quantité	Volumen/Número	Volume/Quantità

Symbol	Portugues	Dansk	Svenska	Ελληνικά
	Consulte as instruções de utilização	Se brugsanvisning	Se bruksanvisningen	Εγχειρίδιο χρήστη
	Conformidade com as normas europeias	Europaeisk overensstemmelse	Europeisk överensstämmelse	Ευρωπαϊκή Συμμόρφωση
	Diagnóstico in vitro	In vitro diagnostik	Diagnostik in vitro	in vitro διαγνωστικό
				
	Catálogo n.º	Katalognummer	Katalog nummer	Αριθμός καταλόγου
	No do lote	Lot nummer	Batch-nummer	Αριθμός Παρτίδος
		Indeholder tilstrækkeligt til "n" test	Innehåller tillräckligt till "n" tester	Περιεχόμενο επαρκές για «n» εξετάσεις
	Temperatura de conservação	Opbevarings-temperatur	Förvaringstemperatur	Θερμοκρασία αποθήκευσης
	Prazo de validade	Udløbsdato	Bäst före datum	Ημερομηνία λήξης
	Fabricante	Producent	Tillverkare	Κατασκευαστής
Distributed by				
Content	Conteúdo	Indhold	Innehåll	Περιεχόμενο
Volume/No.	Volume/Número	Volumen/antal	Volym/antal	Όγκος/αριθ..