Product information

Information about other products is available at: www.demeditec.com

User's Manual

Estriol total ELISA

Enzyme Immunoassay for the quantitative determination of total estriol in human serum and plasma.

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1 INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Total Estriol concentration in human serum or plasma.

Total Estriol kit is intended for laboratory use only.

1.1 CLINICAL SIGNIFICANCE

Estriol (also oestriol) is one of the three main estrogens produced by the human body. It is only produced in significant amounts during pregnancy as it is made by the fetus.

During pregnancy the production of estriol depends on an intact maternal-placental-fetal unit. Fetal-placental production of estriol leads to a progressive rise in maternal circulating levels reaching a late-gestational peak several orders of magnitude greater than non-pregnant levels. In the maternal circulation, estriol undergoes rapid conjugation in the liver followed by urinary excretion with a half-life of ~20 minutes. Since normal estriol production depends on an intact maternal-placental-fetal circulation and functional fetal metabolism, maternal estriol levels have been used to monitor fetal status during pregnancy, particularly during the third trimester.

DHEA is produced by the adrenal cortex of the fetus, this is converted to estriol by the placenta.

If levels are abnormally low in a pregnant woman, this may indicate a problem with the development in the child.

Levels of estriol in non-pregnant women do not change much after menopause, and levels are not significantly different from levels in men.

2 PRINCIPLE

Total Estriol (antigen) in the sample competes with the antigenic estriol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti estriol coated on the microplate (solid phase). After incubation, the bound/free separation is performed by a simple solid-phase washing. Then the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue color that changes into yellow when the Stop Solution is added.

The colour intensity is inversely proportional to the Total Estriol concentration in the sample.

Total Estriol concentration in the sample is calculated based on a series of standards.

3 REAGENT, MATERIAL AND INSTRUMENTATION

3.1 Reagents and materials supplied in the kit

1. Total Estriol **Standards** (4 vials, 1 mL each)
2. Total Estriol **Control** (1 vial, 1 mL)
Concentration of control is Lot-specific and is indicated on Quality Control Report
3. **Incubation Buffer** (1 vial, 30 mL)
Phosphate buffer 50 mM pH 7.5; BSA 1 g/L, stabilisers
4. **Enzyme Conjugate** (1 vial, 1 mL)
Estriol conjugated with horseradish peroxidase (HRP)
5. **Coated Microplate** (1 breakable microplate)
Anti-Estriol antibody adsorbed on microplate
6. **TMB Substrate Solution** (1 vial, 15 mL)
H₂O₂.TMB 0.26 g/L, (avoid any skin contact)
7. **Stop Solution** (1 vial, 15 mL)
Acidic solution 0.3 N, (avoid any skin contact)

3.2 Reagents necessary not supplied

Distilled water.

3.3 Auxiliary materials and instrumentation

Automatic dispenser.

Microplate reader (450 nm)

Note

Store all reagents between 2 °C - 8 °C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close immediately after use.

Do not remove the adhesive sheet from the unused strips.

4 WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300 as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted acidic solution. Acidic solution is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Total Estriol from 2 ng/mL to 200 ng/mL.
- The clinical significance of Estriol determination can be invalidated if the patient was treated with natural or syntetic steroids.

5 PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 °C - 8 °C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22 °C - 28 °C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry dates printed on the labels of the box and of the vials must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemeic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6 PROCEDURE

6.1 Preparation of the Standard (S₁,S₂,S₃,S₄)

Before using, mix for 2 minutes.

The standards are ready to use and have the following concentration of Estriol:

| | S ₁ | S ₂ | S ₃ | S ₄ |
|-------|----------------|----------------|----------------|----------------|
| ng/mL | 2.0 | 20.0 | 80.0 | 200.0 |

Stable until the expiry date of the kit at 2 °C - 8 °C.

Once opened, the standards are stable for 6 months at 2 °C - 8 °C.

6.2 Preparation of diluted Conjugate

Prepare immediately before use.

Add 10 µL Conjugate (reagent 4) to 2.0 mL of Incubation Buffer (reagent 3).

Mix gently for 5 minutes, with a rotating mixer.

Stable for 3 hours at room temperature (22 °C – 28 °C).

6.3 Preparation of Sample

The determination of Total Estriol should be performed in human serum or plasma

Store samples at -20 °C if the determination is not performed on the same day of sample collection.

Avoid repetitive freezing and thawing of samples.

The Control is ready for use.

6.4 Procedure

Allow all reagents to reach room temperature (22 °C – 28 °C).

Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 °C - 8 °C.

To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.

As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for B₀, two for each point of the calibration curve (S₁-S₄), two for each Control, two for each sample, one for Blank.

| Reagent | B ₀ | Standard | Sample/ Control | Blank |
|--|----------------|----------|--------------------|--------|
| Incubation buffer | 20 µL | | | |
| Sample/Control | | | 20 µL | |
| Standard S ₁ -S ₄ | | 20 µL | | |
| Diluted conjugate | 200 µL | 200 µL | 200 µL | |
| Incubate at 37 °C for 1 hour. Remove the content from each well; wash the wells 2 times with 300 µL of distilled water. | | | | |
| TMB Substrate | 100 µL | 100 µL | 100 µL | 100 µL |
| Incubate at 22 °C – 28 °C for 15 minutes in the dark. | | | | |
| Stop Solution | 100 µL | 100 µL | 100 µL | 100 µL |
| Shake gently the microplate. Read the absorbance (E) at 450 nm against Blank. | | | | |

7 QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Total Estriol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8 RESULTS

8.1 Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the standard curve (S1-S4) and of each sample

8.2 Standard Curve

Plot the values of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points (e.g: Four Parameter Logistic).

8.3 Calculation of Results

Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in ng/mL.

9 REFERENCE VALUE

Serum concentrations of Estriol are included in the following ranges:

| weeks | Median | Range (ng/mL) |
|-------|--------|---------------|
| 17° | 18.0 | (10 - 27) |
| 18° | 25.9 | (14 - 51) |
| 19° | 39.5 | (26 - 52) |
| 20° | 40.0 | (27 - 53) |
| 21° | 45.6 | (24 - 66) |
| 22° | 39.2 | (25 - 58) |
| 23° | 56.1 | (27 - 70) |
| 24° | 56.3 | (28 - 75) |
| 25° | 64.3 | (29 - 84) |
| 26° | 68 | (41 - 105) |
| 27° | 57.4 | (41 - 110) |
| 28° | 78.0 | (38 - 127) |
| 29° | 87 | (45 - 146) |
| 30° | 75 | (45 - 160) |
| 31° | 88.0 | (50 - 170) |
| 32° | 90.5 | (46 - 175) |
| 33° | 100 | (60 - 180) |
| 34° | 105.6 | (60- 190) |
| 35° | 114.2 | (65 - 200) |
| 36° | 126.0 | (74 - 210) |
| 37° | 177.0 | (90 - 234) |
| 38° | 190.0 | (101 - 288) |
| 39° | 190.0 | (102 - 306) |
| 40° | 180.0 | (60 - 325) |
| 41° | 177.5 | (95 - 280) |

10 PERFORMANCE AND CHARACTERISTICS

10.1 Precision

10.1.1 Intra Assay Variation

Within run variation was determined by replicate measurements (16x) of two different control sera in one assay. The within assay variability is $\leq 9.7\%$.

10.1.2 Inter Assay Variation

Between run variation was determined by replicate measurements (12x) of three different control sera in different kit lots. The between assay variability is $\leq 10\%$.

10.2 Accuracy

The recovery of 10 – 40 – 100 ng/mL of Total Estriol added to two samples gave an average value (\pm SD) of $94.88\% \pm 4.47\%$ with reference to the original concentrations.

10.3 Sensitivity

The lowest detectable concentration of Total Estriol that can be distinguished from the zero standard is 0.22 ng/mL at the 95 % confidence limit.

10.4 Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

| | |
|---------------------------------|------------------------|
| Estriol | 100 % |
| 16 epi-estriol | 10.5 % |
| 15 α OH-estriol | 7.0 % |
| Estriol 3 Sulphate | 2.0 % |
| Estradiol | 0.1 % |
| 17 epi-estriol | $< 1 \times 10^{-2}\%$ |
| Estriol 3 α Glucuronate | $< 1 \times 10^{-2}\%$ |
| Estriol 16 α Glucuronate | $< 1 \times 10^{-2}\%$ |
| Estrone | $< 1 \times 10^{-4}\%$ |

10.5 Correlation with RIA

The total Estriol ELISA was compared to another commercially available Total Estriol assay. 32 serum samples were analysed.

The linear regression curve was calculated

$$y = 0.86x + 3.85$$

$$r^2 = 0.952$$

$$y = \text{Total Estriol ELISA (DE3717)}$$

$$x = \text{Total Estriol Adaltis RIA Kit}$$

11 WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

12 BIBLIOGRAPHY

1. Fischer-Rasmussen, W., et al Acta Obstet Gynecol. Scand. 60-417 - 420 (1981)
2. Truran, P.L., et al Clin. Chem. 28/12, 2393 (1982)
3. Vining, R. F., et al J. Clin. Endoc. Metab. 56, 454 (1983)
4. Bagger, P.V, et al Acta Obstet Gynecol Scand 60, 187 (1981)
5. Osterman, T.M, et al Clin. Chem. 25(5) 716 (1979)
6. Wisdom, G.B. Clin. Chem. 22 (8) 1243-1255 (1976)

13 TROUBLESHOOTING

ERRORS / POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed)

too high within-run CV%

- reagents and/or strips not pre-warmed to room temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)

too high between-run CV %

- incubation conditions not constant (time, temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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