Product information

Information about other products is available at: www.demeditec.com

User's Manual

Estrone ELISA

DE4174

96

Contents / Inhaltsverzeichnis / Contenuti / Contenido / Contenu

1	INTRODUCTION	4	1	INTRODUCCIÓN	24	
2	PRINCIPLE OF THE TEST	4	2	FUNDAMENTO DEL ENSAYO	24	
3	WARNINGS AND PRECAUTIONS.....	5	3	PRECAUCIONES.....	24	
4	REAGENTS.....	6	4	COMPONENTES DEL KIT.....	25	
5	SPECIMEN COLLECTION AND PREPARATION....	7	5	MUESTRAS	26	
6	ASSAY PROCEDURE.....	7	6	PROCEDIMIENTO DE ENSAYO	26	
7	EXPECTED NORMAL VALUES.....	9	7	VALORES ESPERADOS	28	
8	QUALITY CONTROL.....	9	8	CONTROL DE CALIDAD	28	
9	PERFORMANCE CHARACTERISTICS	9	9	CARACTERÍSTICAS DEL ENSAYO.....	28	
10	LIMITATIONS OF USE.....	10	10	LIMITACIONES DE USO	29	
11	LEGAL ASPECTS	11	11	ASPECTOS LEGALES	29	
12	REFERENCES / LITERATURE.....	11	12	REFERENCIAS / BIBLIOGRAFÍA.....	29	
1	EINLEITUNG	12	1	INTRODUCTION.....	30	
2	TESTPRINZIP	12	2	PRINCIPE DU TEST.....	30	
3	VORSICHTSMAßNAHMEN.....	12	3	PRECAUTIONS D'UTILISATION.....	30	
4	BESTANDTEILE DES KITS	13	4	COMPOSITION DU KIT	31	
5	PROBENVORBEREITUNG.....	14	5	ECHANTILLON	32	
6	TESTDURCHFÜHRUNG.....	14	6	RÉALISATION DU TEST	32	
7	ERWARTETE WERTE	16	7	VALEURS ATTENDUES.....	34	
8	QUALITÄTS-KONTROLLE.....	16	8	CONTROLE DE QUALITE	34	
9	ASSAY CHARACTERISTIKA	16	9	CARACTERISTIQUES DU TEST	34	
10	GRENZEN DES TESTS	17	10	LIMITES D'UTILISATION.....	35	
11	RECHTLICHE GRUNDLAGEN.....	17	11	ASPECTS LEGAUX.....	35	
12	REFERENZEN / LITERATUR	17	12	REFERENCES / BIBLIOGRAPHIE	35	
1	INTRODUZIONE	18	SYMBOLS USED WITH DEMEDITEC ELISAS			36
2	PRINCIPIO DEL TEST	18				
3	PRECAUZIONI	18				
4	COMPONENTI DEL KIT.....	19				
5	CAMPIONI	20				
6	ATTUAZIONE DEL TEST.....	20				
7	VALORI NORMALI	22				
8	CONTROLLO QUALITÀ	22				
9	CARATTERISTICHE DEL TEST	22				
10	LIMITAZIONE DEL TEST	23				
11	ASPETTI LEGALI	23				
12	BIBLIOGRAFIA.....	23				

1 INTRODUCTION

1.1 Intended Use

The **DEMEDIATEC Estrone ELISA** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of Estrone in serum and plasma

1.2 Summary and Explanation

Estrone (3-hydroxy-1,3,5 (10)-estratrien-17-one) is beside estradiol and estriol one of the three major naturally occurring estrogens. The estrogens are involved in the development of female sex organs and secondary sex characteristics. Bioassay data indicate that the estrogenic activity of estrone is considerably lower in comparison to estradiol (1). However, the physiological role of endogenous estrone is not well defined.

Estrone is produced primarily from androstenedione. In premenopausal women, more than 50% of the estrone is secreted by the ovary. In prepubertal children, men and postmenopausal women, the major portion of estrone is derived from peripheral tissue conversion (2). During the follicular phase of the menstrual cycle the estrone level increases with a clear peak around day 13. The peak is of short duration and by day 16 of the cycle levels will be low again. A second peak during the luteal phase occurs around day 21 of the cycle. If fertilization does not occur production of estrone decreases again.

These changes of estrone concentration are in parallel to that of estradiol (3). Until the 4. to 6. week of pregnancy, estrone originates primarily from maternal sources such as the ovaries, adrenals, or peripheral conversion thus remaining within the normal values (4). After the 6. to 10. week of pregnancy the values increase gradually due to placental secretion of estrone. After menopause, estrone levels do not decline as dramatically as estradiol levels. In postmenopausal women estrone is the major estrogen. In males the concentration of E1 has been reported to rise up with age inversely to that of 17-OH-progesterone (5). In premenopausal women excessive estrone levels can result from the conversion of large amounts of androstenedione produced in polycystic ovary syndrom (6) and ovarian tumors.

2 PRINCIPLE OF THE TEST

The DEMEDIATEC Estrone ELISA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. The microtiter wells are coated with streptavidin and biotinylated Estrone molecules. Endogenous Estrone of a patient sample competes with this immobilized Estrone for binding sites of a HRP-conjugated sheep antibody directed towards an antigenic site of estrone.

After incubation the unbound antibody conjugate is washed off. The amount of bound peroxidase conjugate is inversely proportional to the concentration of Estrone in the sample. After addition of the substrate solution, the intensity of colour developed is inversely proportional to the concentration of Estrone in the patient sample.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 1 N acidic solution. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Material Safety Data Sheets.
Material Safety Data Sheets for this product are available upon request directly from DEMEDITEC.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with biotinylated Estrone.
2. **Standard (Standard 0-6)**, 7 vials, 1 mL each, ready to use;
Concentrations: 0 – 13 – 40 – 100 – 220 – 500 - 1000 pg/mL
Conversion: pg/mL x 3.69 = pmol/L
Contain non-mercury preservative.
3. **Control Low & High**, 2 vials, 1 mL each, ready to use;
For control values and ranges please refer to vial label or QC-Datasheet.
Contain non-mercury preservative.
4. **Enzyme Conjugate**, 1 vial, 14 mL, ready to use,
Sheep anti Estrone antibody conjugated to horseradish peroxidase;
Contains non-mercury preservative.
5. **Substrate Solution**, 1 vial, 25 mL, ready to use,
Tetramethylbenzidine (TMB).
6. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 1 N acidic solution,
Avoid contact with the stop solution. It may cause skin irritations and burns.
7. **Wash Solution**, 1 vial, 30 mL (40X concentrated),
see „Preparation of Reagents“.

Note: Additional *Standard 0* for sample dilution is available upon request.

4.2 Materials required but not provided

- A microtiter plate calibrated reader (450 ± 10 nm)
- Calibrated variable precision micropipettes.
- Absorbent paper.
- Distilled or deionized water
- Timer
- Scale paper or semi-logarithmic graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again. Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

Wash Solution

Add deionized water to the 40X concentrated Wash Solution. Dilute 30 mL of concentrated *Wash Solution* with 1170 mL deionized water to a final volume of 1200 mL. *The diluted Wash Solution is stable for 2 weeks at room temperature.*

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheet.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DEMEDITEC has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA-, heparin- or citrate plasma) can be used in this assay. Do not use haemolytic, icteric or lipaemic specimens. Please note: Samples containing sodium azide should not be used in the assay.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 5 days at 2 °C to 8 °C prior to assaying. Specimens held for a longer time should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Standard 0* and reassayed as described in Assay Procedure. For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 µL sample + 90 µL *Standard 0* (mix thoroughly)
- b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Standard 0* (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **100 µL** of each **Standard, Control** and **samples** with new disposable tips into appropriate wells.
3. Dispense **100 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes** at room temperature.
5. Briskly shake out the contents of the wells.
Rinse the wells **4 times** with diluted *Wash Solution* (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.

Important note:

The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!

6. Add **200 µL** of **Substrate Solution** to each well.
7. Incubate for **15 minutes** at room temperature.
8. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
9. Determine the absorbance (OD) of each well at **450 ± 10 nm** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using scale paper or semi-logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4 Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.)
Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 1000 pg/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Standard 0 (0 pg/mL)	2.43
Standard 1 (13 pg/mL)	2.04
Standard 2 (40 pg/mL)	1.69
Standard 3 (100 pg/mL)	1.16
Standard 4 (220 pg/mL)	0.69
Standard 5 (500 pg/mL)	0.36
Standard 6 (1000 pg/mL)	0.18

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently normal healthy adults, using the DEMEDITEC Estrone ELISA the following values are observed:

Population	Valid N	5% Percentile (pg/mL)	95% Percentile (pg/mL)
Males	40	42.28	127.15
Females	64	32.5	154.85

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above mentioned items without finding any error contact your distributor or DEMEDITEC directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 2.21 pg/mL –1000 pg/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross reactivity of the assay:

Compound	Cross reactivity %
Progesterone	0.0
Testosterone	0.0
Estradiol	2.5
Free Estriol	2.1

9.3 Sensitivity

The analytical sensitivity of the DEMEDITEC ELISA was calculated by subtracting 2 standard deviations from the mean of 20 replicate analyses of the Standard 0 (S0) and was found to be 2.21 pg/mL.

9.4 Reproducibility

9.4.1 Intra Assay

The within assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	20	23.4	6.5
2	20	34.0	9.4
3	20	24.9	9.5

9.4.2 Inter Assay

The between assay variability is shown below:

Sample	n	Mean (pg/mL)	CV (%)
1	30	31.4	12.8
2	30	23.7	13.6
3	30	14.7	13.0
4	30	22.9	14.8

9.5 Recovery

Samples have been spiked by adding Estrone solutions with known concentrations

	Sample 1	Sample 2	Sample 3
Concentration (pg/mL)	260.7	51.0	90.5
Average Recovery [%]	95.7	97.5	101.0
Range of Recovery [%]	from	85.2	85.9
	to	101.5	114.6

9.6 Linearity

	Sample 1	Sample 2	Sample 3
Concentration (pg/mL)	340.0	572.0	462.0
Average Recovery [%]	102.9	101.4	104.0
Range of Recovery [%]	from	96.9	92.6
	to	112.4	110.5

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of Estrone in a sample.

10.3 High-Dose-Hook Effect

No hook effect was observed in this test.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DEMEDITEC.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12 REFERENCES / LITERATURE

1. Resnik R, Killam AP, Battaglia FC et al: The stimulation of uterine blood flow by various estrogens. *Endocrinology* 94:1192, 1974.
2. Fayman C, Winter JSD, Reyes FI. Patterns of gonadotropins and gonadal steroids throughout life. *Clin. Obstet. Gynecol.* 3: 467-483, 1976.
3. Baird DT, Fraser IS. Blood production and ovarian secretion rates of estradiol-17 β and estrone in women throughout the menstrual cycle. *J Clin Endocrinol. Metab* 38: 1009-1017. 1974
4. Lindbert BS, Johansson EDB, Nilsson BA: Plasma levels of non-conjugated oestrone, oestradiol-17b and oestriol during uncomplicated pregnancy. *Acta Obstet Gynecol Scand* 32:21, 1974.
5. Drafta D, Schindler AE, Stroe EW, Neacsu E. Age-related changes of plasma steroids in normal adult males. *J. Steroid Biochem.* 17: 683-687, 1982.
6. DeVane GW, Czekala NM, Judd HL, Yen SSC. Circulating gonadotropins, estrogens, and androgens in polycystic ovarian disease. *Am J Obstet Gynecol* 1975; 121:496.

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