# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

## User's Manual

# m/rIGFBP-2 ELISA

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## PACKAGE INSERT ENGLISH

### TECHNICAL FEATURES

- Highly specific and sensitive assay for quantitative detection of IGFBP-2 in mouse and rat serum
- Recombinant mouse IGFBP-2 as standard
- Uses antibodies against complete mouse and rat IGFBP-2
- No sample extraction is required
- Detection limit: 0.01 ng/ml

### INTRODUCTION

Insulin-like growth factors (IGFs) regulate the proliferation, differentiation, apoptosis, cell adhesion and metabolism in various tissues and cell types. The IGF-I, which is produced mainly in liver under the influence of Growth Hormone (GH), regulates as hormone the linear growth of the bones and the process of sexual maturity, while IGF-II is mainly a growth factor of foetal tissue (11-13). The biological actions of IGF over the IGF-Type-I receptor are modulated variably through the IGF binding proteins (IGFBP-1 to-6) (14). IGFBP-2 is, after IGFBP-3, the second most frequent IGFBP in the human blood. IGFs, especially tumor typical pro-IGF-forms and hormones regulate the expression of IGFBP-2, GH effect is thereby inhibiting. At cellular level IGFBP-2 seems to stimulate the proliferation and dissemination of solid tumors via an IGF-independent mechanism (15,16).

### PHYSIOLOGICAL MEANING

IGFBP-2 is an unglycosylated polypeptide of 31.3 kDa, which forms binary IGF-complexes and shows no circadian rhythm in the circulation.

The serum concentration of IGFBP-2 increases in fasting, after major surgery and after trauma, but the increasing of the concentration is most intensive in malignant diseases. The correlation of the IGFBP-2 level to the degree of progression is a striking feature in various tumor types as is the normalization of the IGFBP-serum levels after remission (5-10). During the GH-therapy, e.g. in short stature and in GH-abuse (doping) the IGFBP-2 level decreases. In Trisomy 18 IGFBP-2 in maternal serum is decreased and IGFBP-1 is increased; therefore the ratio IGFBP-2 /IGFBP-1 is a marker for this chromosome abnormality (17).

Transgenic organisms are a good opportunity to investigate the function of genes or proteins. The mouse or rat model is a well-suited system for investigation of the relevance of IGFBP-2 in physiological and pathological processes. Over expression of the IGFBP-2 gene in mice results in a weight reduction of 30% in spleen and moderately reduced weight in other organs (18). Effects of IGFBP-2 on the organism can be compensated through the modified expression of other IGF-Binding proteins.

Especially in tumor biology the mouse and rat systems enable investigation of the systemic relevance of IGFBP-2. IGFBP-2 influences tumor cells as it induces catalase activity in adrenocortical cells (19). Furthermore IGFBP-2 interacts with tumor cells via its RGD-amino acid sequence and seems to stimulate cell invasion of glioma cells (20).

## INTENDED USE

This IGFBP-2 Enzyme Immunoassay-Kit is suited for quantitative determination of IGFBP-2 in mouse and rat serum for scientific purposes.

## METHODOLOGY

### Assay Characteristics and Validation

The ELISA for IGFBP-2 utilizes two different specific high affinity polyclonal antibodies for this protein. The ELISA recognizes quantitatively mouse and rat IGFBP-2 and is unaffected by an excess of IGF-I or IGF-II levels. Related molecules such as IGFBP-3 show no cross-reactions in the assay.

The standards are prepared of recombinant mouse-IGFBP-2 in the range of 0.03125 to 2 ng/ml.

The theoretical sensitivity of the assay is approx. 0.01 ng/ml (2 x SD of zero standard). Intra-assay and inter-assay variation coefficients were found both < 10%. Exemplary determinations are shown in the tables 1 and 2.

**Table 1 :** Inter-Assay-Variation

<b>Sample 1</b>	24.9 ng/ml	n= 6	CV = 6.4 %
<b>Sample 2</b>	105.7 ng/ml	n= 6	CV = 2.4 %
<b>Sample 3</b>	171.0 ng/ml	n= 8	CV = 5.3 %

**Table 2:** Intra-Assay-Variation

<b>Sample</b>	67.87 ng/ml	n = 9	CV = 4.6 %
---------------	-------------	-------	------------

Dilution of samples has been found over a wide range with very good linearity (see table 3).

**Table 3:** Linearity of the sample dilution:

<b>Dilution</b>	<b>Mouse Serum (ng/ml)</b>	<b>Rat Serum (ng/ml)</b>
1:20	>max.	23.06
1:40	>max.	23.81
1:80	109.08	24.10
1:160	118.30	25.73
1:320	125.42	27.73
1:640	127.48	32.14

## Calibration

The assay has been calibrated against the recombinant Mouse-IGFBP-2 of R&D Systems Inc. (Minneapolis, USA; [www.rndsystems.com](http://www.rndsystems.com)).

## Expectation Values

Several commercially available mouse and rat sera have been tested for their IGFBP-2 concentrations, following results were obtained :

	n	Median	min.	max.
Mouse Sera	5	<b>81.0 ng/ml</b>	59.7 ng/ml	105.7 ng/ml
Rat Sera	4	<b>24.2 ng/ml</b>	10.7 ng/ml	38.1 ng/ml

Significant variations of serum values depending on the individual animal or the respective strain or mutant are likely, prior verification is recommended.

## Sample Preparation and Storage

The samples should be stored frozen until measurement. IGFBP-2 levels are influenced by improper handling or storage and do not remain stable over several days at elevated temperatures. Store undiluted samples frozen in a tightly closed plastic vial. **Repeated freezing and thawing of serum/plasma should be avoided**, it seems to have a measurable effect on IGFBP-2 levels.

The high sensitivity of the assay allows measurement of IGFBP-2 in small sample volumes, which is limited by pipetting accuracy rather than the amount of IGFBP-2.

Serum samples should be diluted prior to measurement 1:20 – 1:500-fold with **Dilution Buffer VP**, depending on the expected values (see chapter Expectation Values). In general a dilution of 1:100 should be appropriate ( the recommended minimal essential sample volume is: 10 µl serum).

Sample extraction is not required.

### Suggestion for dilution protocol (double determination):

Mix 10 µl serum manually or with the aid of a dilutor with 990µl **Dilution Buffer VP** (1:100), or, more simple for larger series with 1000 µl Dilution Buffer VP (1:101).

If sample size is limiting, a minimum of 2.5 µl sample might be used alternatively, dilution in 250 µl VP yields a dilution of 1:101 (care should be taken to accuracy of pipetting such low volumes !).

Use 2 x 100 µl of this dilution in the assay.

**REAGENTS PROVIDED**

1)	MTP	<b>Microtiter Plate</b> with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with polyclonal anti-IGFBP-2 antibodies and packed in a laminate bag.
2)	CAL	<b>Standards A-G</b> , 1 ml, lyophilized: Contain recombinant mouse IGFBP-2: Standard values are between 31.25 - 2000 pg/ml (31.25; 62.5; 125; 250; 500; 1000; 2000 pg/ml) IGFBP-2 and have to be reconstituted with <b>1 ml Dilution Buffer VP</b> each.
3)	DILU	<b>Dilution Buffer VP, 120 ml</b> , ready for use.
4)	Control	<b>Control Serum KS, 100 µl</b> , lyophilized: Contains mouse serum and has to be reconstituted with <b>100 µl Dilution Buffer VP</b> . The exact concentration is given on the vial label.
5)	Ab	<b>Antibody Conjugate AK</b> , 120 µl, 100fold concentrated: Contains biotinylated anti-IGFBP-2 antibody and has to be diluted immediately before use <b>1:100 with Dilution Buffer VP</b> .
6)	CONJ	<b>Enzyme Conjugate EK</b> , 120 µl, 100fold concentrated: Contains HRP-labelled Streptavidin and has to be diluted immediately before use <b>1:100 with Dilution Buffer VP</b> .
7)	WASHBUF 20x	<b>Washing Buffer (WP)</b> , 50 ml, 20 X concentrated solution. <b>Washing Buffer (WP)</b> has to be diluted <b>1:20</b> with <b>distilled or demineralised water</b> before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Attention: After dilution the Washing Buffer is only limited stable, dilute only according to requirements.
8)	SUBST	<b>Substrate (S)</b> , 12 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H <sub>2</sub> O <sub>2</sub> Tetramethylbendidine.
9)	H <sub>2</sub> SO <sub>4</sub>	<b>Stopping Solution (SL)</b> , 12 ml, ready for use, 0.2 M sulphuric acid, Caution acid!
10)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

**MATERIALS REQUIRED BUT NOT PROVIDED**

- Precision pipettes (100 and 200µl) Micropipettes and multichannel pipettes with disposable plastic tips
- Distilled or Deionized water for dilution of the Washing Buffer (WP)
- Vortex-mixer
- Device to aspirate the standards and the samples from the wells (recommended because of the potential danger of infection by human samples)
- Timer (120 min. range)
- Reservoirs (disposable)
- Plate washer and plate shaker (recommended)
- Calibrated Micro plate reader ("ELISA-Reader") with filter for 450 and 620nm (or ≥590 nm)
- Foil welding device for laminate bags (recommended)

**TECHNICAL RECOMMENDATIONS**

**Bring all reagents to room temperature (20 - 25°C) before use.**

For the **reconstitution** of the lyophilised components (**Standards A - G** and **Control Serum KS**) the kit **Dilution Buffer VP** should be used. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

**Incubation at room temperature means: Incubation at 20 - 25 °C.**

## STORAGE CONDITIONS

The microtiter plate wells and all undiluted reagents are stable until the expiry date if stored in the dark at 2-8°C.

Store the unused seal strips and microtiter wells together with the desiccant at 2° to 8°C.

The Substrate Solution (S), stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbencidine, is photosensitive – store and incubate in the dark.

Reconstituted components (Standards A – G and Control Serum KS) should be stored at -20°C (or below) for up to 2 months. **Avoid repeated freeze-thaw cycles.**

In case you plan to perform multiple independent m/r-IGFBP-2 determinations over a longer period with one kit, you should aliquote the components prior to freezing into suitable smaller volumes. This is strongly recommended. The 1:20 diluted **Washing Buffer WP** is only limited stable. Please dilute only according to requirements.

## WARNINGS AND PRECAUTIONS

**For in-vitro diagnostic use only. For professional use only.**

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to **room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming. **Temperature WILL affect the absorbance** readings of the assay. However, values for the patient samples will not be affected.

Do not mix reagents of different lots. Do not use expired reagents.

The microplate contains snap-off strips. Unused wells must be stored at 2 - 8°C in the sealed foil pouch and used in the frame provided.

Caution: This kit contains material of animal origin, therefore all components and specimens should be treated as potentially infectious.

Following components contain **0.01% 2-Methyl-4-isothiazolin-3-one Solution** as preservative: AK, EK, VP

R34	Irritating to eyes and skin
R43	Sensibilisation through skin contact possible
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S36/37	Wear suitable protective clothing and gloves
S45	In case of accident or if you feel unwell seek medical advice

Following components contain **0.01% (w/w) 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-isothiazol-3-one** as preservative: AK, EK, VP, WP

R36/38	Irritating to eyes and skin
R43	Sensibilisation through skin contact possible
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1	After contact with skin, wash immediately with plenty of water

### Stop solution contains 0.2 M Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>)

R36/38	Irritating to eyes and skin
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1	After contact with skin, wash immediately with plenty of water
S36/37	Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

**TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.** Store and Incubate in tge dark.

R20/21/R22	Harmful by inhalation, in contact with skin and if swallowed
R36/37/38	Irritating to eyes, respiratory system and skin
S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1	After contact with skin, wash immediately with plenty of water
S36/37	Wear suitable protective clothing and gloves

#### General first aid procedures:

- Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.
- Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.
- Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.
- Do not eat, drink or smoke in these areas.
- Never pipette the materials with the mouth.
- Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

### ASSAY PROCEDURE

NOTES: All determinations (Standards, Control Serum and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Standards, Control Serum and the samples should be pipette as fast as possible (e.g., <15 minutes). To avoid distortions due to differences in incubation times, **Antibody- Conjugate AK, Enzyme Conjugate** as well as the following **Substrate Solution S** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution SL** should be added to the plate in the same order as the Substrate Solution.

- 1) add **100 µl Dilution Buffer VP** in wells A1/2 (blank) and
- 2) pipette in positions B1/2           **100 µl Standard A,**  
     pipette in positions C1/2           **100 µl Standard B,**  
     pipette in positions D1/2           **100 µl Standard C,**  
     pipette in positions E1/2           **100 µl Standard D,**  
     pipette in positions F1/2           **100 µl Standard E,**  
     pipette in positions G1/2           **100 µl Standard F,**  
     pipette in positions H1/2           **100 µl Standard G.**

To control correct accomplishment **100 µl** of the **(1:100)** diluted **Control Serum KS** can be pipetted in positions A3/4.

Pipette **100 µl** of the **diluted sample** in the rest of the wells, according to requirements.

- 3) Cover the wells with sealing tape and incubate the plate for **1 hour** at **room temperature** (shake at 350 rpm).
- 4) After incubation aspirate the contents of the wells and wash the wells 3 times with **250 µl** of **Washing Buffer WP** / well respectively.
- 5) Following the last washing step pipette **100 µl** of the of the **(1:100)** diluted **Antibody Conjugate AK** in each well, and incubate **1 hour** at **room temperature** (shake at 350 rpm).
- 6) After incubation wash the wells 3 times with **Washing Buffer WP** as described above.
- 7) Following the last washing step pipette **100 µl** of the **(1:100)** diluted **Enzyme Conjugate EK** in each well, and incubate **30 min** at **room temperature** (shake at 350 rpm).
- 8) After incubation wash the wells 3 times with **Washing Buffer WP** as described above.
- 9) Pipette **100 µl** of the **TMB-Substrate Solution S** in each well.
- 10) Incubate the plate for **30 minutes** in the dark at **room temperature**.
- 11) Stop the reaction by adding **100 µl** of **Stopping Solution SL** to all wells.  
     Measure the absorbance within **30 minutes** at **450 nm** (reference filter: 620 nm).

## EVALUATION OF RESULTS

### Establishing the Standard Curve

The standards provided contain the following concentrations of recombinant mIGFBP-2 :

Standard	A	B	C	D	E	F	G
pg/ml	31.25	62.5	125	250	500	1000	2000
ng/ml	0.03125	0.0625	0.125	0.25	0.5	1	2

- 1) Calculate the mean absorbance (MA) value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance (MA) of the blank from the mean absorbancies of all other values.
- 3) Plot the **Standard** concentrations **A-G** on the x-axis versus the mean value of the absorbancies of the **Standards** on the y-axis. By using the mean absorbancies of the samples herewith the sample concentrations can be received.
- 4) Recommendation: Calculation of standard curve and sample concentrations should be done by using a computer programme, because the standard curve is in general best described by a non-linear regression or a higher-grade polynom or four parametric (4PL) curve fits.
- 5) The m/rIGFBP-2 concentration of the samples can be calculated with the standard curve equation and by multiplication with the respective dilution factor.

**LITERATUR / LITERATURE**

1. Elmlinger MW, Wimmer K, Biemer E, Blum WF, Ranke MB, Dannecker GE, (1996) Insulin-like growth factor binding protein 2 is differentially expressed in leukaemic T-and B-cell lines. *Growth Regulation* 6, 152-157
2. Hoeflich A, Nedbal S, Blum WF, Erhard M, Lahm H, Brem G, Kolb HJ, Wanke R, Wolf E. (2001) Growth inhibition in giant growth hormone transgenic mice by overexpression of insulin-like growth factor-binding protein-2. *Endocrinology* 142:1889-98
3. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Hormone Research* 54: 60-68
4. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001). Relevance of IGF-I, IGFBP-3, and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Hormone Research* 55: 115-124
5. Crofton PM, Ahmed SF, Wade JC, Elmlinger MW, Ranke MB, Kelnar CJH, Wallace WHB (1999) Effects of a third intensification block of chemotherapy on bone and collagen turnover, insulin-like growth factor I, its binding proteins and short term growth in children with acute lymphoblastic leukemia. *European Journal of Cancer* 35: 960-967
6. Muller HL, Oh Y, Lehrnbecher T, Blum WF, Rosenfeld RG. (1994) Insulin-like growth factor-binding protein-2 concentrations in cerebrospinal fluid and serum of children with malignant solid tumors or acute leukemia. *J Clin Endocrinol Metab.* 79: 428-34
7. Martin W Elmlinger, Martin H Deininger, Burkhardt S Schuett, Richard Meyermann, Frank Duffner, Ernst H Grote, and Michael B Ranke (2001) In vivo expression of the insulin-like growth factor binding protein-2 in human gliomas increases with the tumor grade. *Endocrinology* 142: 1652-1658
8. Cohen P, Peehl DM, Stamey TA, Wilson KF, Clemmons DR, Rosenfeld RG (1993) Elevated levels of insulin-like growth factor-binding protein-2 in the serum of prostate cancer patients. *J Clin Endocrinol Metab.* 76: 1031-1035
9. Boulle N, Baudin E, Gicquel C, Logie A, Bertherat J, Penfornis A, Bertagna X, Luton JP, Schlumberger M, Le Bouc Y. (2001) Evaluation of plasma insulin-like growth factor binding protein-2 as a marker for adrenocortical tumors. *Eur J Endocrinol.* 144: 29-36
10. Ranke MB, Maier KP, Schweizer R, Stadler B, Schleicher s, Elmlinger MW, Flehmig B (2003), Pilot study of elevated levels of insulin-like growth factor-binding protein-2 as indicators of hepatocellular carcinoma. *Horm Res* 60:174-180
11. Rosenfeld RG, Roberts CT Jr. (eds.) (1999) *The IGF system: Contemporary Endocrinology Series; Humana Press*
12. Jones JL, Clemmons DR. (1995) Insulin-like growth factors and their binding proteins: biological actions. *Endocr Rev.* 16: 3-34
13. Chard T (1994) Insulin-like growth factors and their binding proteins in normal and abnormal foetal growth. *Growth Reg* 4: 91-100
14. Ranke MB, Elmlinger MW (1997) Functional role of insulin-like growth factor binding proteins. *Hormone Research* 48 (suppl 4): 9-15
15. Elmlinger MW, Bell M, Schütt, BS, Langkamp M, Kutoh E, Ranke MB (2001) Transactivation of the IGFBP-2 promoter in human tumor cell lines. *Molecular and Cellular Endocrinology* 175: 211-218
16. Hoeflich A, Reisinger R, Lahm H, Kiess W, Blum WF, Kolb HJ, Weber MM, Wolf E. (2001) Insulin-like growth factor-binding protein 2 in tumorigenesis: protector or promoter? *Cancer Res.* 15: 8601-8610
17. Miell JP, Langford KS, Jones JS, Noble P, Westwood M, White A, Nicolaidis KH. (1997) The maternal insulin-like growth factor (IGF) and IGF-binding protein response to trisomic pregnancy during the first trimester: a possible diagnostic tool for trisomy 18 pregnancies. *J Clin Endocrinol Metab.* 82: 287-292
18. Hoeflich A, Wu M, Mohan S, Foll J, Wanke R, Froehlich T, Arnold GJ, Lahm H, Kolb HJ, Wolf E. (1999) Overexpression of insulin-like growth factor-binding protein-2 in transgenic mice reduces postnatal body weight gain. *Endocrinology* 140(12):5488-5496.
19. Hoeflich A, Fettscher O, Preta G, Lahm H, Kolb HJ, Wolf E, Weber MM (2003) Increased activity of catalase in tumor cells overexpressing IGFBP-2. *Horm Metab Res* 35(11-12):816-821.
20. Song SW, Fuller GN, Khan A, Kong S, Shen W, Taylor E, Ramdas L, Lang FF, Zhang W. (2003) lip45, an insulin-like growth factor binding protein 2 (IGFBP-2) binding protein, antagonizes IGFBP-2 stimulation of glioma cell invasion. *Proc Natl Acad Sci USA* 100(24):13970-13975.

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## Summary of the Assay

Reagent preparation:	Reconstitution:	Dilution:
Standards A-G	in 1 ml Dilution Buffer VP	
Control Serum KS	in 100 µl Dilution Buffer VP	1:100 with Dilution Buffer VP
Antibody Conjugate AK		1:100 with Dilution Buffer VP
Enzyme Conjugate EK		1:100 with Dilution Buffer VP
Washing Buffer WP		1:20 with Aqua. dest. (e.g., add the complete contents of the flask (50 ml) into a graduated flask and fill with A.dest. to 1000 ml).
<b>Sample dilution:</b> e.g. 1:100 (e.g. Mix 10 µl Serum with 990 µl Dilution Buffer VP).		

## Assay Procedure for double determination

Pipette	Reagents	Well positions
100 µl	Dilution Buffer VP (Blank)	A1/2
100 µl	Standard A (31.25 pg/ml)	B1/2
100 µl	Standard B (62.5 pg/ml)	C1/2
100 µl	Standard C (125 pg/ml)	D1/2
100 µl	Standard D (250 pg/ml)	E1/2
100 µl	Standard E (500 pg/ml)	F1/2
100 µl	Standard F (1000 pg/ml)	G1/2
100 µl	Standard G (2000 pg/ml)	H1/2
100 µl	Control Serum KS	A3/4
100 µl	Sample dilution	following wells
Cover the wells with the sealing tape.		
<b>Incubation: 1 h</b> at RT, ≥350 rpm		
3x 250 µl	Aspirate the contents of the wells and <b>wash</b> 3x with <b>250 µl Wash Buffer WP</b>	each well
100 µl	1:100 diluted Antibody Conjugate AK	each well
<b>Incubation: 1 h</b> at RT, ≥350 rpm		
3x 250 µl	Aspirate the contents of the wells and <b>wash</b> 3x with <b>250 µl Wash Buffer WP</b>	each well
100 µl	1:100 diluted Enzyme Conjugate EK	each well
<b>Incubation: 30 min</b> at RT, ≥350 rpm		
3x 250 µl	Aspirate the contents of the wells and <b>wash</b> 3x with <b>250 µl Wash Buffer WP</b>	each well
100 µl	Substrate Solution S	each well
<b>Incubation: 30 min</b> in the <b>dark</b> at RT		
100 µl	Stop Solution SL	each well
Measure the absorbance within 30 min at <b>450 nm</b> with <b>620 nm</b> as reference wavelength.		

REF DEE008

## International Test description

<b>CAL</b> <b>A-G</b>	A-G	<b>Rec in</b> 1 ml <b>VP</b>	
<b>Control</b>	KS	<b>Rec in</b> 100 µl <b>VP</b>	1:100 <b>DILU</b> <b>VP</b>
<b>Ab</b>	AK		1:100 <b>DILU</b> <b>VP</b>
<b>CONJ</b>	EK		1:100 <b>DILU</b> <b>VP</b>
<b>WASHBUF</b> <b>20x</b>	WP		1:20 <b>DILU</b> <b>A. dest.</b>
<b>SPE</b>			1:100 <b>DILU</b> <b>VP</b>
<b>°C</b> 20-28 °C			

100 µl	<b>VP</b>	A1/2
100 µl	<b>CAL</b> <b>A</b> (31.25 pg/ml)	B1/2
100 µl	<b>CAL</b> <b>B</b> (62.5 pg/ml)	C1/2
100 µl	<b>CAL</b> <b>C</b> (125 pg/ml)	D1/2
100 µl	<b>CAL</b> <b>D</b> (250 pg/ml)	E1/2
100 µl	<b>CAL</b> <b>E</b> (500 pg/ml)	F1/2
100 µl	<b>CAL</b> <b>F</b> (1000 pg/ml)	G1/2
100 µl	<b>CAL</b> <b>G</b> (2000 pg/ml)	H1/2
100 µl	<b>CONTROL</b> <b>KS</b> 1:100 <b>DILU</b> <b>VP</b>	A3/4
100 µl	<b>SPE</b> 1:100 <b>DILU</b> <b>VP</b>	B3/4...
<b>TAPE</b>		

1 h **°C** 20-25 ≥ 350 rpm

3x 250 µl	3x <b>WASHBUF</b> <b>WP</b>
100 µl	<b>Ab</b> <b>AK</b>
<b>TAPE</b>	

1 h **°C** 20-25 ≥ 350 rpm

3x 250 µl	3x <b>WASHBUF</b> <b>WP</b>
100 µl	<b>CONJ</b> <b>EK</b>
<b>TAPE</b>	

0.5 h **°C** 20-25 ≥ 350 rpm

3x 250 µl	3x <b>WASHBUF</b> <b>WP</b>
100 µl	<b>SUBST</b> <b>TMB</b> <b>S</b>

0.5 h **°C** 20-25

100 µl	<b>H<sub>2</sub>SO<sub>4</sub></b> <b>SL</b>
<b>MEASURE</b>	