# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

User's Manual

# IGF-I ELISA

DEE020

96 wells

EN/ DE/ FR/ IT/ ES/ PT/ NL/ DK/ SE/ PL/ HU/ SK/ CZ/ BG/ EE/ GR/ RO/ SL/ FI

**Symbols/ Symbole /Symboles/ Simboli/ Símbolos/ Símbolos/ Symbolen/ Symboler/ Symboler/ Symbole/ Szimbólumok/ Symboly/ Symboly/ Символи/ Sümbolid/ Σύμβολα/ Simboluri/ Simboli/ Symbolit**

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## PACKAGE INSERT ENGLISH

### TECHNICAL FEATURES+APPLICATIONS

- High Specificity for IGF-I
- Correct measurement of IGF-I in **non-extracted** samples. No physical separation of IGF-I from IGF-binding proteins required
- Elimination of **interference by IGF-binding proteins** through excess IGF-II
- Calibration against the WHO International Reference Standard preparation of IGF-I,  
**WHO NIBSC 02/254**
- 98,7 % recovery of recombinant IGF-I leads to correct absolute values
- Precise measurement of very low IGF-I levels: **sensitivity of 0.09 ng/ml**
- Small sample volume requirement, thus ideal for young patients.
- For clinical application: **only one dilution** necessary
- Inter- and Intra-assay variance: 6.8 and 6.7%


### Intended Use

Measurement of human IGF-I in human serum and plasma samples.

### INTRODUCTION

Insulin-like growth factors (IGF) I and II play a pivotal role in regulating the proliferation, differentiation and specific functions of many cell types (1-3). IGF-I is identical with Somatomedin C (Sm-C) (4) and has a molecular weight of 7649 daltons (5). Its major regulators are growth hormone (GH) and nutrition (6), although its production in specific tissues is affected by a multitude of tropic hormones and other peptide growth factors. In contrast to many other peptide hormones, IGFs are avidly bound to specific binding proteins (IGFBP). The seven classes of IGFBPs which are known at present (7,8,22) either bind IGF-I and IGF-II with similar affinities or show a preference for IGF-II (9,10).

A major problem of IGF-I measurement results from the interference of IGFBPs in the assay. Direct determinations in untreated serum samples (11) give false values because of the extremely slow dissociation of the IGF-I/IGFBP-3 complexes during the assay incubation. Depending on the ratio IGF-I to IGFBP in the sample interference comes up (see example Figure 1):

**Figure 1.** Interference of IGFBP in IGF-I measurements. Known concentrations of IGF-I were assayed in the presence of 0.5 ng (left) or 5 ng (right) hIGFBP-3 by a conventional (□) and by the IGFBP-blocked assay (\*).

Therefore, various techniques were applied to physically separate IGF-I from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2,12,13). These techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries. The most widely used method is the acid-ethanol extraction (13,14) with a recovery of only 70-80 % of IGFBP-bound IGF-I as a result of co-precipitation. The absolute results of such an extraction are therefore false low (15). The extraction removes the IGFBPs only insufficiently and leads to reduction in sensitivity of the assay due to pre-dilution of the samples by the extraction procedure. Furthermore, the remaining IGFBP may still interfere in the assay. In addition, the acid-ethanol extraction is ineffective in specimens other than serum or plasma (e.g. cell culture media), in which determination of IGF-I is already difficult enough due to the fact that IGFBPs are frequently present at large excess. To avoid these difficulties, an uncomplicated assay was developed, in which special sample preparation is not required before measurement.

### Clinical Significance

There are apart from GH, a number of variables that influence serum IGF-I. Decreased levels are found in states of malnutrition/malabsorption, hypothyroidism, liver disease, untreated diabetes mellitus, chronic inflammatory disease (1,6), malignant disease or polytrauma. High levels, on the other hand, are likely to be present in precocious puberty or obesity. Crucially important to the correct interpretation of IGF-I measurements is the relationship between age and IGF-I levels. It is certainly inadequate to use a common cut-off point to define "normal" levels for all age groups, particularly in children and adolescents.

Due to its GH-dependence, determination of serum IGF-I was shown to be a useful tool in diagnosis of growth disorders, especially with regard to GH deficiency (GHD) or acromegaly (6,16-19,23,24). The major advantage of IGF-I determination compared to GH determination is its stable circadian concentration; therefore a single measurement is sufficient. Hence IGF-I determination should be the first in a series of laboratory test. Clearly normal levels would then rule out disturbances of the GH-IGF-

I-axis. Low levels, i.e. close to or below the age-related 5th percentile would indicate the necessity of further diagnostic efforts. Subnormal levels of IGF-I would be evidence for reduced GH secretion, if other causes of low serum IGF-I (e.g. malnutrition or impaired liver function) can be ruled out. For differentiation of healthy short children without GH deficiency and children with "classical" GH deficiency, the 0.1st percentile proved to be an appropriate cut-off point, especially after the age of eight. However, IGF-I levels of short children not suffering from GHD may nevertheless lay between the 0.1st and 5th percentile (19). In contrast, acromegaly is characterized by pathologically elevated IGF-I levels, which apparently reflect the severity of the disease better than GH-levels (17,18,20).

### **Scientific Use**


IGF-I is present in low concentrations in various body fluids and in conditioned cell culture media of many cell lines. However, the determination of IGF-I in these specimens is particularly difficult due to the presence of IGFBPs usually in excessive amounts. This explains why conventional assays, in which IGFBPs are not removed, result in incorrect IGF-I values, which reflect more the present amount of IGFBP rather than the exact concentration of IGF-I (Figure 1.) (15,21). The low IGF-I concentrations require often additional efforts after the extraction procedure to concentrate the extract for obtaining a satisfactory sensitivity. The IGFBP-blocked IGF-I ELISA avoids these problems and allows the simple, correct and sensitive IGF-I determination in numerous samples at the same time.

### **INTENDED USE**

This ELISA kit is suitable for the scientific and diagnostic measurement of IGF-I in human serum or plasma, in cerebrospinal fluid and other human body fluids or conditioned media of human cell lines. Due to the high cross-reactivity with IGF-I from other mammalian species, it can also be used as a **assay** for determination of IGF-I in **primates, cattle, pig, sheep, horse, donkey, goat, dog, cat, rabbit and guinea pig**, however for rat, mouse and chicken derived samples the kit is not suited.

## PRINCIPLE

The Demeditec ELISA for IGF-I DEE020 is a so-called Sandwich-Assay. It utilizes two specific and high affinity antibodies for this protein. The IGF-I in the sample binds to the immobilized first antibody on the microtiter plate, the biotinylated and Streptavidin-Peroxidase conjugated second specific anti-IGF-I-Antibody binds in turn to the immobilized IGF-I. In the closing

substrate reaction the turn of the colour will be high specific catalysed, quantitatively depending on the IGF-I-level of the samples.

In order to dissociate IGF-I from the IGFBPs, the samples must be diluted in an acidic buffer (**Sample Buffer PP**) (Figure 2). The diluted samples are then pipetted into the wells, by this the pH-value will be neutralized. After neutralization of the samples, the excess IGF-II occupies the IGF-binding sites of the binding proteins, thus allowing the measurement of resulting free IGF-I. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized. Due to the extremely low cross-reactivity of the IGF-I antibody with IGF-II, the excess of IGF-II does not disturb the interaction with IGF-I. The test runs like a conventional ELISA using a Streptavidin-Peroxidase-Enzyme Conjugate.

Figure 2: Principle of the IGFBP blocked IGF-I ELISA

## PERFORMANCE CHARACTERISTICS AND VALIDATION

The standards of the ELISA DEE020 are **recombinant IGF-I** in concentrations of **2,5,15,30 and 50 ng/ml**, respectively the assay range covers –at recommended normal sample dilution- the **range from 42 to 1050 ng/ml**. By varying the sample dilution this can be adapted to the special individual requirements.

### Sensitivity

The **analytical sensitivity** of the ELISA DEE020 yields **0.09 ng/ml** (2 SD of zero standard in 17fold determination).

The **Inter-** and **Intra-Assay variation** coefficients are less than 6.8% and 6.7 % respectively. Exemplary determinations are represented in the Table 1 and 2.

**Table 1: Inter-Assay-Variation**

	Mean value (ng/ml)	Standard deviation	VC %
Sample 1	174	11.79	6.79
Sample 2	494	11.11	2.25
Sample 3	142	8.68	6.11

**Table 2: Intra-Assay-Variation (n=18)**

	Mean value (ng/ml)	Standard deviation	VC %
Sample 1	144.8	9.63	6.65
Sample 2	140.79	7.15	5.08
Sample 3	138.02	7.86	5.69

**Linearity**

Dilution:	Sample 1 (recalculated, ng/ml)	Dilution:	Sample 2 (recalculated, ng/ml)
1:10	137.2	1:10	439.1
1:20	133.5	1:20	500.2
1:40	133.6	1:40	499.2
1:80	134.6	1:80	490.5
1:160	134.4	1:160	494.5
1:320	135.7	1:320	526.4
		1:640	463.69
<b>AV / SD / VC%</b>	134.8 / 1.4 / 1.04	<b>AV / SD / VC%</b>	487.6 / 28.2 / 5.79

AV = Average Value , SD = Standard Deviation; VC = Coefficient of Variation

**Recovery**

In different human sera the **recovery** was on average **98.67%** of the hypothetical expected amount.

The Demeditec IGF-I ELISA is calibrated against the WHO International Reference Standard preparation of IGF-I, **WHO NIBSC 02/254** (25-26).

**Validation**

The Demeditec IGF-I ELISA was developed in accordance with the Demeditec IGF-I RIA. The clinical validation of the radioimmunoassay was achieved by determining the IGF-I levels in a large number of normal children and adults, normal short statured children without GH deficiency, girls with Ullrich-Turner syndrome, children with Silver-Russell syndrome, patients with GH deficiency, children with familial tall stature, Sotos syndrome, patients with acromegaly, and children with precocious puberty

**SPECIMEN COLLECTION, PREPARATION AND STORAGE**

Serum samples as well as Heparin-, EDTA- and Citrat-Plasma samples are suited. Possible dilution of the sample by the anticoagulant must be considered.

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen -20°C or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although IGF-I levels were found to be unaffected by few cycles in our experiments.

The high sensitivity of the assays allows IGF-I determinations in small sample volumes. For most of the determinations (Serum- or Plasma-Samples and no extreme values expected) the **dilutions of 1:10 to 1:50 in Sample Buffer PP** should be suited. Generally, a dilution of **1:21 is well suited** for serum or plasma samples.

Suggestion for dilution protocol:

Please pipette 200 µl **Sample Buffer PP** in PE/PP-Tubes (application of a multi-stepper is recommended in larger series); subsequently add 10 µl serum or plasma samples (dilution 1:21). After mixing use 20 µl of this dilution per determination within max. 2 h.

If necessary – depending on the expected IGF-I level- it is possible to use higher or lower dilution in **Sample Buffer PP**. Attention: Serum or Plasma samples must be diluted at least 1:10 in Sample Buffer PP in order to achieve the sufficient sample acidification.

IGF-I concentration in other body fluids or cell culture supernatants however could differ strongly.

**REAGENTS PROVIDED**

1)	<b>SORB MTP</b>	<b>Microtiter plate</b> , ready for use: Microtiter plate with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with anti-human IGF-I Antibody, packed in a laminate bag.
2)	<b>CAL</b>	<b>Standards A-E</b> , lyophilised, contain recombinant human IGF-I. Standard values are <b>between 2 50 ng/ml</b> (2, 5, 15, 30 and 50 ng/ml) IGF-I and have to be reconstituted in <b>500 µl (each) in Sample Buffer PP</b> . After using store the reconstituted standards in the original flasks as soon as possible at –20°C. When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. 20 µl per well are used in the assay
3)	<b>SAM BUF</b>	<b>Sample Buffer PP</b> 25 ml, ready for use, please use for reconstitution of Standards and Controls and for dilution of samples and Controls
4)	<b>Control</b>	<b>Control Serum KS1, KS2</b> , 500 µl, lyophilised, contain human serum and has to be reconstituted in <b>500 µl Sample Buffer PP</b> . The reconstituted Control Sera must be stored in the original flask as soon as possible at –20°C after using. When using anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. The IGF-I target value concentrations and the respective ranges are given on the QC-Certificate. The dilution of the Control Sera should be according to the dilution of the respective samples.
5)	<b>Ab CONJ</b>	<b>Antibody Conjugate AK</b> , 9 ml, ready for use, contains the biotinylated anti-IGF-I antibody. use 80 µl for each well in the assay. <b>ATTENTION: READY FOR USE!</b>
6)	<b>ENZ CONJ</b>	<b>Enzyme Conjugate EK</b> , 12 ml, ready for use, contains horseradisch-peroxidase conjugate to streptavidin, use 100 µl for each well in the assay. <b>ATTENTION: READY FOR USE!</b>
7)	<b>WASHBUF 20x</b>	<b>Washing Buffer (WP)</b> , 50 ml, <b>20X concentrated</b> solution. <b>Washing Buffer (WP)</b> has to be diluted 1:20 with distilled or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Attention: After dilution the Washing Buffer is only 4 weeks stable, dilute only according to requirements.
8)	<b>SUB TMB</b>	<b>Substrate (S)</b> , 12 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H <sub>2</sub> O <sub>2</sub> Tetramethylbencidine.
9)	<b>STOP SOLN</b>	<b>Stopping Solution (SL)</b> , 12 ml, ready for use, 0.4 N acididc solution, Caution acid!
10)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

## MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes and multichannel pipettes with disposable plastic tips  
Distilled or deionized water for dilution of the Washing Buffer (WP)  
Vortex-mixer  
Microtiter plate shaker (350 rpm)  
Microtiter plate washer (recommended)  
Micro plate reader ("ELISA-Reader") with filter for 450 and  $\geq 590$  nm  
Polyethylene PE/Polypropylene PP tubes for dilution of samples

## TECHNICAL NOTES

### Room temperature incubation means: Incubation at 20 - 25°C.

The microtiter plate and all reagents are stable unopened until the expiry date, if stored in the dark at 2° - 8°C (see label).

The Standards **A – E** and **Control Sera** are reconstituted with the **Sample Buffer PP** provided in the Kit. It is recommended to keep the reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The shelf life of the components after opening is not affected, if used appropriately. Store the unused seal stripes of the microtiter plate together with the desiccant at 2-8°C. Reconstituted Components (**Standards A – E** and **Control Sera**) should be stored at -20°C (or below). Freezing extends the expiry at least 3 months. When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay.

The 1:20 diluted Washing Buffer **WP** is only 4 weeks stable. Please dilute only according to daily requirements.

Before use, all kit components should be brought to room temperature. **Precipitates, possible in buffers, should be dissolved before use through mixing and warming.**

The **Substrate Solution S**, stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbenzidine, is photosensitive – store and incubate in the dark.

When performing the assay, the Standards **A-E**, Control Sera **KS** and the samples should be pipetted as fast as possible (e.g., 15 minutes). To avoid distortions due to differences in incubation times the Enzyme Conjugate **EK** as well as the succeeding **Substrate Solution S** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution SL** should be added to the plate in the same order as the Substrate Solution **S**.

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtitre plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must become adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/or false values, excessive shaking may result in high optical densities and/or false values.

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided washing buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µl at least.

The danger of handling with potentially infectious material must be taken into account.

When using an automatic microtitre plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

---

Demeditec Diagnostics GmbH • Lise-Meitner-Straße 2 • D-24145 Kiel (Germany)

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtitre plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

## STORAGE CONDITIONS

The microtiter plate wells and all undiluted reagents are stable until the expiry date, if stored in the dark at 2-8°C.

Store the unused seal strips and microtiter wells together with the desiccant at 2° to 8°C.

The Substrate Solution (S), stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbenzidine, is photosensitive – store and incubate in the dark.

Store the unused seal stripes of the microtiter plate together with the desiccant at 2-8°C. Reconstituted components (**Standards A – E** and **Control Sera**) should be stored at -20°C. Freezing extends the expiry at least 2 months. When using the Standards or Control Serum KS anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay.

## WARNINGS AND PRECAUTIONS

### For In Vitro Diagnostic Use.

#### For professional use only.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to **room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming. **Temperature WILL affect the absorbance** readings of the assay. However, values for the patient samples will not be affected.

Reagents from Kits with different lot numbers should not be mixed. Do not use expired reagents.

The microplate contains snap-off strips. Unused wells must be stored at 2 - 8°C in the sealed foil pouch and used in the frame provided.

Caution: This kit contains material of human and/or animal origin. Source of human serum for the Control Serum provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

### Stop Solution contains 0.2 M Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>)

R36/38 Irritating to eyes and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

Following components contain < 0.01% 2-Methyl-4-isothiazolin-3-one solution as preservative **AK, EK, PP**

R34 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

---

Demeditec Diagnostics GmbH • Lise-Meitner-Straße 2 • D-24145 Kiel (Germany)

Phone: +49 (0)431/71922-0 • Fax. +49 (0)431/71922-55

Email: info@demeditec.de • http://www.demeditec.com

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice  
S36/37 Wear suitable protective clothing and gloves  
S45 In case of accident or if you feel unwell seek medical advice

Following components contain < 0.01% (w/w) **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** as preservative: **AK, EK, PP, WP**

R36/38 Irritating to eyes and skin  
R43 Sensibilisation through skin contact possible  
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice  
S28.1 S28.1 After contact with skin, wash immediately with plenty of water

**TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.**

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed  
R36/37/38 Irritating to eyes, respiratory system and skin  
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice  
S28.1 After contact with skin, wash immediately with plenty of water  
S36/37 Wear suitable protective clothing and gloves

**General first aid procedures:**

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

## ASSAY PROCEDURE

NOTES: All determinations (Standards, Control Sera and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

- 1) Add **80 µl Antibody Conjugate AK** in all wells used
- 2) Pipette in positions A1/2 **20 µl Sample Buffer PP**
- 3) Pipette in positions B1/2 **20 µl of the Standard A (2 ng/ml)**,  
pipette in positions C1/2 **20 µl of the Standard B (5 ng/ml)**,  
pipette in positions D1/2 **20 µl of the Standard C (15 ng/ml)**,  
pipette in positions E1/2 **20 µl of the Standard D (30 ng/ml)**,  
pipette in positions F1/2 **20 µl of the Standard E (50 ng/ml)**.  
To control the correct accomplishment of the assay **20 µl** of the 1:21 (or in respective dilution ratio of the samples) in Sample Buffer diluted **Control Sera KS1/KS2** can be pipetted in positions G1/2 and H1/2.  
Pipette **20 µl** each of the diluted samples (e.g. dilute 1:21 with Sample Buffer **PP**)  
In the rest of wells, according to your requirements.
- 4) Cover the wells with sealing tape and incubate the plate for **1 hour** at **room temperature** (shake at 350 rpm)
- 5) After incubation aspirate the contents of the wells and wash the wells **5 times 300 µl Washing Buffer WP** / well.
- 6) Following the last washing step pipette **100 µl** of the **Enzyme Conjugate EK** in each well.
- 7) Cover the wells with sealing tape and incubate the plate for **30 Minutes** at **room temperature** (shake at 350 rpm).
- 8) After incubation wash the wells **5 times** with Washing Buffer **WP** as described in step 5.
- 9) Pipette **100 µl** of the **Substrate Solution S** in each well.
- 10) Incubate the microtiter plate for **15 minutes in the dark** at **room temperature**.
- 11) Stop the reaction by adding **100 µl Stopping Solution SL** to all wells.
- 12) Measure the absorbance within **30 minutes at 450 nm (Reference filter ≥ 590 nm)**.

## CALCULATION OF RESULTS

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of standard E should be above 1.00. Samples, which yield higher absorbance values than **Standard E**, are beyond the standard curve, for reliable determinations such samples should be retested at a higher dilution.

## Establishing the Standard Curve

The standards provided contain the following concentration of recombinant hIGF-I:

Standard	A	B	C	D	E
ng/ml	2	5	15	30	50
nmol/L	0.26	0.66	1.96	3.92	6.54

The Conversion factor of ng/mL in nmol/L is 0.13074

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other values.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The **concentration in ng/ml** of the samples can be calculated by **multiplication** with the respective **dilution factor**.

## EXPECTED NORMAL VALUES

IGF-I levels are highly age-dependent in children, less so in adults until the age of about 60. The normal ranges in various age groups, which are log-normally distributed, are given in Table 5 by percentiles. Between 8 and 19 years of age, values are given for boys and girls separately, because the pubertal peak usually occurs approximately 2 years earlier in girls. A graphic presentation is shown in Figures 3, 4 and 5. A major problem for the interpretation of IGF-I values arises from the fact that short stature is often due to developmental delay rather than any metabolic or endocrine disorder (constitutional delay of growth and adolescence). The sharp rise in IGF-I levels during puberty may therefore cause some uncertainty as to whether or not it would be appropriate to relate measured values to chronological age. It is recommended to take the pubertal stage into account (Table 4 and Figure 6) to get a more complete picture of this situation.

## LIMITATIONS OF PROCEDURE

IGF-I levels depend to a great degree on GH secretion. Diminished IGF-I values, however, do not prove GH deficiency, because a number of other factors can influence the plasma concentration of IGF-I and must therefore be taken into account in order to make a correct interpretation. IGF-I levels decrease during fasting (more than 1 day), as a result of malnutrition, malabsorption, cachexia, impaired hepatic function, or in hypothyroidism and untreated diabetes mellitus. They may also be low in chronic inflammatory disease and malignancies. IGF-I levels are high in states of accelerated sexual development. In clinical situations with hyperprolactinemia or in patients with craniopharyngioma, normal levels may be observed despite GH deficiency. In late pregnancy, IGF-I levels are moderately elevated.

## Appendix /Anhang

**Table 4:** Normal range of serum IGF-I levels (ng/ml) at different pubertal stages according to Tanner. Because no significant difference between boys and girls is observed, both sexes are combined. Only children and adolescents between 7 and 17 years of age are included.

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**Table 5:** Serum levels of IGF-I in healthy subjects at various ages. Individuals between 8 and 19 years of age were classified according to gender, as the pubertal peak occurs almost 2 years earlier in girls than in boys

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Serum concentrations are given in ng/mL. Determined with IGFBP-blocked IGF-I RIA without extraction step (Blum and Breier 1994),27.


Fig. 3: Age-dependent normal range of serum IGF-I levels in girls

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
Fig. 5: Age-dependent normal range of serum IGF-I levels in adults

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
Fig 6.: Serum IGF-I levels in normal children and adolescents (7 to 17 years) according to pubertal stages. Both sexes were included.

## LITREATURE / LITERATUR

- 1) Baxter RC. 1986 The somatomedins: insulin-like growth factors. *Adv Clin Chem.*25:49-115
- 2) Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
- 3) Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors.* New York: Elsevier.
- 4) Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology.* 112:2215-2217.
- 5) Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem.* 253:2769-2276.
- 6) Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab.* 13:113-143.
- 7) Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh).* 121:751-752.
- 8) Drop SLS. 1992 Report on the nomenclature of the IGF binding proteins. *J. Clin Endocrinol Metab.* 74:1215-1216.
- 9) Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem.* 261:8754-8760.
- 10) Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J.* 8:2497-2502.
- 11) Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest.* 60:648-657.
- 12) Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med.* 109:355-363.
- 13) Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol.* 128:347-357.
- 14) Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab.* 51:781-788.
- 15) Ranke MB (ed): *Diagnostics of Endocrine Function in Children and Adolescents,* Basel, Karger, 2003, pp 166-199
- 16) Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr.* 109:428-433.
- 17) Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N.Engl J Med.* 301:1138-1142
- 18) Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest.* 68:1321-1330.
- 19) Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinologic Diagnostics in Children and Adolescence.* Mannheim: J + J Verlag; 102-117.
- 20) Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab.* 55:147-153.
- 21) Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh).* 118:374-380.
- 22) Wilson EM, Oh Y, Rosenfeld RG (1997) Generation and characterization of an IGFBP-7 antibody: Identification of 31 kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media. *J Clin Endocrinol Metab* Vol 82, 4:1301-1303
- 23) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Horm Res* 54:60-68
- 24) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001) Relevance of IGF-I, IGFBP-3 and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Horm Res* 55:155-124
- 25) Address NIBSC:Blance Lane, South Mimms, Potters Bar, Hertford EN 6 3 QG, Great Britain

- 26) Burns C, Rigsby P, Moore M, Rafferty B. (2009) The first International Standard for Insulin-like Growth Factor-1 (IGF-1) for immunoassay: Preparation and calibration in an international collaborative study. *Growth Hormone & IGF Research* (2009)
- 27) Blum WH, Breier BH (1994) Radioimmunoassays for IGFs and IGF-BPs. *Growth Regulation* 4 (Suppl. 1):11-19
- 28) Ranke MB, Osterziel KJ, Schweitzer R, Schuett B., Weber K., Röbbel P, Vornwald A, Blumenstock G, Elmlinger MB (2003) Reference Levels of Insulin-Like-Growth Factor I in the Serum of Healthy Adults: Comparison of Four Immunoassays. *Clin Chem Lab Med* Vol 41(10):1329-1334

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## SUMMARY – DEMEDITEC IGF-I ELISA DEE020

Reconstitution / Dilution of Reagents		
Standards A-E	Reconstitution in <b>Sample Buffer PP</b>	<b>500 µl</b>
Control Serum KS1	Reconstitution in <b>Sample Buffer PP</b>	<b>500 µl</b>
Control Serum KS2	Reconstitution in <b>Sample Buffer PP</b>	<b>500 µl</b>
Wash Buffer WP	dilute in <b>A. dest.</b> (eg. total volume of 50 ml in a graduated flask and fill up to 1000 ml)	<b>1:20</b>
<b>Sample + Control Sera KS1 and KS2: dilute 1:21 in Sample Buffer PP, mix immediately, incubate max. 2h. Use 20 µl for each well in the assay.</b>		
Before conducting the assay equilibrate all <b>reagents to room temperature.</b>		

### Assay Procedure for Double Determinations:

Pipette	Reagent	Position
80 µl	Antibody Conjugate <b>AK</b>	in <b>all</b> wells used
20 µl	Sample Buffer <b>PP</b> (blank)	A1 and A2
20 µl	Standard <b>A (2 ng/ml)</b>	B1 and B2
20 µl	Standard <b>B (5 ng/ml)</b>	C1 and C2
20 µl	Standard <b>C (15 ng/ml)</b>	D1 and D2
20 µl	Standard <b>D (30 ng/ml)</b>	E1 and E2
20 µl	Standard <b>E (50 ng/ml)</b>	F1 and F2
20 µl	Control Serum <b>KS1</b>	G1 and G2
20 µl	Control Serum <b>KS2</b>	H1 and H2
20 µl	Diluted Samples	following wells
Cover the wells with the sealing tape.		

### Incubation: 1 h at RT, 350 rpm

5x 300 µl	Aspirate the contents of the wells and wash <b>5x</b> with <b>300 µl Wash Buffer WP</b>	each well
100 µl	Enzyme Conjugate <b>EK</b>	each well

### Incubation: 30 min at RT, 350 rpm

5x 300 µl	Aspirate the contents of the wells and wash <b>5x</b> with <b>300 µl Wash Buffer WP</b>	each well
100 µl	Substrate <b>S</b>	each well

### Incubation: 15 min in the dark RT

100 µl	Stop Solution <b>SL</b>	each well
Measure the absorbance within <b>30 min</b> at <b>450 nm</b> with <b>≥ 590 nm</b> as reference wavelength.		

<b>CAL</b> A-E	A -E	<b>Rec in</b> 500 µl <b>BUF</b> PP	
<b>Control</b>	KS1&KS2	<b>Rec in</b> 500 µl <b>BUF</b> PP	
<b>WASHBUF</b> 20x	WP		1:20 <b>DILU</b> A. dest.

**SPE** + **Control** 1:21 **SAM** **BUF** PP ↔ ⌚ max. 2 h

°C 20-25 °C

80 µl	<b>Ab CONJ</b>	A1 - End
20 µl	<b>SAM</b> <b>BUF</b> PP	A1/2
20 µl	<b>CAL</b> <b>A</b> Std A (2 ng/ml)	B1/2
20 µl	<b>CAL</b> <b>B</b> Std B (5 ng/ml)	C1/2
20 µl	<b>CAL</b> <b>C</b> Std C (15 ng/ml)	D1/2
20 µl	<b>CAL</b> <b>D</b> Std D (30 ng/ml)	E1/2
20 µl	<b>CAL</b> <b>E</b> Std E (50 ng/ml)	F1/2
20 µl	<b>CONTROL</b> KS 1	G1/2
20 µl	<b>CONTROL</b> KS 2	H1/2
20 µl	<b>SPE</b>	
<b>TAPE</b>		

⌚ 1 h °C 20-25 ↔ 350 rpm

5x 300 µl	5x <b>WASHBUF</b> WP
100 µl	<b>ENZ CONJ</b>
<b>TAPE</b>	

⌚ 0.5 h °C 20-25 ↔ 350 rpm

5 x 300 µl	5x <b>WASHBUF</b> WP
100 µl	<b>SUB</b> <b>TMB</b> S

⌚ 15 min °C 20-25

100 µl	<b>STOP SOLN</b> SL
<b>MEASURE</b>	