# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

User's Manual

# Adiponectin mouse ELISA

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## PACKAGE INSERT ENGLISH

### Demeditec Maus Adiponectin ELISA DEE091M

- is suited for Adiponectin determination in **Mouse Serum** and **Plasma** samples,
- is **fast**: incubation time a total of 2 hours and 30 minutes
- Single Standards with **0.025, 0.075, 0.15, 0.3, 0.65 and 1 ng/ml** native Mouse-Adiponectin are provided in the Kit
- 2 Control Sera are provided for quality control
- uses **high affinity monoclonal antibodies** against mouse Adiponectin
- Microtiter plates are separately breakapart

### Intended Use

Measurement of Adiponectin in Mouse serum and plasma samples.

### Introduction

Adiponectin was described for the first time in the early 90th of the last century as an endocrine factor produced by adipocytes. Adiponectin is involved in regulation of energy- and fat metabolism. So its concentration in the circulation is said to reflect the risk of atherosclerosis and the degree of insulin resistance. Based on the high incidence of these diseases, adiponectin was and still is object of intensive research regarding the underlying biological mechanisms and regarding its value as biomarker. Beside different cell culture models and studies with human patients, mice and rats are suitable model organisms for basic research and pre-clinical studies.

Therefore we developed and validated this testsystem as a tool for adiponectin measurements in mice usable in research and pre-clinical studies.

Even if the comparability of mice and humans is limited we offer some background information on *human* adiponectin physiology in the following section:

Adiponectin is a 30kDa protein which percentage in serum proteins is 0.01%. It is mainly synthesized by Adipocytes, but also muscle cells and hepatocytes have the ability to synthesize Adiponectin. Until now, IGF-I is the only known natural inductor of the synthesis. It consists of a Collagen-like N-terminal and a globular C-terminal domain [1]. In vivo Adiponectin appears with different oligomers. Beside the trimer and dimer also high molecular multimers exist [1-3]. Up to now two different receptors are known, both receptors are ubiquitarily expressed, though the distribution in the tissues varies. The Adiponectin Receptor 1 (AdipoR1) is especially in muscle- and AdipoR2 in liver tissue synthesized [4].

The significance for the organism is not clear until now. First studies show, that adiponectin correlates negatively with BMI and thus it could have relevance for the energy metabolism for example through the regulation of fatty acid oxidation. Beside the correlation with BMI, Adiponectin level is associated with the Insulin-Resistance [5-7] and so also linked with Type II Diabetes. Adiponectin is associated also with glucose- und lipometabolism [8, 9].

The formerly proposed diagnostic value of the high molecular weight form of adiponectin was not verified using a commercially available testsystem for the determination of HMW adiponectin [10]. Blueher et al. clearly demonstrate that regarding the diagnosis of insulin resistance, measured by whole body glucose uptake below 40  $\mu\text{mol}/\text{kg}\cdot\text{min}$ , total adiponectin as determined with the Demeditec DEE09, is with an area of 0.92 under the receiver-operating curve, of greater diagnostic value [10].

Furthermore it is involved in inflammatory processes [11-15] and therewith it is of importance for appearance of arteriosclerosis [4, 5, 16] and coronary disease [17, 18], thus the determination of Adiponectin level in plasma could serve to estimate the risk of coronary disease [19, 20]. Beside this Adiponectin influences further physiological processes as for example the angiogenesis [21, 22].

## Reagents Provided

1)	MTP	<b>Microtiter plate</b> , ready for use: Microtiter plate with 96 wells, dived up in 12 stripes à 8 wells (separately breakapart), coated with anti-Maus Adiponectin antibody.
2)	CAL	<b>Standards A-F</b> , lyophilised, contain native Mouse Adiponectin. Standard values are between <b>0.025 - 1 ng/ml</b> (0.025, 0.075, 0.15, 0.3, 0.65 und 1 ng/ml) Adiponectin and have to be reconstituted in <b>1 ml (each) Dilution Buffer VP</b> . 100 µl per well are used in the assay. If the standards are required for more than one assay process we recommend to store the reconstituted Standards frozen at -20°C. <b>Attention:</b> Standards should be thawed only once – where required please store aliquoted in adequate volumes.
3)	DILU	<b>Dilution buffer VP, 125 ml</b> , ready for use, please use for the reconstitution of Standards A-F, Control Sera KS1 & KS2 and for the serum dilution.
4)	Control	<b>Control Serum KS1 &amp; KS2</b> , each <b>100 µl</b> lyophilised: Contains mouse Serum and has to be reconstituted in <b>100 µl Dilution Buffer VP</b> . The Adiponectin target value concentration and the respective range is given the vial label. The dilution of the Control Sera KS 1&2 should be according to the dilution of the respective samples, the target value concentration should be obtained by <b>multiplication</b> with the respective <b>dilution factor</b> .
5)	Ab CONJ	<b>Antibody-POD-Conjugate AK, 12 ml, ready for use</b> , contains a mixture of biotinylated anti-Adiponectin antibody and HRP (Horseradish Peroxidase)-labelled Streptavidin. Use 100 µl/well in the assay.
6)	WASHBUF 20x	<b>Washing Buffer (WP), 50 ml, 20 X concentrated</b> solution. <b>Washing Buffer (WP)</b> has to be diluted 1:20 with distilled or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Attention: After dilution the Washing Buffer is only 4 weeks stable, dilute only according to requirements.
7)	SUBST	<b>Substrate (S), 12 ml</b> , ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H <sub>2</sub> O <sub>2</sub> Tetramethylbencidine.
8)	H <sub>2</sub> SO <sub>4</sub>	<b>Stopping Solution (SL), 12 ml</b> , ready for use, 0.2 M sulphuric acid, Caution acid!
9)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

## Materials Required but not provided

- Precision pipettes (100 and 200µl) Micropipettes and multichannel pipettes with disposable plastic tips
- Distilled or Deionized water for dilution of the Washing Buffer (WP)
- Vortex-mixer
- Device to aspirate the standards and the samples from the wells
- Timer (120 min. range)
- Reservoirs (disposable)
- Plate washer and plate shaker (recommended)
- Calibrated Micro plate reader ("ELISA-Reader") with filter for 450 and 620nm (or ≥590 nm)
- Foil welding device for laminate bags (recommended)

## WARNINGS AND PRECAUTIONS

### For in-vitro diagnostic use only. For professional use only.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought **to room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming.

Do not mix reagents of different lots. Do not use expired reagents.

The microplate contains snap-off strips. Unused wells must be stored at 2 - 8°C in the sealed foil pouch and used in the frame provided.

Caution: This kit contains material of human and/or animal origin.

No known test methods can offer total assurance of the absence of infectious agents; therefore all components and specimens should be treated as potentially infectious.

Following components contain < 0.01% 2-Methyl-4-isothiazolin-3-one solution as preservative : **A-F, AK, VP**

R34 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S36/37 Wear suitable protective clothing and gloves

S45 In case of accident or if you feel unwell seek medical advice

Following components contain < 0.01%(w/w) 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-isothiazol-3-one as preservative: **A-F, AK, VP, WP**

R36/38 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice S28.1 After contact with skin, wash immediately with plenty of water

### Stop solution contains 0.2 M Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>)

R36/38 Irritating to eyes and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

**TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.** Store and Incubate in the dark.

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed  
R36/37/38 Irritating to eyes, respiratory system and skin  
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice  
S28.1 After contact with skin, wash immediately with plenty of water  
S36/37 Wear suitable protective clothing and gloves

#### **General first aid procedures:**

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

#### **Method**

The Demeditec ELISA for Adiponectin DEE091M is a so-called Sandwich-Assay using two specific and high affinity antibodies. The Adiponectin in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-Adiponectin-Antibody binds in turn to the immobilised Adiponectin. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction the turn of the colour will be catalysed quantitatively depending on the Adiponectin-level of the samples.

#### **Specimen**

Serum and plasma samples of mice can be used in this assay.

Influence of Heparin (30IE/mL), EDTA (6,8mM) and NaCitrat (0,015M) on the measurement of Adiponectin has been investigated by recovery experiments. PBS was enriched with recombinant mouse Adiponectin and the above mentioned substances. No significant influence on the recovery of adiponectin was detected.

Hemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation.

#### **Storage of the samples**

Storage at RT max. 2 days

Storage at -20 °C max. 2 years

More than five freeze/thaw cycles are not possible.

#### **Sample Preparation**

Samples have to be diluted in Dilution Buffer (VP).

A sample dilution of 1:10 000 is in general suitable. However, the Adiponectin levels can vary individually significant, we would therefore recommend to check this and adjust the dilution respectively.

## Technical Recommendations

The assay has to be conducted strictly according the test protocol herein.

Reagents with different lot numbers cannot be mixed. The microtiterplate and reagents are stable until the indicated expiry if stored unopened and protected from sunlight at 2 – 8°C.

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming.

Incubation at room temperature means: 20-25°C

## Standards and Controls

For the reconstitution of the lyophilised components (Standards A - F and Control Sera KS1 &KS2) the kit Dilution Buffer VP has to be used. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam!) with a Vortex mixer.

The reconstituted standard and controls can be stored for 2 month at –20°C. Repeated freeze/thaw cycles have to be avoided.

## Washing Buffer

The required volume of washing buffer is prepared by 1:20 dilution of the provided 20fold concentrate with deionised water. The diluted Washing Buffer is stable for 4 weeks at 2-8°C. It has to be at room temperature for usage!

## Microtiterplate

Unused microtiterplate stripes have to be stored airtight together with the desiccant bag at 2-8°C. The labelled expiry is not influenced in case of proper storage.

**Assay Procedure**

All determinations (Standards, Control Sera KS1 & KS2 and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Standards, Control Sera and the samples should be pipette as fast as possible (e.g., <15 minutes).

All incubations have to be conducted at room temperature (20-25°C)

To avoid distortions due to differences in incubation times, Antibody-POD-Conjugate AK as well as the following Substrate Solution S should be added to the plate in the same order and in the same time interval as the samples. Stop Solution SL should be added to the plate in the same order as the Substrate Solution.

- 1) Add **100 µl** Dilution Buffer VP in the first wells. Subsequently add 100 µl Standard or 100 µl of diluted Control Sera or diluted samples.
- 2) Cover the wells with sealing tape and incubate the plate for **1 hour at room temperature** (if possible, shake at  $\geq 350$  rpm )
- 3) After incubation aspirate the contents of the wells and wash the wells 3 times **250 µl Washing Buffer WP** / well. The washing buffer should incubate for at least for 15 seconds/cycle.
- 4) Following the last washing step pipette **100 µl** of the **Antibody-POD-Conjugate AK** in each well.
- 5) Cover the wells with sealing tape and incubate the plate for **1 hour at room temperature** (if possible shake  $\geq 350$  rpm).
- 6) After incubation wash the wells 3 times with Washing Buffer as described in step 3
- 7) Pipette **100 µl of the TMB Substrate** Solution in each well.
- 8) Incubate the plate for **30 minutes in the dark at room temperature (20 - 25°C)**.
- 9) Stop the reaction by adding **100 µl of Stopping Solution**.
- 10) Measure the colour reaction within 30 minutes at 450nm (reference filter  $\geq 590$  nm).

## Calculation of Results

### Establishing the Standard Curve

For the evaluation of the assay it is preconditioned that the absorbance values of the blank should be below 0.25, these of standard F should exceed 1.0.

Samples, which yield higher absorbance values than Standard F are beyond the standard curve, for reliable determinations these samples should be tested anew with a higher dilution.

Standards are provided in the following concentrations (use the concentration unit as preferred):

Standard	A	B	C	D	E	F
ng/ml	0.025	0.075	0.15	0.3	0.65	1

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other values.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program. **Non-linear regression**, a higher-grade polynomial or four parametric logistic (4-PL) curve fit are in general suited for the evaluation.
- 6) The Adiponectin concentration of the diluted sample or the diluted control sera KS1&2 in ng/ml (or µg/ml according the chosen unit for the standards) is calculated in this way, the Adiponectin concentration of the **undiluted sample** and of KS1 & KS2 is calculated **by multiplication** with the respective dilution factor.

## Performance Characteristics

### Standards

The Standards of the ELISA DEE09 are prepared from **native Adiponectin** in concentrations of **0.025, 0.075, 0.15, 0.3, 0.65 and 1 ng/ml**. The native Adiponectin was calibrated against the recombinant Protein (Manufacturer: R&D Systems, Wiesbaden).

### Sensitivity

The analytical sensitivity of the ELISA DEE091M yields 0.008 ng/mL (equal to < 0.0008 ng per well; as 2xSD of zero standard in 22fold determination).

### Interassay Variability

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6	Sample 7
Mean	1.12	6.12	7.65	8.04	2.15	7.54	6.87
SD	0.09	1.19	0.60	0.56	0.25	0.55	0.42
VC%	7.66	3.79	7.80	6.98	11.80	7.28	6.08
n	8	8	9	9	9	9	6

### Intraassay Variability

	Sample 1	Sample 2
Mean [µg/mL]	7.395	1.712
SD	159	37
VC%	2.15	2.16
n	21	21

### Linearity

Linearity				
Dilution	Sample 1	Dilution	Sample 2	Sample 3
1:14000	11.673	1:10000	9.59	
1:28000	10.928	1:12000	11.229	11.886
1:56000	10.67	1:24000	9.856	10.822
1:112000	10.421	1:48000	9.915	10.222
		1:96000	9.684	9.999

recalculated Values in ng/mL

## Literatur / Literature

1. Nakano, Y., et al., *Isolation and characterization of GBP28, a novel gelatin-binding protein purified from human plasma*. J Biochem (Tokyo), 1996. **120**(4): p. 803-12.
2. Pajvani, U.B., et al., *Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity*. J Biol Chem, 2003. **278**(11): p. 9073-85.
3. Tsao, T.S., et al., *Role of disulfide bonds in Acrp30/adiponectin structure and signaling specificity. Different oligomers activate different signal transduction pathways*. J Biol Chem, 2003. **278**(50): p. 50810-7.
4. Shimada, K., T. Miyazaki, and H. Daida, *Adiponectin and atherosclerotic disease*. Clin Chim Acta, 2004. **344**(1-2): p. 1-12.
5. Higashiura, K., et al., *Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations*. Clin Endocrinol (Oxf), 2004. **61**(6): p. 753-9.
6. Spranger, J., et al., *Adiponectin is independently associated with insulin sensitivity in women with polycystic ovary syndrome*. Clin Endocrinol (Oxf), 2004. **61**(6): p. 738-46.
7. Zoico, E., et al., *Adipocytokines, fat distribution, and insulin resistance in elderly men and women*. J Gerontol A Biol Sci Med Sci, 2004. **59**(9): p. M935-9.
8. Ye, J.M., et al., *PPARalpha /gamma ragaglitazar eliminates fatty liver and enhances insulin action in fat-fed rats in the absence of hepatomegaly*. Am J Physiol Endocrinol Metab, 2003. **284**(3): p. E531-40.
9. Yamauchi, T., et al., *Adiponectin stimulates glucose utilization and fatty-acid oxidation by activating AMP-activated protein kinase*. Nat Med, 2002. **8**(11): p. 1288-95.
10. Blüher (Blueher), M., et al., *Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays*. Diabetes Care, 2007. **30**(2): p. 280-5.
11. Delaigle, A.M., et al., *Induction of adiponectin in skeletal muscle by inflammatory cytokines: in vivo and in vitro studies*. Endocrinology, 2004. **145**(12): p. 5589-97.
12. Winzer, C., et al., *Plasma adiponectin, insulin sensitivity, and subclinical inflammation in women with prior gestational diabetes mellitus*. Diabetes Care, 2004. **27**(7): p. 1721-7.
13. Xydakis, A.M., et al., *Adiponectin, inflammation, and the expression of the metabolic syndrome in obese individuals: the impact of rapid weight loss through caloric restriction*. J Clin Endocrinol Metab, 2004. **89**(6): p. 2697-703.
14. Motoshima, H., et al., *Adiponectin suppresses proliferation and superoxide generation and enhances eNOS activity in endothelial cells treated with oxidized LDL*. Biochem Biophys Res Commun, 2004. **315**(2): p. 264-71.
15. Wolf, A.M., et al., *Adiponectin induces the anti-inflammatory cytokines IL-10 and IL-1RA in human leukocytes*. Biochem Biophys Res Commun, 2004. **323**(2): p. 630-5.
16. Okamoto, Y., et al., *Adiponectin reduces atherosclerosis in apolipoprotein E-deficient mice*. Circulation, 2002. **106**(22): p. 2767-70.
17. Schlegel, A., *Adiponectin and risk of coronary heart disease*. Jama, 2004. **292**(1): p. 40; author reply 40.
18. Choi, K.M., et al., *Inflammation, Insulin Resistance, and Glucose Intolerance in Acute Myocardial Infarction Patients without a Previous Diagnosis of Diabetes Mellitus*. J Clin Endocrinol Metab, 2004.
19. Nakamura, Y., et al., *Implications of plasma concentrations of adiponectin in patients with coronary artery disease*. Heart, 2004. **90**(5): p. 528-33.
20. Pischon, T., et al., *Plasma adiponectin levels and risk of myocardial infarction in men*. Jama, 2004. **291**(14): p. 1730-7.
21. Shibata, R., et al., *Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling*. J Biol Chem, 2004. **279**(27): p. 28670-4.
22. Fernandez-Real, J.M., et al., *Adiponectin is associated with vascular function independent of insulin sensitivity*. Diabetes Care, 2004. **27**(3): p. 739-45.

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**SUMMARY – DEMEDITEC Mouse Adiponectin ELISA DEE091M**

Reagent preparation:	Reconstitution:	Dilution
<b>Standards A-F</b>	in 1 ml Dilution Buffer <b>VP</b>	
<b>Control Sera KS1 &amp; KS2</b>	in 100 µl Dilution Buffer <b>VP</b>	<b>1:10 000 with</b> Dilution Buffer <b>VP</b>
<b>Washing Buffer WP</b>		<b>1:20 with Aqua. dest.</b> (e.g., add the complete contents of the flask <b>(50 ml)</b> into a graduated flask and fill with A.dest. to 1000 ml).
<b>Sample Dilution:</b> e.g. 1 :10 000		

**Assay Procedure for Double Determination:**

Pipette	Reagents	Position
100 µl	Dilution Buffer <b>VP</b>	A1/2
100 µl	Standard <b>A (0.025 ng/ml)</b>	B1/2
100 µl	Standard <b>B (0.075 ng/ml)</b>	C1/2
100 µl	Standard <b>C (0.15 ng/ml)</b>	D1/2
100 µl	Standard <b>D (0.3 ng/ml)</b>	E1/2
100 µl	Standard <b>E (0.65 ng/ml)</b>	F1/2
100 µl	Standard <b>F (1 ng/ml)</b>	G1/2
100 µl	Control Serum <b>KS1</b>	H1/2
100 µl	Control Serum <b>KS2</b>	A3/4
100 µl	Sample dilution	following wells
Cover the wells with the sealing tape.		
<b>Incubation: 1 h at RT, ≥ 350 rpm</b>		
3x 250 µl	Aspirate the contents of the wells and wash 3x with <b>250 µl Wash Buffer WP</b>	each well
100 µl	Antibody-POD-Conjugate <b>AK</b>	each well
<b>Incubation: 1 h at RT, ≥350 rpm</b>		
3x 250 µl	Aspirate the contents of the wells and wash 3x with <b>250 µl Wash Buffer WP.</b>	each well
100 µl	Substrate Solution <b>S</b>	each well
<b>Incubation: 30 min in the Dark at RT</b>		
100 µl	Stopping Solution <b>SL</b>	each well
Measure the absorbance within <b>30 min</b> at <b>450 nm</b> with <b>≥590 nm</b> as reference wavelength.		

REF DEE091M

## International Test description

<b>CAL</b> A-F	A -F	<b>Rec in</b> 1 ml VP	
<b>Control</b>	KS1 & KS2	<b>Rec in</b> 100 µl VP	1:10 000 <b>DILU</b> VP
<b>WASHBUF</b> 20x	WP		1:20 <b>DILU</b> A. dest.

<b>SPE</b>		1:10 000 <b>DILU</b> VP
<b>°C</b> 20-25 °C		
100 µl	<b>VP</b>	A1/2
100 µl	<b>CAL</b> A (0.025 ng/ml)	B1/2
100 µl	<b>CAL</b> B (0.075 ng/ml)	C1/2
100 µl	<b>CAL</b> C (0.15 ng/ml)	D1/2
100 µl	<b>CAL</b> D (0.3 ng/ml)	E1/2
100 µl	<b>CAL</b> E (0.65 ng/ml)	F1/2
100 µl	<b>CAL</b> F (1 ng/ml)	G1/2
100 µl	<b>CONTROL</b> KS1 1:10 000 <b>DILU</b> VP	H1/2
100 µl	<b>CONTROL</b> KS2 1:10 000 <b>DILU</b> VP	A3/4
100 µl	<b>SPE</b> 1:10 000 <b>DILU</b> VP	
<b>TAPE</b>		

1 h **°C** 20-25 350 rpm

3x 250 µl	3x <b>WASHBUF</b> WP
100 µl	<b>AbCONJ</b> AK
<b>TAPE</b>	

1 h **°C** 20-25 ≥ 350 rpm

3x 250 µl	3x <b>WASHBUF</b> WP
100 µl	<b>SUBST</b> <b>TMB</b> S

0.5 h **°C** 20-25

100 µl	<b>H<sub>2</sub>SO<sub>4</sub></b> SL
<b>MEASURE</b>	