Instructions for use
Adrenaline Plasma ELISA

REF

DEE4100

IVD

1. Introduction

1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of adrenaline (epinephrine) in plasma.

Adrenaline (epinephrine) is extracted from a plasma sample*) by using a cis-diol-specific affinity gel, acylated and then modified enzymatically.

The competitive ELISA kit uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The derivatized standards, controls and samples and the solid phase bound analytes compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a reference curve prepared with known standard concentrations.

*) Flexible sample volumes between 100 – 600 µl can be used with this assay.

1.2 Clinical application

In humans the catecholamines adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine are neurotransmitter of the sympathetic nervous system and are involved in many physiological processes. The sympathetic nervous system sets the body to a heightened state of alert, also called as the body's fight-or-flight response.

In the human body the catecholamines and their metabolites indicate the adaption of the body to acute and chronic stress.

Next to the metanephrine/normetanephrine the catecholamines are important for the diagnosis and the follow-up of tumors of the sympathicoadrenal system like the pheochromocytomas. The quantitative determination of catecholamines in urine is preferred for the diagnosis of these tumors, whereas the determination of catecholamines in plasma is medically sensible for the localization of the tumor and for function testing. Values above the cut-off can provide an indication for neuroendocrine tumors.

However, in literature various diseases like hypertension, cardiovascular diseases, schizophrenia and manic depression are described with different levels of catecholamines

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

2. Procedural cautions, guidelines, warnings and limitations

2.1 Procedural cautions, guidelines and warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for certain types of samples as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) Reagents of this kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- (4) The principles of Good Laboratory Practice (GLP) have to be followed.
- (5) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (6) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (7) For dilution or reconstitution purposes, use deionized, distilled or ultra-pure water.
- (8) The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided. Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.

- (9) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (10) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (11) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (12) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (13) A standard curve must be established for each run.
- (14) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (15) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (16) Avoid contact with Stop Solution containing 0.25 M H₂SO₄. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (17) Some reagents contain sodium azide (NaN₃) as preservatives. In case of contact with eyes or skin, rinse off immediately with water. NaN₃ may react with lead and copper plumbing to form explosive metal azides. When disposing reagents, flush with a large volume of water to avoid azide build-up.
- (18) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (19) For information on hazardous substances included in the kit please refer to Material Safety Data Sheet (MSDS). The Material Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (20) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (21) The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence but have to be correlated to other diagnostic tests and clinical observations.
- (22) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.

2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

2.2.1 Interfering substances

Plasma

Samples containing precipitates or fibrin strands or which are haemolytic or lipemic might cause inaccurate results.

2.2.2 Drug interferences

There are no known substances (drugs) which ingestion interferes with the measurement of adrenaline level in the sample.

2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 - 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

4. Materials

4.1 Contents of the kit

BA D-0090  **Adhesive Foil** - Ready to use

Contents: Adhesive Foils in a resealable pouch

Volume: 1 x 4 foils

BA D-0032  **Microtiter Plate** - Ready to use

Contents: 2 x 96 wells, empty in a resealable pouch

- BA E-0030** **WASH-CONC 50x** **Wash Buffer Concentrate** - Concentrated 50x
 Contents: Buffer with a non-ionic detergent and physiological pH
 Volume: 1 x 20 ml/vial, light purple cap
- BA E-0040** **CONJUGATE** **Enzyme Conjugate** - Ready to use
 Contents: Goat anti-rabbit immunoglobulins, conjugated with peroxidase
 Volume: 1 x 12 ml/vial, red cap
- BA E-0055** **SUBSTRATE** **Substrate** - Ready to use
 Contents: Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide
 Volume: 1 x 12 ml/black vial, black cap
- BA E-0080** **STOP-SOLN** **Stop Solution** - Ready to use
 Contents: 0.25 M sulfuric acid
 Volume: 1 x 12 ml/vial, light grey cap
- BA E-0131** **ADR MN** **Adrenaline Microtiter Strips** - Ready to use
 Contents: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable blue foil pouch with desiccant, blue coloured
- BA E-4110** **ADR-AS** **Adrenaline Antiserum** - Ready to use
 Contents: Rabbit anti-adrenaline antibody, blue coloured
 Volume: 1 x 6 ml/vial, blue cap

Standards and Controls - Ready to use

Cat. no.	Component	Colour/ Cap	Concentration	Concentration	Volume/ Vial	
			pg/ml	pmol/l		
			ADR	ADR		
BA R-4601	STANDARD A	white	0	0	4 ml	
BA R-4602	STANDARD B	light yellow	20	109	4 ml	
BA R-4603	STANDARD C	orange	60	328	4 ml	
BA R-4604	STANDARD D	dark blue	200	1 092	4 ml	
BA R-4605	STANDARD E	light grey	800	4 368	4 ml	
BA R-4606	STANDARD F	black	3 200	17 472	4 ml	
BA R-4651	CONTROL 1	light green	Refer to QC-Report for expected value and acceptable range!		4 ml	
BA R-4652	CONTROL 2	dark red			4 ml	

Conversion: Adrenaline (pg/ml) x 5.46 = Adrenaline (pmol/l)

Contents: Acidic buffer with non-mercury stabilizer, spiked with defined quantity of adrenaline

BA R-0050 **ADJUST-BUFF** **Adjustment Buffer** - Ready to use

Contents: TRIS buffer

Volume: 1 x 4 ml/vial, green cap

BA R-4617 **TE-BUFF** **TE Buffer** - Ready to use

Contents: TRIS-EDTA buffer

Volume: 1 x 4 ml/vial, brown cap

BA R-6611 **ACYL-BUFF** **Acylation Buffer** - Ready to use

Contents: Buffer with light alkaline pH for the acylation

Volume: 1 x 20 ml/vial, white cap

BA R-6612 ACYL-REAG **Acylation Reagent** - Ready to use

Contents: Acylation reagent in DMF and DMSO

Volume: 1 x 3 ml/vial, light red cap

Hazards

identification:

H225 Highly flammable liquid and vapour.
H360 May damage fertility or the unborn child.
H319 Causes serious eye irritation.

BA R-6614 COENZYME **Coenzyme** - Ready to use

Contents: S-adenosyl-L-methionine

Volume: 1 x 4 ml/vial, purple cap

BA R-6615 ENZYME **Enzyme** - Lyophilized

Contents: Catechol-O-methyltransferase

Volume: 4 vials, pink cap

BA R-6618 EXTRACT-PLATE 48 **Extraction Plate** - Ready to use

Contents: 2 x 48 well plates coated with boronate affinity gel in a resealable pouch

BA R-6619 HCL **Hydrochloric Acid** - Ready to use

Contents: 0.025 M Hydrochloric Acid, yellow coloured

Volume: 1 x 20 ml/vial, dark green cap

4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 25 – 700 µl; 1 ml
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 – 650 nm
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Temperature controlled incubator (37 °C) or similar heating device
- Absorbent material (paper towel)
- Water (deionized, distilled or ultra-pure)
- Vortex mixer

5. Sample collection and storage

Plasma

Whole blood should be collected by venipuncture into centrifuge tubes containing EDTA as anti-coagulant (Plasma Monovette™ or Vacuette™) and centrifuged at room temperature immediately after collection.

Haemolytic and especially lipemic samples should not be used for the assay.

Storage: up to 6 hours at 2 - 8 °C, for longer period (up to 6 month) at -20 °C.

Repeated freezing and thawing should be avoided.

6. Test procedure

A plasma volume between **100 µl-600 µl** is needed per single determination.

If a plasma volume **< 600 µl** is used, water (deionized, distilled or ultra-pure) has to be added to a **final volume of 600 µl** and this **prediluted sample** has to be used for the extraction procedure (please refer to point 6.2 of this protocol).

This sample predilution has to be considered in the calculation of results (please refer to point 7 of this protocol).

Allow all reagents to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent, and the extinction values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. Corresponding variations also apply to the incubation times. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

 In case of overflow, read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 405 nm

6.1 Preparation of reagents

Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled or ultra-pure) to a final volume of 1000 ml.

Storage: 1 month at 2 – 8 °C

Enzyme Solution

Reconstitute the content of the vial labelled 'Enzyme' with 1 ml water (deionized, distilled or ultra-pure) and mix thoroughly. Add 0.3 ml of Coenzyme followed by 0.7 ml of Adjustment Buffer. The total volume of the Enzyme Solution is 2.0 ml.

 The Enzyme Solution has to be prepared freshly prior to the assay (not longer than 10 - 15 minutes in advance). Discard after use!

6.2 Sample preparation, extraction and acylation

1.	Pipette 30 µl of standards, controls and 600 µl of plasma samples into the respective wells of the Extraction Plate .
2.	Add 500 µl of water (deionized, distilled or ultra-pure) to the wells with standards and controls .
3.	Pipette 25 µl of TE Buffer into all wells.
4.	Cover the plate with Adhesive Foil . Incubate 60 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
5.	Remove the foil and empty the plate. Blot dry by tapping the inverted plate on absorbent material.
6.	Pipette 1 ml of Wash Buffer into all wells.
7.	Incubate 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
8.	Blot dry by tapping the inverted plate on absorbent material.
9.	Wash one more time as described (step 6, 7 and 8)!
10.	Pipette 150 µl of Acylation Buffer into all wells.
11.	Pipette 25 µl of Acylation Reagent into all wells.
12.	Incubate 20 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
13.	Empty the plate and blot dry by tapping the inverted plate on absorbent material.
14.	Pipette 1 ml of Wash Buffer into all wells.
15.	Incubate 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
16.	Blot dry by tapping the inverted plate on absorbent material.
17.	Wash one more time as described (step 14, 15, 16).
18.	Pipette 200 µl of Hydrochloric Acid into all wells.
19.	Cover plate with Adhesive Foil . Incubate 10 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
	Do not decant the supernatant thereafter!
	190 µl of the supernatant is needed for the subsequent enzymatic conversion

6.3 Enzymatic conversion

1.	Pipette 190 µl of the extracted standards, controls and samples into the respective wells of the Microtiter Plate .		
2.	Add 50 µl of Enzyme Solution (refer to 6.1) to all wells.		
3.	Cover plate with Adhesive Foil . Incubate 1 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).		
4.	Incubate for 2 h at 37 °C . The following volumes of the supernatants are needed for the subsequent ELISA: <table border="1" data-bbox="199 2004 571 2038"><tr><td>Adrenaline</td><td>100 µl</td></tr></table>	Adrenaline	100 µl
Adrenaline	100 µl		

6.4 Adrenaline ELISA

1.	Pipette 100 µl of standards, controls and samples from the Enzyme Plate (refer to 6.3) into the respective pre-coated Adrenaline Microtiter Strips .
2.	Pipette 50 µl of the respective Adrenaline Antiserum into all wells.
3.	Cover the plate with Adhesive Foil . Incubate 1 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
4.	Incubate for 15 – 20 h (overnight) at 2 – 8 °C .
5.	Remove the foil. Discard or aspirate the contents of the wells. Wash the plate 4 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
6.	Pipette 100 µl of Enzyme Conjugate into all wells.
7.	Cover the plate with Adhesive Foil . Incubate 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
8.	Remove the foil. Discard or aspirate the contents of the wells. Wash the plate 4 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
9.	Pipette 100 µl of Substrate into all wells.
10.	Incubate 20 - 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). ⚠ Avoid exposure to direct sunlight!
11.	Pipette 100 µl of Stop Solution into all wells.
12.	Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended).

7. Calculation of results

Measuring range	Adrenaline
	10 – 3 200 pg/ml

The calibration curves are obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

The concentrations of the **undiluted plasma samples** and the **controls** can be read directly from the standard curve.

Concentration of diluted plasma samples:

If only a plasma volume **< 600 µl** was used for the extraction, the concentration read from the standard curve has to be **multiplied** with a **volume-factor**:

$$\text{Volume-factor} = \frac{600 \mu\text{l}}{\text{used plasma volume } (\mu\text{l})}$$

Conversion

Adrenaline (pg/ml) x 5.46 = Adrenaline (pmol/l)

Expected reference value

It is strongly recommended that each laboratory should determine its own reference value.

In a study conducted with apparently normal healthy adults, using the Adrenaline Plasma ELISA ^{High Sensitive} the following value is observed:


Expected Reference Value	Adrenaline
	< 100 pg/ml

7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. The kit controls or other commercial controls should fall within established confidence limits. The confidence limits of the kit controls are indicated on the QC-Report.

7.2 Typical standard curve

 Example, do not use for calculation!

8. Assay characteristics

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
	Derivatized Adrenaline	100
	Derivatized Noradrenaline	0.20
	Derivatized Dopamine	< 0.01
	Metanephrine	0.64
	Normetanephrine	< 0.01
	3-Methoxytyramine	< 0.01
	3-Methoxy-4-hydroxyphenylglycol	0.03
	Tyramine	< 0.01
	Phenylalanine, Caffeinic acid, L-Dopa, Homovanillic acid, Tyrosine, 3-Methoxy-4-hydroxymandelic acid	< 0.01

Analytical Sensitivity (600 µl undiluted sample) Limit of Detection	Adrenaline
	5.2 pg/ml

Precision Intra-Assay CV				
	Sample	Mean (pg/ml)	SD (pg/ml)	CV (%)
Adrenaline	low	54.0	9.1	16.8
	medium	113	13.2	11.7
	high	447	48.2	10.8

Recovery	Mean (%)	Range (%)	Measurement Range (pg/ml)
Adrenaline	110	93 - 121	28.9 - 2669


Linearity		Serial dilution up to	Mean (%)	Range (%)
	Adrenaline		1:512	90

9. References/Literature

- (1) Kim et al. Vitamin C prevents stress-induced damage on the heart caused by the death of cardiomyocytes, through the down-regulation of the excessive production of catecholamine, TNF- α , and ROS production in GULO(-I-) Vit C-Insufficient mice. Free Radical Biology and Medicine, 65:573-583 (2013)
- (2) Bada et al. Peripheral vasodilatation determines cardiac output in exercising humans: insight from atrial pacing. The Journal of Physiology, 590(8):2051-2060 (2012)
- (3) Parks et al. Employment and work schedule are related to telomere length in women. Occupational & Environmental Medicine 68(8):582-589 (2011)

 **For updated literature or any other information please contact your local supplier.**

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	LOT	Batch code	IVD	For in-vitro diagnostic use only!
	Consult instructions for use	CONT	Content	CE	CE labelled
	Caution	REF	Catalogue number	RUO	For research use only!

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