---

**Instructions for use**  
**2-CAT Plasma ELISA**

**REF**

**DEE4500**

  
2 x 96

  
+2 / +8 °C

**IVD**

**CE**

## **1. Introduction**

### **1.1 Intended use and principle of the test**

Enzyme Immunoassay for the quantitative determination of adrenaline (epinephrine) and noradrenaline (norepinephrine) in plasma.

Adrenaline (epinephrine) and noradrenaline (norepinephrine) are extracted from a plasma sample\*) by using a cis-diol-specific affinity gel, acylated and then modified enzymatically.

The competitive ELISA kit uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The derivatized standards, controls and samples and the solid phase bound analytes compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a reference curve prepared with known standard concentrations.

\*) Flexible sample volumes between 100 – 600 µl can be used with this assay.

### **1.2 Clinical application**

In humans the catecholamines adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine are neurotransmitter of the sympathetic nervous system and are involved in many physiological processes. The sympathetic nervous system sets the body to a heightened state of alert, also called as the body's fight-or-flight response.

In the human body the catecholamines and their metabolites indicate the adaption of the body to acute and chronic stress.

Next to the metanephrine/normetanephrine the catecholamines are important for the diagnosis and the follow-up of tumors of the sympathicoadrenal system like the pheochromocytomas. The quantitative determination of catecholamines in urine is preferred for the diagnosis of these tumors, whereas the determination of catecholamines in plasma is medically sensible for the localization of the tumor and for function testing. Values above the cut-off can provide an indication for neuroendocrine tumors.

However, in literature various diseases like hypertension, cardiovascular diseases, schizophrenia and manic depression are described with different levels of catecholamines

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

## **2. Procedural cautions, guidelines, warnings and limitations**

### **2.1 Procedural cautions, guidelines and warnings**

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for certain types of samples as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) Reagents of this kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- (4) The principles of Good Laboratory Practice (GLP) have to be followed.
- (5) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (6) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (7) For dilution or reconstitution purposes, use deionized, distilled or ultra-pure water.
- (8) The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided.

- (9) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (10) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (11) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (12) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (13) A standard curve must be established for each run.
- (14) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (15) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (16) Avoid contact with Stop Solution containing 0.25 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (17) Some reagents contain sodium azide (NaN<sub>3</sub>) as preservatives. In case of contact with eyes or skin, rinse off immediately with water. NaN<sub>3</sub> may react with lead and copper plumbing to form explosive metal azides. When disposing reagents, flush with a large volume of water to avoid azide build-up.
- (18) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (19) For information on hazardous substances included in the kit please refer to Material Safety Data Sheet (MSDS). The Material Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (20) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (21) The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence but have to be correlated to other diagnostic tests and clinical observations.
- (22) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.

## 2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

### 2.2.1 Interfering substances

#### Plasma

Samples containing precipitates or fibrin strands or which are haemolytic or lipemic might cause inaccurate results.

### 2.2.2 Drug interferences

There are no known substances (drugs) which ingestion interferes with the measurement of catecholamine level in the sample.

### 2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

## 3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 - 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

## 4. Materials

### 4.1 Contents of the kit

<b>BA D-0090</b>	<b>FOILS</b>	<b>Adhesive Foil</b> - Ready to use
Contents:	Adhesive Foils in a resealable pouch	
Volume:	2 x 4 foils	
<b>BA E-0030</b>	<b>WASH-CONC 50x</b>	<b>Wash Buffer Concentrate</b> - Concentrated 50x
Contents:	Buffer with a non-ionic detergent and physiological pH	
Volume:	2 x 20 ml/vial, light purple cap	
<b>BA R-6618</b>	<b>EXTRACT-PLATE 48</b>	<b>Extraction Plate</b> - Ready to use
Contents:	2 x 48 well plates coated with boronate affinity gel in a resealable pouch	

- BA E-0040** **CONJUGATE** **Enzyme Conjugate** - Ready to use  
 Contents: Goat anti-rabbit immunoglobulins, conjugated with peroxidase  
 Volume: 2 x 12 ml/vial, red cap
- BA E-0055** **SUBSTRATE** **Substrate** - Ready to use  
 Contents: Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide  
 Volume: 2 x 12 ml/black vial, black cap
- BA E-0080** **STOP-SOLN** **Stop Solution** - Ready to use  
 Contents: 0.25 M sulfuric acid  
 Volume: 2 x 12 ml/vial, light grey cap
- BA E-0131** **ADR MN** **Adrenaline Microtiter Strips** - Ready to use  
 Contents: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable blue foil pouch with desiccant, blue coloured
- BA E-0231** **NAD MNM** **Noradrenaline Microtiter Strips** - Ready to use  
 Contents: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable yellow foil pouch with desiccant, yellow coloured
- BA E-4110** **ADR-AS** **Adrenaline Antiserum** - Ready to use  
 Contents: Rabbit anti-adrenaline antibody, blue coloured  
 Volume: 1 x 6 ml/vial, blue cap
- BA E-4210** **NAD-AS** **Noradrenaline Antiserum** - Ready to use  
 Contents: Rabbit anti-noradrenaline antibody, yellow coloured  
 Volume: 1 x 6 ml/vial, yellow cap

**Standards and Controls** - Ready to use

Cat. no.	Component	Colour/ Cap	Concentration pg/ml		Concentration pmol/l		Volume/ Vial
			ADR	NAD	ADR	NAD	
<b>BA R-4601</b>	STANDARD A	white	0	0	0	0	4 ml
<b>BA R-4602</b>	STANDARD B	light yellow	20	80	109	473	4 ml
<b>BA R-4603</b>	STANDARD C	orange	60	240	328	1 418	4 ml
<b>BA R-4604</b>	STANDARD D	dark blue	200	800	1 092	4 728	4 ml
<b>BA R-4605</b>	STANDARD E	light grey	800	3 200	4 368	18 912	4 ml
<b>BA R-4606</b>	STANDARD F	black	3 200	12 800	17 472	75 648	4 ml
<b>BA R-4651</b>	CONTROL 1	light green	Refer to QC-Report for expected value and acceptable range!				4 ml
<b>BA R-4652</b>	CONTROL 2	dark red					4 ml

Conversion: Adrenaline (pg/ml) x 5.46 = Adrenaline (pmol/l)  
 Noradrenaline (pg/ml) x 5.91 = Noradrenaline (pmol/l)

Contents: Acidic buffer with non-mercury stabilizer, spiked with defined quantity of adrenaline and noradrenaline

- BA R-0050** **ADJUST-BUFF** **Adjustment Buffer** - Ready to use  
 Contents: TRIS buffer  
 Volume: 1 x 4 ml/vial, green cap

- BA D-0032** **96** **Microtiter Plate** - Ready to use  
 Contents: 1 x 96 wells, empty in a resealable pouch

- BA R-4617** TE-BUFF **TE Buffer** - Ready to use  
 Contents: TRIS-EDTA buffer  
 Volume: 1 x 4 ml/vial, brown cap
- BA R-6611** ACYL-BUFF **Acylation Buffer** - Ready to use  
 Contents: Buffer with light alkaline pH for the acylation  
 Volume: 1 x 20 ml/vial, white cap
- BA R-6612** ACYL-REAG **Acylation Reagent** - Ready to use  
 Contents: Acylation reagent in DMF and DMSO  
 Volume: 1 x 3 ml/vial, light red cap  
 Hazards identification:   
 H225 Highly flammable liquid and vapour.  
 H360 May damage fertility or the unborn child.  
 H319 Causes serious eye irritation.
- BA R-6614** COENZYME **Coenzyme** - Ready to use  
 Contents: S-adenosyl-L-methionine  
 Volume: 1 x 4 ml/vial, purple cap
- BA R-6615** ENZYME **Enzyme** - Lyophilized  
 Contents: Catechol-O-methyltransferase  
 Volume: 4 vials, pink cap
- BA R-6619** HCL **Hydrochloric Acid** - Ready to use  
 Contents: 0.025 M Hydrochloric Acid, yellow coloured  
 Volume: 1 x 20 ml/vial, dark green cap

#### 4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 25 – 700 µl; 1 ml
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 – 650 nm
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Temperature controlled incubator (37 °C) or similar heating device
- Absorbent material (paper towel)
- Water (deionized, distilled or ultra-pure)
- Vortex mixer

#### 5. Sample collection and storage

##### Plasma

Whole blood should be collected by venipuncture into centrifuge tubes containing EDTA as anti-coagulant (Monovette™ or Vacuette™ for plasma) and centrifuged at room temperature immediately after collection.

Haemolytic and especially lipemic samples should not be used for the assay.

Storage: up to 6 hours at 2 - 8 °C, for longer period (up to 6 month) at -20 °C.

Repeated freezing and thawing should be avoided.

#### 6. Test procedure

A plasma volume between **100 µl-600 µl** is needed per single determination.

If a plasma volume **< 600 µl** is used, water (deionized, distilled or ultra-pure) has to be added to a **final volume of 600 µl** and this **prediluted sample** has to be used for the extraction procedure (please refer to point 6.2 of this protocol).

This sample predilution has to be considered in the calculation of results (please refer to point 7 of this protocol).

Allow all reagents to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent, and the extinction values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. Corresponding variations also apply to the incubation times. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

*In case of overflow, read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 405 nm*

## 6.1 Preparation of reagents

### Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled or ultra-pure) to a final volume of 1000 ml.

Storage: 1 month at 2 – 8 °C

### Enzyme Solution

Reconstitute the content of the vial labelled 'Enzyme' with 1 ml water (deionized, distilled or ultra-pure) and mix thoroughly. Add 0.3 ml of Coenzyme followed by 0.7 ml of Adjustment Buffer. The total volume of the Enzyme Solution is 2.0 ml.

*The Enzyme Solution has to be prepared freshly prior to the assay (not longer than 10 - 15 minutes in advance). Discard after use!*

## 6.2 Sample preparation, extraction and acylation

1.	Pipette <b>30 µl</b> of <b>standards, controls</b> and <b>600 µl</b> of <b>plasma samples</b> into the respective wells of the <b>Extraction Plate</b> .
2.	Add <b>500 µl</b> of <b>water</b> (deionized, distilled or ultra-pure) to the wells with <b>standards</b> and <b>controls</b> .
3.	Pipette <b>25 µl</b> of <b>TE Buffer</b> into all wells.
4.	Cover the plate with <b>Adhesive Foil</b> . Incubate <b>60 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
5.	Remove the foil and empty the plate. Blot dry by tapping the inverted plate on absorbent material.
6.	Pipette <b>1 ml</b> of <b>Wash Buffer</b> into all wells.
7.	Incubate <b>5 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
8.	Blot dry by tapping the inverted plate on absorbent material.
9.	<b>Wash one more time</b> as described (step 6, 7 and 8)!
10.	Pipette <b>150 µl</b> of <b>Acylation Buffer</b> into all wells.
11.	Pipette <b>25 µl</b> of <b>Acylation Reagent</b> into all wells.
12.	Incubate <b>20 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
13.	Empty the plate and blot dry by tapping the inverted plate on absorbent material.
14.	Pipette <b>1 ml</b> of <b>Wash Buffer</b> into all wells.
15.	Incubate <b>5 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
16.	Blot dry by tapping the inverted plate on absorbent material.
17.	<b>Wash one more time</b> as described (step 14, 15, 16).
18.	Pipette <b>200 µl</b> of <b>Hydrochloric Acid</b> into all wells.
19.	Cover plate with <b>Adhesive Foil</b> . Incubate <b>10 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
	<b>Do not decant the supernatant thereafter!</b>
	<b>190 µl of the supernatant is needed for the subsequent enzymatic conversion</b>

### 6.3 Enzymatic conversion

1.	Pipette <b>190 µl</b> of the <b>extracted standards, controls</b> and <b>samples</b> into the respective wells of the <b>Microtiter Plate</b> .				
2.	Add <b>50 µl</b> of <b>Enzyme Solution</b> (refer to 6.1) to all wells.				
3.	Cover plate with <b>Adhesive Foil</b> . Incubate <b>1 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).				
4.	Incubate for <b>2 h</b> at <b>37°C</b> . The following volumes of the supernatants are needed for the subsequent ELISA: <table border="1" data-bbox="199 403 1029 448"><tr><td><b>Adrenaline</b></td><td><b>100 µl</b></td><td><b>Noradrenaline</b></td><td><b>10 µl</b></td></tr></table>	<b>Adrenaline</b>	<b>100 µl</b>	<b>Noradrenaline</b>	<b>10 µl</b>
<b>Adrenaline</b>	<b>100 µl</b>	<b>Noradrenaline</b>	<b>10 µl</b>		

### 6.4 Adrenaline ELISA

1.	Pipette <b>100 µl</b> of <b>standards, controls</b> and <b>samples</b> from the <b>Enzyme Plate</b> (refer to 6.3) into the respective pre-coated <b>Adrenaline Microtiter Strips</b> .
2.	Pipette <b>50 µl</b> of the respective <b>Adrenaline Antiserum</b> into all wells.
3.	Cover the plate with <b>Adhesive Foil</b> . Incubate <b>1 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
4.	Incubate for <b>15 – 20 h</b> (overnight) at <b>2 – 8 °C</b> .
5.	Remove the foil. Discard or aspirate the contents of the wells. Wash the plate <b>4 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
6.	Pipette <b>100 µl</b> of <b>Enzyme Conjugate</b> into all wells.
7.	Cover the plate with <b>Adhesive Foil</b> . Incubate <b>30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
8.	Remove the foil. Discard or aspirate the contents of the wells. Wash the plate <b>4 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
9.	Pipette <b>100 µl</b> of <b>Substrate</b> into all wells.
10.	Incubate <b>20 - 30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).  <b>Avoid exposure to direct sunlight!</b>
11.	Pipette <b>100 µl</b> of <b>Stop Solution</b> into all wells.
12.	<b>Read</b> the <b>absorbance</b> of the solution in the wells within 10 minutes, using a microplate reader set to <b>450 nm</b> (if available a reference wavelength between 620 nm and 650 nm is recommended).

### 6.5 Noradrenaline ELISA

1.	Pipette <b>10 µl</b> of <b>standards, controls</b> and <b>samples</b> from the <b>Enzyme Plate</b> (refer to 6.3) into the respective pre-coated <b>Noradrenaline Microtiter Strips</b> .
2.	Pipette <b>50 µl</b> of the respective <b>Noradrenaline Antiserum</b> into all wells.
3.	Cover the plate with <b>Adhesive Foil</b> . Incubate <b>1 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
4.	Incubate for <b>15 – 20 h</b> (overnight) at <b>2 – 8 °C</b> .
5.	Remove the foil. Discard or aspirate the contents of the wells. Wash the plate <b>4 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
6.	Pipette <b>100 µl</b> of <b>Enzyme Conjugate</b> into all wells.
7.	Cover the plate with <b>Adhesive Foil</b> . Incubate <b>30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).
8.	Remove the foil. Discard or aspirate the contents of the wells. Wash the plate <b>4 x</b> by adding <b>300 µl</b> of <b>Wash Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.
9.	Pipette <b>100 µl</b> of <b>Substrate</b> into all wells.
10.	Incubate <b>20 - 30 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm).  <b>Avoid exposure to direct sunlight!</b>
11.	Pipette <b>100 µl</b> of <b>Stop Solution</b> into all wells.
12.	<b>Read</b> the <b>absorbance</b> of the solution in the wells within 10 minutes, using a microplate reader set to <b>450 nm</b> (if available a reference wavelength between 620 nm and 650 nm is recommended).

## 7. Calculation of results

Measuring range	Adrenaline	Noradrenaline
	10 – 3 200 pg/ml	20 – 12 800 pg/ml

The calibration curves are obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

The concentrations of the **undiluted plasma samples** and the **controls** can be read directly from the standard curve.

### Concentration of diluted plasma samples:

If only a plasma volume **< 600 µl** was used for the extraction, the concentration read from the standard curve has to be **multiplied** with a **volume-factor**:

$$\text{Volume-factor} = \frac{600 \mu\text{l}}{\text{used plasma volume } (\mu\text{l})}$$

### Conversion

Adrenaline (pg/ml) x 5.46 = Adrenaline (pmol/l)

Noradrenaline (pg/ml) x 5.91 = Noradrenaline (pmol/l)

### Expected reference values

It is strongly recommended that each laboratory should determine its own reference values.

In a study conducted with apparently normal healthy adults, using the 2-CAT Plasma ELISA <sup>High Sensitive</sup> the following values are observed:

Expected Reference Values	Adrenaline	Noradrenaline
	< 100 pg/ml	< 600 pg/ml


### 7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. The kit controls or other commercial controls should fall within established confidence limits. The confidence limits of the kit controls are indicated on the QC-Report.

### 7.2 Typical standard curves

*Examples, do not use for calculation!*

## 8. Assay characteristics

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)	
		Adrenaline	Noradrenaline
	Derivatized Adrenaline	100	0.14
	Derivatized Noradrenaline	0.20	100
	Derivatized Dopamine	< 0.01	0.2
	Metanephrine	0.64	< 0.01
	Normetanephrine	< 0.01	0.48
	3-Methoxytyramine	< 0.01	< 0.01
	3-Methoxy-4-hydroxyphenylglycol	0.03	0.01
	Tyramine	< 0.01	< 0.01
	Phenylalanine, Caffeinic acid, L-Dopa, Homovanillic acid, Tyrosine, 3-Methoxy-4-hydroxymandelic acid	< 0.01	< 0.01

Analytical Sensitivity (600 µl undiluted sample) Limit of Detection	Adrenaline	Noradrenaline
	5.2 pg/ml	20 pg/ml

Precision				
Intra-Assay CV				
	Sample	Mean (pg/ml)	SD (pg/ml)	CV (%)
Adrenaline	low	54.0	9.1	16.8
	medium	113	13.2	11.7
	high	447	48.2	10.8
Noradrenaline	low	535	59.9	11.2
	medium	679	75.9	11.2
	high	2123	344	16.2

Recovery	Mean (%)	Range (%)	Measurement Range (pg/ml)
Adrenaline	110	93 - 121	38.9 - 2669
Noradrenaline	96	81 - 112	614 - 9606

Linearity		Serial dilution up to	Mean (%)	Range (%)
	Adrenaline	1:512	90	73 - 105
	Noradrenaline	1:512	96	76 - 107

## 9. References/Literature

- (1) Kim et al. Vitamin C prevents stress-induced damage on the heart caused by the death of cardiomyocytes, through the down-regulation of the excessive production of catecholamine, TNF- $\alpha$ , and ROS production in GULO(-I-) Vit C-Insufficient mice. *Free Radical Biology and Medicine*, 65:573-583 (2013)
- (2) Bada et al. Peripheral vasodilatation determines cardiac output in exercising humans: insight from atrial pacing. *The Journal of Physiology*, 590(8):2051-2060 (2012)
- (3) Parks et al. Employment and work schedule are related to telomere length in women. *Occupational & Environmental Medicine* 68(8):582-589 (2011)

 **For updated literature or any other information please contact your local supplier.**

**Symbols:**

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date	<b>LOT</b>	Batch code	<b>IVD</b>	For in-vitro diagnostic use only!
	Consult instructions for use	<b>CONT</b>	Content	<b>CE</b>	CE labelled
	Caution	<b>REF</b>	Catalogue number	<b>RUO</b>	For research use only!

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## 9. Referenzen/Literatur

- (1) Kim et al. Vitamin C prevents stress-induced damage on the heart caused by the death of cardiomyocytes, through the down-regulation of the excessive production of catecholamine, TNF- $\alpha$ , and ROS production in GULO(-/-) Vit C-Insufficient mice. Free Radical Biology and Medicine, 65:573-583 (2013)
- (2) Bada et al. Peripheral vasodilatation determines cardiac output in exercising humans: insight from atrial pacing. The Journal of Physiology, 590(8):2051-2060 (2012)
- (3) Parks et al. Employment and work schedule are related to telomere length in women. Occupational & Environmental Medicine 68(8):582-589 (2011)

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