Product information

Information about other products is available at: www.demeditec.com

User's Manual

Progesterone free in Saliva ELISA

IVD

REF

DES6633

96 Wells

CONTENTS / INHALTSVERZEICHNIS

1	INTRODUCTION	3
2	PRINCIPLE.....	3
3	WARNINGS AND PRECAUTIONS.....	3
4	REAGENTS.....	3
5	SPECIMEN COLLECTION AND PREPARATION.....	5
6	ASSAY PROCEDURE.....	6
7	EXPECTED NORMAL VALUES.....	7
8	QUALITY CONTROL.....	7
9	PERFORMANCE CHARACTERISTICS.....	8
10	LIMITATIONS OF PROCEDURE.....	10
11	LEGAL ASPECTS.....	10
12	REFERENCES.....	10
1	EINLEITUNG.....	13
2	TESTPRINZIP.....	13
3	VORSICHTSMAßNAHMEN.....	13
4	BESTANDTEILE DES KITS.....	13
5	PROBENVORBEREITUNG.....	15
6	TESTDURCHFÜHRUNG.....	16
7	ERWARTETE WERTE.....	17
8	QUALITÄTSKONTROLLE.....	18
9	ASSAY CHARACTERISTIKA.....	18
10	GRENZEN DES TESTS.....	18
11	RECHTLICHE GRUNDLAGEN.....	19
12	REFERENZEN.....	19
	SYMBOLS USED WITH DEMEDITEC ELISA.....	20

1 INTRODUCTION

1.1 Intended Use

Enzyme immunoassay for the *in vitro diagnostic* quantitative measurement of active free progesterone in saliva. Measurements obtained by this device may be used in the diagnosis and treatment of disorders of the ovaries or placenta and as an aid to confirm that ovulation takes place.

1.2 Summary and explanation

Progesterone (4-pregnene-3, 20-dione) is a C21 steroid hormone containing a keto-group (at C-3) and a double bond between C-4 and C-5. Like other steroids, it is synthesized from cholesterol via a series of enzyme-mediated steps (1). Progesterone is a female sex hormone of primary importance in ovulation, fertility and menopause. It is particularly important in preparing the endometrium for the implantation of the blastocyte and in maintaining pregnancy (2). In the follicular phase of menstrual cycle progesterone is produced in low levels. It increases to the LH peak and then sharply rises to high levels. Next there is a sharp decline to low levels of follicular phase. In non-pregnant women progesterone is mainly secreted by the corpus luteum whereas in pregnancy the placenta becomes the major source (3,4). Minor sources for progesterone are the adrenal cortex for both sexes and the testes for males.

The determination of progesterone in saliva combines a highly sensitive technique and non-invasive collection and represents the concentration of the metabolic active free progesterone.

2 PRINCIPLE

The **DEMEDIATEC Progesterone free in Saliva ELISA** Kit is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the principle of competitive binding. An unknown amount of free progesterone present in the sample and a defined amount of progesterone conjugated to horseradish peroxidase compete for the binding sites of rabbit polyclonal progesterone antiserum coated to the wells of a microplate. After one-hour incubation on a shaker the microplate is washed four times. After addition of the substrate solution the concentration of progesterone is inversely proportional to the optical density measured.

3 WARNINGS AND PRECAUTIONS

1. This kit is for *in vitro* diagnostic use only. For professional use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
4. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
5. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
6. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
7. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
8. Allow the reagents to reach room temperature (21-26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
9. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
10. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
11. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
12. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
13. Do not use reagents beyond expiry date as shown on the kit labels.

14. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiterplate readers.
15. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
16. Avoid contact with Stop Solution. It may cause skin irritation and burns.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
18. For information please refer to Material Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from Demeditec Diagnostics GmbH.

4 REAGENTS

4.1 Reagents provided

1. **SORB MT Microtiterplate**, 12 x 8 (break apart) strips with 96 wells;
Wells coated with an anti-Progesterone antibody (rabbit polyclonal antibody).
2. **CAL 0 Calibrator 0**, 1 vial, 3.0 ml, ready to use
3. **CAL 1-6 Calibrator (Calibrator 1-6)**, 6 vials, 1 ml each, ready to use;
Concentrations: 10 – 30 – 100 – 300 – 1000 – 5000 pg/ml
Conversion: Progesterone (pg/ml) x 3.18 = pmol/l
4. **CONTROL 1-2 Control low / Control high**, 2 vials, 1.0 ml each, ready to use;
For control values and ranges please refer to QC-Datasheet.
5. **ENZ CONJ Enzyme Conjugate**, 1 vial, 7 ml, ready to use;
Progesterone conjugated to horseradish peroxidase
6. **SUB TMB Substrate Solution**, 1 vial, 22 ml, ready to use;
contains tetramethylbenzidine (TMB).
7. **STOP SOLN Stop Solution**, 1 vial, 7 ml, ready to use;
contains 2 N Hydrochloric Acid solution.
8. **WASH SOLN 10x Wash Solution**, 1 vial, 50 ml (10X concentrated);
see „Preparation of Reagents“.

Note: Additional Calibrator 0 for sample dilution is available upon request.

4.2 Materials required but not provided

- A microtiter plate reader capable for endpoint measurement at 450nm
- Calibrated variable precision micropipettes (50 µl, 100 µl, 200 µl).
- Microplate mixer operating more than 600 rpm
- Microcentrifuge
- Vortex mixer
- Absorbent paper.
- Distilled or deionized water
- Timer
- Semi logarithmic graph paper or software for data reduction

4.3 Storage conditions

When stored at 2°C to 8°C unopened reagents will be stable until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2°C to 8°C. Microtiter wells must be stored at 2°C to 8°C. Once the foil bag has been opened, care should be taken to close it tightly again.

4.4 Reagent preparation

All reagents should be at room temperature before use.

Wash Solution:

Dilute 50 ml of 10X concentrated *Wash Solution* with 450 ml deionized water to a final volume of 500 ml.
The diluted Wash Solution is stable for at least 3 months at room temperature.

4.5 Disposal of the kits

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheet.

4.6 Damaged test kits

In case of any severe damage of the test kit or components, DEMEDITEC have to be informed written, latest one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Samples containing sodium azide should not be used in the assay. The saliva samples should be completely colorless. Even the slightest red color shows blood contamination resulting in falsely elevated concentration values. In case of visible blood contamination the patient should discard the sample, rinse the sampling device with water, wait for 10 minutes and take a new sample.

5.1 Specimen Collection

For the correct collection of saliva we are recommending to only use appropriate devices made from ultra-pure polypropylene. Do not use any PE devices or cotton based Salivettes for sampling. False readings will result. Glass tubes can be used as well, but in this case special attention is necessary for excluding any interference caused by the stopper. Please contact Demeditec Diagnostics for more details.

As food might contain significant amounts of steroid hormones samples preferably should be taken while fasting. If fasting should be a problem at least any food of animal origin (meat or dairy products) should be avoided prior to finalizing the collection. In the morning breakfast should be done only after finalizing the collection procedure. During the day the collection period should be timed just before an anticipated meal. As the steroid hormone secretion in saliva as well in serum shows an obvious dynamic secretion pattern throughout the day it is important to always collect 5 samples during a 2 hour period; this means every 30 minutes one sample. If possible the volume of each single sample should be a minimum of 0.5 ml (better 1 ml). Saliva flow may be stimulated by drinking water. This is allowed and even recommended before and during the collection period. Drinking of water is not allowed during the last 5 minutes before taking the single samples.

5.2 Specimen Storage and Preparation

Saliva samples in general are stable at ambient temperature for up to seven days. Therefore mailing of such samples by ordinary mail without cooling will not create any problem. Storage at 4°C can be done for a period of up to one month. Whenever possible samples should preferably be kept at a temperature of -20°C. Even repeated thawing and freezing is not a problem. Each sample has to be frozen, thawed, and centrifuged at least once anyhow in order to precipitate the mucins by centrifugation. Upon arrival of the samples at the lab the samples have to be kept frozen at least overnight. Next morning the samples are thawed and mixed carefully. The samples have to be centrifuged for 5 to 10 minutes. The clear colorless supernatant is easy to pipette. If the sample should show even a slightly red color it should be discarded. Otherwise the value most probably will be falsely elevated. Due to the episodic variations of the steroid secretion we highly recommend the strategy of multiple sampling. If such a set of multiple samples has to be tested the staff of lab (after at least one freezing, thawing, and centrifugation cycle) should mix aliquots of the 5 single samples and perform the determination using the mixture.

5.3 Specimen Dilution

Samples expected to contain progesterone concentrations higher than the highest calibrator (5000 pg/ml) should be diluted with the zero calibrator before assay. The additional dilution step has to be taken into account for the calculation of the result.

Example:

- a) Dilution 1:10: 10 µl saliva + 90 µl Calibrator 0 (mix thoroughly)
b) Dilution 1:100: 10 µl of dilution a) + 90 µl Calibrator 0 (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Respect the incubation times as stated in this instructions for use.

6.2 Assay procedure

Each run must include a standard curve.

1. Prepare a sufficient number of microplate wells to accommodate calibrators, controls and patient samples.
2. Dispense **100 µl** of each **calibrator, control and sample** with new disposable tips into appropriate wells.
3. Dispense **50 µl** of **Enzyme Conjugate** into each well.
4. Incubate for **60 minutes** at room temperature on a Microplate mixer.

Important Note:

Optimal reaction in this assay is markedly dependent on shaking of the microplate!

5. Discard the content of the wells and rinse the wells **4 times** with diluted **Wash Solution** (300 µl per well). Remove as much Wash Solution as possible by beating the microplate on absorbent paper.
6. Add **200 µl** of **Substrate Solution** to each well.
7. Incubate without shaking for **30 minutes** in the dark.
8. Stop the reaction by adding **50 µl** of **Stop Solution** to each well.
9. Determine the absorbance of each well at 450 nm. It is recommended to read the wells within 15 minutes.

6.3 Calculation of results

1. Calculate the average absorbance values for each set of calibrators, controls and patient samples.
2. Using semi logarithmic graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample, determine the corresponding concentration from the calibration curve.
4. Automated method: Computer programs using cubic spline, 4 PL (4 Parameter Logistics) or Logit-Log are recommended.
5. The concentration of the samples can be determined directly from this calibrator curve. Samples with concentrations higher than that of the highest calibrator have to be further diluted. For the calculation of the concentrations, this dilution factor has to be taken into account.

Conversion to SI units:

Progesterone (pg/ml) x 3.18 = pmol/l

6.3.1 Example of typical calibrator curve

Following data are intended for illustration only and should not be used to calculate results from another run.

Standard	Absorbance Units
Calibrator 0 (0 pg/ml)	2.326
Calibrator 1 (10 pg/ml)	2.103
Calibrator 2 (30 pg/ml)	1.912
Calibrator 3 (100 pg/ml)	1.548
Calibrator 4 (300 pg/ml)	1.257
Calibrator 5 (1000 pg/ml)	0.831
Calibrator 6 (5000 pg/ml)	0.455

7 EXPECTED NORMAL VALUES

In order to determine the normal range of salivary Progesterone, saliva samples from 101 adult male and 268 adult female apparently healthy subjects, ages 15 to 75 years, were collected in the morning and analyzed using the DEMEDITEC Progesterone free in Saliva ELISA kit. The following ranges are calculated with the results of this study.

	Age group	Salivary progesterone pg/ml (5 – 95% Percentile)
Women	15 - 55 yrs. Follicular phase n = 38	3.7 – 81.4 pg/ml
	15 - 55 yrs. Luteal phase n = 116	73.1 – 271.5 pg/ml
	> 55 yrs. Postmenopausal n = 114	5.5 – 98.9 pg/ml
Men Children (boys)	15 - 75 yrs. n = 101	8.6 – 107.0 pg/ml
	< 12 yrs. n = 7	14.5 – 43.4 pg/ml

Therapy should not be decided based on results alone. The results should be correlated to other clinical observations and diagnostic tests.

Furthermore, we recommend that each laboratory establish its own range for the population tested, because the values differ between age, new born, children, adolescents and adults.

8 QUALITY CONTROL

Good laboratory practice requires that controls should be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to national regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated at the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analyzing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DEMEDITEC directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Analytical Sensitivity

The lowest analytical detectable level of progesterone that can be distinguished from the Zero Calibrator is 5.0 pg/ml at the 2SD confidence limit.

9.2 Specificity (Cross Reactivity)

The following materials have been evaluated for cross reactivity. The percentage indicates cross reactivity at 50% displacement compared to Progesterone.

Steroid	% Cross reaction
Androstenedione	4.0
Androsterone	< 0.1
Cortisol	< 0.1
Cortison	< 0.1
Corticosterone	2.1
Danazol	< 0.1
11-Desoxycorticosterone	4.4
11-Desoxycortisol	< 0.1
Dexamethasone	< 0.1
5 α -Dihydrotestosteron	< 0.1
Estradiol	< 0.1
Estriol	< 0.1
Estrone	< 0.1
17 α -Hydroxyprogesterone	1.6
Prednisolon	< 0.1
Prednison	< 0.1
Pregnenolone	< 0.1
Testosterone	< 0.1

9.3 Assay dynamic range

The range of the assay is between 10 – 5000 pg/ml.

9.4 Reproducibility

9.4.1 Intra-Assay

The intra-assay variation was determined by 20 replicate measurements of 3 saliva samples within one run. The within-assay variability is shown below:

Mean (pg/ml)	129.58	865.37	2203.2
SD	7.72	55.75	212.28
CV (%)	5.96	6.44	9.63
n =	20	20	20

9.4.2 Inter-Assay

The inter-assay (between-run) variation was determined by duplicate measurements of 3 saliva samples over 10 days.

Mean (pg/ml)	111.33	1191.83	2285.5
SD	9.53	111.50	231.36
CV (%)	8.6	9.4	10.1
n =	10	10	10

9.5 Recovery

Using the Calibrator Matrix four spiking solutions were prepared (A = 2000 pg/ml, B = 5000 pg/ml, C = 10000 pg/ml, D = 20000 pg/ml). A 25 µl aliquot of each solution was spiked into 475 µl of different salivas, for a spiking ratio of 1 to 20, leaving the saliva matrix of the spiked samples relatively intact. All samples were then measured by Salivary Progesterone procedure. To calculate expected values 95% of the unspiked values were added to 5% of the spiking solution concentrations.

Sample	Measured (pg/ml)	Expected (pg/ml)	Recovery (%)	Sample	Measured (pg/ml)	Expected (pg/ml)	Recovery (%)
1	204	-	-	4	271	-	-
	302	294	102,7%		439	507	86,6%
	414	444	93,2%		593	757	78,3%
	558	694	80,4%		1089	1257	86,6%
2	50	-	-	5	954	-	-
	151	148	102,0%		1218	1156	105,4%
	249	298	83,6%		1299	1406	92,4%
	470	548	85,8%		2365	1906	124,1%
3	1852	-	-	6	122	-	-
	2061	2009	102,6%		386	366	105,5%
	2255	2259	99,8%		582	616	94,5%
	2734	2759	99,1%		1098	1116	98,4%

9.6 Linearity

Six saliva samples were assayed undiluted and diluted with the calibrator matrix.

Saliva	Dilution	Observed (O)	Expected (E)	O/E %
1	native	1242	-	-
	1 in 2	642	621	103,4%
	1 in 4	354	310	114,2%
	1 in 8	182	155	117,4%
2	native	862	-	-
	1 in 2	440	431	102,1%
	1 in 4	223	216	103,2%
	1 in 8	101	108	93,5%
3	native	1704	-	-
	1 in 2	958	852	112,4%
	1 in 4	470	426	110,3%
	1 in 8	247	213	115,9%
4	native	2527	-	-
	1 in 2	1290	1264	102,1%
	1 in 4	580	632	91,8%
	1 in 8	269	316	85,1%
5	native	324	-	-
	1 in 2	151	162	93,2%
	1 in 4	86	81	106,2%
	1 in 8	40	41	97,5%
6	native	1724	-	-
	1 in 2	856	862	99,3%
	1 in 4	354	431	82,1%
	1 in 8	190	216	88,0%

10 LIMITATIONS OF PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

10.1 Drug Interferences

Until now no substances (drugs) are known influencing the measurement of Progesterone in a saliva sample.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DEMEDITEC.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12 REFERENCES

1. Authors: Charles D. West, Damodar K. Mahajan, Virginia J. Chavré et.al.
Title: Simultaneous Measurement of Multiple Plasma Steroids by Radioimmunoassay Demonstrating Episodic Secretion.
Published in: Journal of Clinical Endocrinology and Metabolism **1973** (36 No.6) pages 1230 – 1236.
2. Authors: Ross F Vining, Robynne A McGingley, Richard G. Symons
Title: Hormones in Saliva: Mode of Entry and Consequent Implications for Clinical Interpretation
Published in: Clinical Chemistry **1983**, Vol 20 (10), pages 1752 - 1756
3. Authors: Jung K. Choe, Firyal S. Khan-Dawood, M. Yusoff Dawood
Title: Progesterone and Estradiol in the saliva and plasma during the menstrual cycle.
Published in: Amer. J. of Obstetrics & Gynecology **1983**, Vol 147 (5), pages 557 - 562
4. Authors: L. Cedard, A. Guichard, Y. Janssens, G. Tanguy, P. Boyer, JR Zorn
Title: Progesterone and Estradiol in saliva after in vitro fertilization and embryo transfer.
Published in: Fertility and Sterility **1987**, Vol 47 (2), pages 278 - 283
5. Authors: Ross F. Vining, Robynne A. McGinley
Title: The Measurement of Hormones in Saliva: Possibilities and Pitfalls
Published in: J. Steroid Biochem. **1987** (27 no.1-3) pages 81-94

6. Authors: Marguerite J.Mc.Neely, Michael R.Soules
Title: The diagnosis of luteal phase deficiency: a critical review
Published in: Fertility and Sterility, **1988** (50 No.1) pages 1 – 15
7. Authors: Michael R. Soules, Donald K. Clifton, Robert A. Steiner et.al.
Title: The Corpus Luteum: Determinants of Progesterone Secretion in the Normal Menstrual Cycle.
Published in: Obstetrics and Gynaecology **1988** (71 no.5) pages 659 – 666
8. Authors: Taru Vuorento, Ilpo Huhtaniemi
Title: Daily levels of salivary Progesterone during menstrual cycle in adolescent girls.
Published in: Fertility and Sterility **1992** (58:4) pages 685 - 690.
9. Authors: M.M. Finn, J.P. Gosling, D.F. Tallon, S. Baynes, F.P. Meehan et.al.
Title: The frequency of salivary Progesterone sampling and the diagnosis of luteal phase insufficiency.
Published in: Gynecol. Endocrinol. **1992** 6 pages 127 - 134
10. Authors: T. Blom, A. Ojanotko-Harri, M. Laine, I. Huhtaniemi
Title: Metabolism of Progesterone and Testosterone in human parotid and submandibular salivary glands *in vitro*.
Published in: J Steroid Biochem Molec Biol **1993**, Vol. 44 No.1, pages 69-76
11. Authors: Torsten M.Delfs, Olaf G.Naether, Susann Klein, Freimut A..Leidenberger, Patrick Fottrell, Ralf C.Zimmermann
Title: 24-Hour profiles of salivary Progesterone
Published in: Fertility and Sterility **1994** (62 no.5) pages 960 – 966
12. Authors: Peter O'Leary, Peter Feddema, Katharine Chan, Mario Taranto, Sharon Evans
Title: Salivary, but not serum or urinary levels of Progesterone are elevated after topical application of Progesterone cream to pre- and postmenopausal women.
Published in: Clinical Endocrinology **2000**, Vol 53, pages 615 - 620.
13. Authors: M. Gröschl, M. Rauh, P. Schmid, H-G Dörr
Title: Relationship between salivary Progesterone, 17OH-Progesterone, and Cortisol levels throughout the normal menstrual cycle of healthy postmenarche girls.
Published in: Fertility & Sterility **2001**, Vol 76 (3), pages 615 - 617
14. Authors: M. Gröschl, M. Rauh, P. Schmid, HG Dörr
Title: Relationship between salivary Progesterone, 17OH-Progetserone, and Cortisol levels throughout the normal menstrual cycle of healthy postmenarcheal girls.
Published in: Fertility and Sterility **2001**, Vol 76 (3), pages 615 - 617
15. Authors: T. Vimpeli, H. Tinkanen, H. Huhtala, L. Rönnerberg, E. Kujansuu
Title: Salivary and serum Progesterone concentrations during two luteal support regimens used in in-vitro fertilization treatments
Published in: Fertility & Sterility **2001**, Vol.4, pages 847-848
16. Authors: Nike Stikkelbroeck, C Sweep, Didi Braat, Ad Hermus, Barto Otten
Title: Monitoring of menstrual cycles, ovulation, and adrenal suppression by saliva sampling in female patients with 21-hydroxylase deficiency.
Published in: Fertility and Sterility **2003**, Vol 80 (4), pages 1030 - 1036
17. Authors: M Gröschl, M Rauh, HG Dörr
Title: Circadian rhythm of salivary Cortisol, 17OHP, and Progesterone in healthy children.
Published in: Clin Chem. **2003**, Vol. 49 (10), pages 1688-1691
18. Authors: RT Chatterton, ET Mateo, N Hou, AW Rademaker, S Acharya, VC Jordan
Title: Characteristics of salivary profiles of Estradiol and Progesterone in postmenopausal women.
Published in: Journal of Endocrinology **2005**, Vol. 186, pages 77-84
19. Author: John G. Lewis
Title: Steroid Analysis in Saliva: An overview (Review article)
Published in: Clin Biochem Rev., August **2006**, Vol.27, pages 139-146
20. Authors: L.J. Bond, E.T. Vella, Y. Kiparissis, K.E. Wynne-Edwards
Title: Anthropometry and body composition do not predict bioavailable Androgen or Progesterone concentration in adolescent girls.
Published in: Am J Hum Biol **2006**, Vol.18, pages 639-653
21. Authors: V.J. Vitzthum, H. Spielvogel, J. Thornburg, B. West
Title: A prospective study of early pregnancy loss in humans
Published in: Fertil Steril **2006**, Vol.86, pages 373-379
22. Authors: MM Wirth, EA Meier, BL Fredrickson, OC Schultheiss
Title: Relationship between salivary Cortisol and Progesterone levels in humans.
Published in: Biological Psychology **2007**, Vol.74, pages 104-107

23. Authors: A. Nunez-de-la-Mora, R.T. Chatterton, O.A. Choudhury, G.R. Bently
Title: Childhood Conditions Influence Adult Progesterone Levels
Published in: www.PLOSmedicine.org **2007**, Vol.4 (5) pages 1-8
24. Authors: B.K. Gandara, L. Leresche, L. Mancl
Title: Patterns of salivary Estradiol and Progesterone across the menstrual cycle.
Published in: Ann. NY Acad. Sci. **2007**, Vol. 1098, pages 446-450

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