
Product information

Information about other products is available at: www.demeditec.com

User's Manual

TPA ELISA

Enzyme Immunoassay for the quantitative determination of TPA (Tissue Polypeptide Antigen) in human serum

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1 INTENDED USE

TPA ELISA is a one step *in vitro* diagnostic assay for the quantitative determination of cytokeratin 8 and 18 in serum. The assay is a sensitive indicator of tumor cell activity and is useful in the management of patients with carcinomas of epithelial origin.

2 PRINCIPLE OF THE ASSAY

TPA ELISA is a solid phase sandwich assay based on immunochemical technique. Standards, controls and samples react simultaneously with solid phase catcher antibodies (6D7 and 3F3) and the HRP conjugated detector antibody during incubation in the microstrip wells. After washing, the TMB substrate is added. Subsequently the reaction is stopped and the absorbance is read. The developed color is directly proportional to the concentration of the analyte.

3 ASSAY SPECIFICITY

TPA measures key epitopes of **TPA (Tissue Polypeptide Antigen)** fragments. The monoclonal 6D7 and 3F3 antibodies used in the test are specific for cytokeratin 8 and 18, with no detectable cross reactivity to other tumor associated antigens that may be present in patient samples.

4 SPECIMENS

Serum samples are recommended. Enough blood should be collected to be sufficient for 2 x 100 µL serum (duplicates) at each analysis. If the analysis will be performed within 24 h, the serum should be refrigerated (2 - 8 °C). If delayed analysis, serum should be frozen in aliquots (<-18 °C). Avoid repeated thawing and freezing. Do not use serum samples that are grossly lipemic, hemolysed or contaminated.

5 PRECAUTIONS FOR USERS

1. TPA ELISA is for *in vitro* diagnostic use only.
2. Wear protective gloves and protective goggles.
3. Do not use the kit after expiry date.
4. Do not mix reagents from different lots.
5. All patient specimens should be regarded as contagious and handled and disposed of according to appropriate regulations.
6. Avoid microbiological contamination of reagents.
7. Analysis should be performed according to GLP.
8. The accuracy of the test is related to adherence to the assay procedure and accurate volume pipetting.
9. The Stop Solution contains 1 N acidic solution, which might cause irritation on skin, and is harmful to eyes. In case of contact flush with plenty of water and seek medical advice.
10. The TMB Substrate might cause irritation on skin and eyes. In case of contact flush with plenty of water and seek medical advice.
11. ProClin 300 (60 ppm) used as preservative in this product might be allergenic. In case of contact flush with plenty of water and seek medical advice.
12. Material Safety Data Sheets are available on request.

6 MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader (wavelength 450 nm).
- Microplate shaker (oscillation ~450 rpm).
- Microplate wash equipment.
- Routine laboratory equipment, *e.g.* precision pipettes and vortex.

7 COMPONENTS IN THE TPA ELISA

Materials supplied for 96 determinations.

1. **TPA Coated Microstrips:** 1 plate, 96 dry wells (12 x 8), coated with monoclonal anti-cytokeratin 8/18 antibodies (6D7 and 3F3). Packed in aluminum bag with desiccating device. Ready for use.
2. **TPA ELISA HRP Conjugate:** 2 vials, 0.5 mL/vial, conjugated antibody in protein stabilized buffer, pH 7.5 (11 x conc.). Should be diluted with TPA Diluent (Standard 0 ng/mL). Blue colored. Preservative added.
3. **TPA Diluent (Standard 0 ng/mL):** 1 vial, 15 mL, sample diluent and standard 0 ng/mL, protein stabilized buffer, pH 7.5. Preservative added. Ready for use.
4. **TPA ELISA Standard (1, 2, 5, 10, 15 ng/mL):** 5 vials standard, 1 mL/vial, TPA standard material in protein stabilized buffer, pH 7.5. Concentrations as stated on the QC data sheet. Preservative added. Ready for use.
5. **TPA Control (Low, High):** 2 vials control (1 mL/vial), Lyophilized Controls, TPA standard material in protein stabilized buffer, pH 7.5. Yellow colored. Preservative added.
6. **Wash Tablet:** 1 package, 1 tablet/package. The tablet should be dissolved in 500 mL fresh deionized water.
7. **TMB Substrate:** 1 vial, 22 mL. Protect from light and keep lid tightly closed. Do not sample more than what is needed for the analysis. Ready for use.
8. **Stop Solution:** 1 vial, 12 mL, 1 N acidic solution. Ready for use.
9. **Sealing Tape:** 1 sheet, Sealing tape for microstrips.
10. **TPA ELISA Certificate:** 1 protocol. Certificate of lot content.

8 ASSAY PROCEDURE

8.1 Reconstitution of Controls

Reconstitute the lyophilized kit controls the day before you run the test.

Add 1 mL fresh deionized water to each vial and keep the vials overnight at 2 °C to 8 °C.

8.2 Procedure

The assay should be performed at room temperature (RT; 22 °C ± 6 °C).

1. Allow all reagents and samples to adjust to RT. Vortex all reagents prior to use.
2. Based on the number of strips needed, dilute TPA ELISA HRP Conjugate (11 x conc.) with TPA Diluent (5.0 mL/vial). Mix thoroughly.
3. Pipette 100 µL standards, controls or samples per well (duplicates). Start with two empty wells for background absorbance measurement (blank).
4. Add 100 µL diluted HRP Conjugate to each well, except the two empty wells. Cover the strips with the supplied Sealing Tape. *NB! Steps 4 and 5 should be performed sequentially without interruption.*
5. Incubate for 2 h ± 2 min on shaker (~450 rpm).
6. Prepare the wash solution. Dissolve Wash Tablet in 500 mL fresh deionized water.
7. Aspirate and wash the wells 3 x 300 µL with wash solution.
8. Add 200 µL TMB Substrate to each well, including the two empty wells. Incubate in darkness for 15 ± 1 min.
9. Add 100 µL Stop Solution to each well. Agitate on shaker 1 min (~450 rpm).
10. Read the absorbance at 450 nm, within 30 min after addition of the Stop Solution.
11. Calculate the cytokeratin 8 and 18 concentration (ng/mL) of the samples. Samples showing concentrations >15 ng/mL value should be suitably diluted with TPA Diluent (Standard 0 ng/mL) before repeated analysis.

9 PROCESSING OF RESULTS

Manual calculation or by using a computer software for handling ELISA-type data (curve fitting - Spline smoothed). For generation of valid data, ensure that included controls are within range.

Manual processing of results: Correct each absorbance value by subtracting the background absorbance (blank). Estimate the mean value for each duplicate. Construct a standard curve by plotting the mean absorbance value for each standard (y-axis) against the corresponding concentration (x-axis). Determine the concentrations of the samples against the standard curve.

10 REAGENT STORAGE

The kit should be stored at 2 - 8 °C. **Do not freeze!** Store reagents in their original containers if not used at once. Reseal the Microstrip bag, including the desiccating device, if not all strips are used at once. The wash solution is stable for 4 weeks when stored at 2 - 8 °C. The diluted TPA ELISA HRP Conjugate is stable for 4 weeks when stored at 2 - 8 °C. The reconstituted controls are stable for 4 weeks when stored at 2 - 8 °C.

11 EXPECTED VALUES

We recommend that each laboratory establish their own normal range (cut-off value), due to differences between laboratories and locales with respect to population, laboratory techniques and selection of reference groups. The following values should only be a guideline. The normal values of the DEMEDITEC TPA ELISA were determined by measuring the values of apparently healthy individuals:

0 – 1.46 ng/mL (5 – 95 % Percentile)

12 LIMITATIONS OF THE PROCEDURE

The assay values should be interpreted in conjunction with all available clinical information. The TPA ELISA is NOT to be used as a cancer screening test.

Has any of the patients tested undergone any kind of treatment for cancer?

Different kind of medications may affect the test result.

Is it known if any of the patients tested have been taking other medications?

The test results may be affected (increased) by general infections like a cold or influenza. Furthermore, increased values may also be found in case of pregnancy, liver disease and renal failure. If a temporary infection is suspected, it may be necessary to re-test in two weeks.

13 ASSAY CHARACTERISTICS

13.1 Measuring range

The measuring range is 0-15 ng/mL. The assay has no "high-dose hook effect" up to 500 ng/mL.

13.2 Reproducibility

The intra-assay reproducibility of the standard curve has a typical CV of <5 % (in duplicates). Inter-assay imprecision is characteristically of 5-10 % CV.

13.3 Sensitivity

The minimal detectable concentration in TPA ELISA is 0.1 ng/mL, defined as the concentration of TPA that corresponds to the absorbance being three standard deviations from the absorbance of standard 0 ng/mL.

13.4 Recovery

Specified quantities of TPA were added to human serum specimens. The recovery was 88-98 %.

14 WARRANTY

The performance data presented here were obtained using the procedure indicated. Any change or modification in the procedure, not recommended by DEMEDITEC, may affect the results. In such event DEMEDITEC disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and the fitness for use.

15 REFERENCES / LITERATURE

1. Stigbrand, T. *et al.* Epitope specificity of 30 monoclonal antibodies against cytokeratin antigens: The ISOBM TD5-1 Workshop. *Tumor Biol* 1998; 19:132-152.
2. Bergqvist, M. *et al.* Cytokeratin 8 and 18 fragments measured in serum and their relation to survival in patients with non-small cell lung cancer. *Anticancer Res* 1999; 19:1833-1836.
3. Sütterlin, M. *et al.* Clinical value of CYFRA 8/18 and TPS in the diagnosis and follow up of invasive breast cancer. *Anticancer Res* 1997; 17:2963-2966.

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15 BIBLIOGRAFIA

1. Stigbrand, T. et al. Epitope specificity of 30 monoclonal antibodies against cytokeratin antigens: The ISOBM TD5-1 Workshop. *Tumor Biol* 1998; 19:132-152.
2. Bergqvist, M. et al. Cytokeratin 8 and 18 fragments measured in serum and their relation to survival in patients with non-small cell lung cancer. *Anticancer Res* 1999; 19:1833-1836.
3. Sütterlin, M. et al. Clinical value of CYFRA 8/18 and TPS in the diagnosis and follow up of invasive breast cancer. *Anticancer Res* 1997; 17:2963-2966.

SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Français	Español	Italiano
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	European Conformity	CE-Konformitätskennzeichnung	Conformité aux normes européennes	Conformidad europea	Conformità europea
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Numéro de catalogue	Número de catálogo	Numero di Catalogo
	Lot. No. / Batch code	Chargen-Nr.	Numéro de lot	Número de lote	Numero di lotto
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore
Content	Content	Inhalt	Conditionnement	Contenido	Contenuto
Volume/No.	Volume / No.	Volumen/Anzahl	Volume/Quantité	Volumen/Número	Volume/Quantità