
Product information

Information about other products is available at: www.demeditec.com

User's Manual

Homocysteine ELISA

Enzyme Immunoassay for the quantitative determination of total L homocysteine in human serum or plasma

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1 INTENDED USE

The Homocysteine ELISA is intended for the quantitative determination of total L-homocysteine in human serum or plasma. The device can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia and homocystinuria.

2 SUMMARY AND EXPLANATION OF THE TEST

Homocysteine (Hcy) is a thiol-containing amino acid produced by the intracellular demethylation of methionine. Hcy is exported into plasma where it circulates mostly in its oxidized forms bound to plasma proteins.[1, 2, 3, 4] Smaller amounts of reduced homocysteine and disulfide homocystin (Hcy-SS-Hcy) are present. Total homocysteine represents the sum of all Hcy species found in plasma and serum (free plus protein-bound).

Hcy is either metabolised to cysteine or to methionine. In the vitamin B₆ dependent transsulphuration pathway Hcy is irreversibly catabolised to cysteine. A major part of Hcy is remethylated to methionine, mainly by the folate and cobalamin- dependent enzyme methionine synthase. Hcy accumulates and is excreted into the blood when these reactions are impaired. [2,4]

Severely elevated concentrations of Hcy are found in subjects with homocystinuria, a rare genetic disorder of the enzymes involved in the metabolism of Hcy. Patients with homocystinuria exhibit mental retardation, early arteriosclerosis and arterial and venous thromboembolism.[1, 5] Hcy reducing therapy improves the prognosis for this disease.5 Other less severe genetic defects which lead to moderately elevated levels of Hcy are also found.[6, 7,8]

Epidemiological studies have investigated the relationship between Hcy levels in blood and cardiovascular disease (CVD). A meta analysis of 27 epidemiological studies, including more than 4000 patients, estimated that a 5 µmol/L increase in Hcy was associated with an odds ratio for coronary artery disease (CAD) of 1.6 for men and 1.8 for women, or the same that is associated with 0.5 mmol/L (20 mg/dL) increase in cholesterol. Peripheral arterial disease also showed a strong association.[9]

Certain patient groups with anemia and/or asthenia also demonstrate increased levels of plasma- or serum Hcy. [10,11]

Patients with chronic renal disease experience an excess morbidity and mortality due to arteriosclerotic CVD. Elevated concentration of Hcy is a frequently observed finding in the blood of these patients. Although such patients may lack some of the vitamins involved in the metabolism of Hcy, the increased levels of Hcy are mainly due to impaired removal of Hcy from the blood by the kidney.[12, 13]

Drugs such as methotrexate, carbamazepine, phenytoin, nitrous oxide and penicillamine interfere with the Hcy metabolism and may give elevated levels of Hcy.[14, 15]

3 ASSAY PRINCIPLE

The Homocysteine ELISA is an enzyme immunoassay for the determination of Hcy in blood.[16] Protein-bound Hcy is reduced to free Hcy and enzymatically converted to S-adenosyl-L-homocysteine (SAH) in a separate procedure prior to the immunoassay.[17] The enzyme is specific for the L-form of homocysteine, which is the only form present in the blood.

Reduction

Hcy, mixed disulfide and protein-bound forms of Hcy in the sample are reduced to free Hcy by use of dithiothreitol (DTT).

Enzymatic conversion

Hcy in the test sample is converted to S-adenosyl-L-homocysteine by the use of SAH hydrolase and excess adenosine (Ad).

The following solid-phase enzyme immunoassay is based on competition between SAH in the sample and immobilised SAH bound to the walls of the microtitre plate for binding sites on a monoclonal anti-SAH antibody. After removal of unbound anti-SAH antibody, a secondary rabbit anti-mouse antibody labelled with the enzyme horse radish peroxidase (HRP) is added. The peroxidase activity is measured spectrophotometrically after addition of substrate, and the absorbance is inversely related to the concentration of Hcy in the sample.

4 WARNINGS AND PRECAUTIONS

1. For in-vitro Diagnostic Use Only.
2. Reagent D contains 0.15% merthiolate ($\leq 0.074\%$ mercury), and is classified as "Harmful". Please handle and dispose of properly (See section "Product Safety Information").
3. 0.01% merthiolate is used as preservative in some reagents. Each kit contains less than 0.028% mercury. Please handle and dispose of appropriately.
4. Reagent F contains mouse antibodies and Reagent G contains rabbit antibodies.
5. Reagent S contains 1.6 N acidic solution, and is classified as "Irritant". Please handle and dispose of properly (See section "Product Safety Information").
6. Calibrators, Controls, Reagent A and Reagent E contain less than 0.10% sodium azide as preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. Controls contain sera originating from human blood samples. The source materials have been tested and found to be negative for Hepatitis B Surface Antigen (HBsAG), HIV-1 Antigen (HIVAg), HCV antibody, HIV-1/2 antibody, HTLV-1/2 antibody and Hepatitis B core Antibody (HBc). However, blood derivatives should be handled according to recommended procedures for handling infectious material. HHS publication no. (CDC) 93-8395 [18] or local/national guidelines on laboratory safety procedures should be consulted.
8. Reagents with different lot numbers must not be interchanged.
9. Do not use the kit after the expiration date on the outer box.

5 KIT COMPONENTS

Homocysteine ELISA Kit, 96 wells

Kit Components	Solution	Component Description	Volume
[REAG A] / BUF	Assay Buffer	Phosphate buffer, sodium azide	54 mL
[REAG B] / ADENO DTT	Adenosine/DTT	Adenosine / dithiothreitol, citric acid	3.5 mL
[REAG C] / SAH HYDROL	SAH-hydrolase	Recombinant S-adenosyl-L-homocysteine hydro- lase, trisbuffer, glycerol, methylparaben	3.5 mL
[REAG D] / ENZ INH	Enzyme Inhibi- tor	Merthiolate, phosphate buffer	55 mL
[REAG E] / ADENO DEAM	Adenosine de- aminase	Adenosine deaminase, phosphate buffer, sodium azide, BSA, phenol-red dye	55 mL
[REAG F] / SAH Ab	a-SAH Antibody	Monoclonal mouse-anti-S-adenosyl-L- homocysteine antibody, BSA, merthiolate	25 mL
[REAG G] / ENZ CONJ	Enzyme Conju- gate	Rabbit anti-mouse-antibody enzyme conjugate, BSA, horse radish peroxidase, blue dye	15 mL
[REAG H] / SUB TMB	Substrate Solu- tion	N-methyl-2-pyrrolidon, propyleneglycol	15 mL
[REAG S] / STOP SOLN	Stop Solution	1.6 N acidic solution	20 mL
BUF WASH 10x	Wash Buffer	Phosphate buffer, merthiolate, Tween 20, BSA	60 mL
CAL 1 - 6	Standards	S-adenosyl-L-homocysteine (2, 4, 8, 15, 30, 50 $\mu\text{mol/L}$) in buffer with preservative	6 x 1.5 mL
SORB MT	Microtiter wells	Coated with S-adenosyl-L-homocysteine	12 x 8 wells

BUF WASH 10x is concentrated and must be diluted (1+9) with purified water before use.
All other components are ready to use.

Homocysteine ELISA Control kit (DE3329)

Kit Components	Component Description	Volume
[CONTROL L]	7.0 µmol/L homocysteine in diluted serum samples of human origin, phosphate buffer and preservative	1.5 mL
[CONTROL M]	12.5 µmol/L homocysteine in diluted serum samples of human origin, phosphate buffer and preservative	1.5 mL
[CONTROL H]	25.0 µmol/L homocysteine in diluted serum samples of human origin, phosphate buffer and preservative	1.5 mL

All controls are ready to use.

Homocysteine ELISA Wash Buffer

Kit Components	Component Description	Volume
BUF WASH 10x	Phosphate buffer, merthiolate, Tween 20, BSA	1000 mL

BUF WASH 10x is concentrated and must be diluted (1+9) with purified water before use.

5.1 Materials required but not provided in the kit:

- Homocysteine controls (see section “Quality Controls” for more information)
- Plastic or glass tubes for pre-treatment of samples
- Pipettes / multipipettes 25 µL, 100 µL, 200 µL and 500 µL or 8 channel multipipette for 100 µL and 200 µL
- Volumetric flask 50 mL and 600 mL
- Incubator, 37 °C
- Washer and reader (450 nm) for microtiter plates

6 PREPARATION AND STORAGE OF KIT COMPONENTS

- Components should be stored refrigerated (2 - 8 °C). Store all bottles upright and tightly capped. The components are stable until the stated expiration date when stored and handled as directed. Once the components in the Homocysteine ELISA Kit are opened they are stable for 12 weeks when stored at 2-8 °C.
- The sample pre-treatment solution has to be made by mixing Reagent A, B and C (see Section “Procedure”). The solution is stable for one hour and has to be freshly made for each run.
- The Wash Buffer must be diluted (1+9) with distilled water before use. The prepared Wash buffer is stable for 4 weeks when stored at room temperature (18-25 °C).
- Reagent D and H are stored in dark bottles to avoid exposure to light.
- It is important that the microtiter strips are kept dry, i.e. in the sealed bag with drying capsules, and stored refrigerated. Equilibration for a minimum of two hours is required to reach room temperature (18-25 °C). Leave the strips in the bag during equilibration.
- Only the necessary number of microtiter strips should be kept in the frame during the run. Unused strips should be kept in the sealed bag with drying capsules.
- Avoid exposure of the kit to temperatures exceeding 37 °C as this may denature the enzymes.

7 SPECIMEN COLLECTION AND PREPARATION

EDTA-plasma or serum may be used with the Homocysteine ELISA assay.

As synthesis of Hcy continues in red blood cells after drawing, it is very important to prepare specimens as follows:

- Serum samples should be allowed to clot for no more than 30 minutes before centrifugation and separation of serum. Serum samples should be kept on ice prior to separation.
- EDTA-plasma samples must be centrifuged or put on ice immediately after drawing. EDTA-plasma samples may be kept on ice for up to 6 hours prior to separation by centrifugation.

Food consumption can affect circulating homocysteine levels. Protein rich meals give higher homocysteine values and should be avoided late in the day before sampling.[19, 20]

Standardised sampling procedures are crucial due to the above mentioned influencing factors. Complete mixing of thawed samples is required before use.

Plasma or serum samples may be stored for 12 weeks at 2 - 8 °C, for up to 3 weeks at room temperature (18 - 25 °C) and have been shown to be stable for at least 8 months if frozen at minus 20 °C.

8 LIMITATIONS

- If an automatic pipetting station is used, thorough washing of the tubing after addition of the blue coloured Reagent G may be needed - preferably with diluted acid followed by water. Any remaining solution in the tubing will interfere with the next assay step; i.e. addition of Reagent H.
- The washing procedure is critical for obtaining good precision. If manual washing is required, use 4 times 350 µL instead of 3 times 400 µL. After washing, empty the wells on paper towels.
- Avoid exposure of the kit to temperatures exceeding 37 °C as this may denature the enzymes.
- Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibody (HAMA). HAMA, present in serum or plasma specimens, may interfere in immunoassays which utilise mouse monoclonal antibodies. These specimens should not be assayed with the Homocysteine ELISA assay.
- Specimens from patients taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anti-convulsants or 6-azauridine triacetate, may have elevated levels of homocysteine due to metabolic interference with the homocysteine metabolism.

9 PROCEDURE

Make sure all solutions and microtiter strips are equilibrated to room temperature before use. Leaving the kit at room temperature over night is recommended. We recommend running the standards in duplicate and to performing a new calibration curve for each run to avoid run-to-run variations using coated microtiter plates.

9.1 Sample pre-treatment procedure

1. Sample pre-treatment solution must be made up no more than 1 hour prior to the start of the assay. Volume needed per 10 samples (no dead volume calculated):
4.5 mL REAG A
0.25 mL REAG B
0.25 mL REAG C
Mix.
 2. Dilute calibrators and samples/controls in plastic or glass tubes as follows:
25 µL standard/sample/control
+ 500 µL sample pre-treatment solution
Mix well.
Incubate for 30 minutes at 37 °C (Cap the tubes or cover with parafilm during incubation).
- Note:** Proceed with step 3 before the samples have cooled.
3. Add 500 µL REAG D
Mix well.
 4. Incubate for 15 minutes at 18-25 °C.
 5. Add 500 µL REAG E
Mix well.
 6. Incubate for 5 minutes at 18-25 °C.

9.2 Microtitre plate procedure

7. Pipette 25 µL diluted calibrator / sample / control from step 4 into the wells of the SAH-coated microtitre strips.
8. Add 200 µL REAG F to each well.
9. Incubate for 30 min at 18-25 °C.
Use the enclosed lid during all incubations.
10. Wash with diluted Wash buffer.
Use 3 x 400 µL. If manual washing is required, use 4 times 350 µL instead of 3 times 400 µL. After washing, empty the wells on paper towels.
11. Add 100 µL REAG G to each well.
12. Incubate for 20 min at 18-25 °C.
13. Wash with diluted Wash buffer.
Use 3 x 400 µL. If manual washing is required, use 4 times 350 µL instead of 3 times 400 µL. After washing, empty the wells on paper towels.
14. Add 100 µL REAG H to each well.
15. Incubate for 10 min at 18-25 °C.
16. Add 100 µL REAG S to each well.
17. Shake and read at 450 nm within 15 minutes (Automatic plate shaker is preferred to ensure proper mixing).

10 INTERPRETATION OF RESULTS

Results should be interpreted considering all other test results and the clinical status of the patient. We recommend that a four parameter logistic curve fit is used for preparing the calibration curve and calculation of unknown samples.

11 QUALITY CONTROL

We recommend that each laboratory use a homocysteine control with known value. DEMEDITEC provides a set of low, medium and high controls, (REF DE3329). The controls contain L-homocysteine in processed human serum at the following concentrations:

Control	Mean value Hcy (µmol/L)	Range Hcy (µmol/L)
Low	7.0	5.6 - 8.4
Medium	12.5	10.0 – 15.0
High	25.0	20.0 – 30.0

12 REFERENCE RANGE

The reference range should be determined by each laboratory to confirm the characteristics of the population being tested. As a point of reference, the following data may be used until the laboratory has analysed a sufficient number of samples to determine its own reference range.

The Hcy concentration in plasma or serum of healthy individuals varies with age, gender, geographical area and genetic factors. Scientific literature reports reference values for adult male and females between 5 and 15 µmol/L, men having higher values than women, and post menopausal women having higher homocysteine values than pre menopausal women.[21, 22, 23] Hcy values will normally increase with age, giving a reference range among an elderly population (> 60 years) of 5 - 20 µmol/L.[24] In countries with folic acid fortification programmes, reduced levels of Hcy may be observed.[25, 26]

Samples from 382 males and females (100 Scandinavians; 54 males aged 30 - 60, 46 females aged 29 - 70. 185 Hispanics; predominantly males aged 20 - 65. 97 Americans; 54 males aged 16 - 74, 43 females aged 15 - 79), apparently healthy, without information on current medications, disease states or known risk conditions for elevated homocysteine, were tested using the Homocysteine ELISA.

The median value of the homocysteine concentration among Scandinavians was 8.4 µmol/L, among Hispanics 8.9 µmol/L and among Americans 9.3 µmol/L.

The homocysteine reference range was established based on 95% confidence limits as 5 - 15 µmol/L for the Scandinavian population, 3.6 - 15.0 µmol/L for the American population and 2.9 - 16.0 µmol/L for the Hispanic population.

13 MEASURING RANGE

The calibrator range is from 2 to 50 µmol/L.

14 PERFORMANCE DATA**14.1 Assay Precision**

Precision of the Homocysteine ELISA was evaluated according to NCCLS Protocol EP5-T2. Three levels of controls were assayed for 20 days with 4 replicates per run at each level. Precision data are summarised in

Samples	Average Hcy µmol/L	Within run precision	Total precision
Control Low	6.1	8 %	10 %
Control Medium	10.5	7 %	9 %
Control High	20.6	8 %	10 %

14.2 Limit of quantification

The quantification limit (CV < 20%) is 1.0 µmol/L.

14.3 Linearity of diluted plasma samples

If the homocysteine concentration of a sample exceeds the range of the calibration curve, the sample should be diluted with Reagent A and reanalysed.

The linearity was evaluated by diluting four high patient samples with varying amounts of Reagent A as diluent.

Linear regression analysis gave: Slope: 0.98
Intercept: - 0.4 µmol/L
Correlation coefficient r^2 : 0.99

14.4 Method Comparison

The Homocysteine ELISA was compared to the University of Bergen HPLC method.[27]

A comparison of 164 patient samples ranging from 3 - 37 µmol/L homocysteine gave the linear regression:

Slope: 0.94
Intercept: - 0.09 µmol/L
Correlation coefficient r^2 : 0.94

14.5 Interfering Substances

Bilirubin, haemoglobin, lipids, red blood cells, protein and sodium fluoride were spiked into plasma samples and tested for interference by the Homocysteine ELISA.

The assay demonstrated less than 10% interference in the presence of: bilirubin (0.5 g/L), haemoglobin (10 g/L), triglycerides (10 g/L), red blood cells (5.0% v/v), protein (80 g/L) and sodium-fluoride (10 g/L).

14.6 Cross-reactivity

Cross-reactivity was tested for compounds that may interfere with the Homocysteine ELISA.

The assay demonstrated 16% cross-reactivity in the presence of S-adenosyl-L-methionine (0.5 mmol/L) and less than 1% cross-reactivity in the presence of: Adenosine (5.0 mmol/L), cystathionine (0.5 mmol/L), L-cysteine (100 mmol/L), glutathione (100 mmol/L) and thiolactone (0.5 mmol/L).

15 PRODUCT SAFETY INFORMATION

Reagent D R-20/21/22 Harmful by inhalation, in contact with skin and if swallowed.
R-33 Danger of cumulative effects.
S-24 Avoid contact with skin.
S-28 After contact with skin, wash immediately with plenty of water.
S-36/37/39 Wear suitable protective clothing, gloves and eye/face protection.

Reagent S R-36/38 Irritating to eyes and skin.
S-25 Avoid contact with eyes.
S-26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S-37/39 Wear suitable gloves and eye/face protection.

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16 REFERENCES / LITERATURE

1. Malinow, MR. Plasma homocyst(e)ine and arterial occlusive diseases: A mini-review. *Clin Chem* 1995;41:173-176.
2. Ueland PM. Homocysteine species as components of plasma redox thiol status. *Clin Chem* 1995;41:340-342.
3. Perry IJ, Refsum H, Morris RW, Ebrahim SB, Ueland PM, Shaper AG. Prospective study of serum total homocysteine concentration and risk of stroke in middle-aged British men. *The Lancet* 1995;346:1395-1398.
4. Finkelstein JD. Methionine metabolism in mammals. *J Nutr Biochem* 1990;1:228-237.
5. Mudd SH, Levy HL, Skovby F. Disorder of transsulfuration. In: C.R. Scriver, A.L. Beaudet, W.S. Sly, et al. eds. *The metabolic basis of inherited disease* New York: McGraw-Hill; 1995:1297-1327.
6. Clarke R, Daly L, Robinson K, Naughten E, Cahalane S, Fowler B, Graham I. Hyperhomocysteinemia: An independent risk factor for vascular disease. *N Engl J Med* 1991;324:1149-1155.
7. Deloughery TG, Evans A, Sadeghi A, McWilliams J, Henner WD, Taylor LM Jr, Press RD. Common mutation in methylenetetrahydrofolate reductase. Correlation with homocysteine metabolism and late-onset vascular disease. *Circulation* 1996;94:3074-3078.
8. Schmitz C, Lindpaintner K, Verhoef P, Gaziano JM, Buring J. Genetic polymorphism of methylenetetrahydrofolate reductase and myocardial infarction. *Circulation* 1996;94:1812-1814.
9. Boushey CJ, Beresford SAA, Omenn GS, Motulsky AG. A Quantitative Assessment of Plasma Homocysteine as a Risk Factor for Vascular Disease, *JAMA* 1995;274:1049-1 057.
10. Allen, RH, Stabler SP, Savage DG, Lindenbaum J. Diagnosis of cobalamin deficiency I: Usefulness of serum methylmalonic acid and total homocysteine concentrations, *Am J Hematol* 1 990;34:90-98
11. Stabler, SP, Lindenbaum J, Allen RH. The Use of homocysteine and other metabolites in the specific diagnosis of vitamin B-12 deficiency, *J Nutr* 1996;126: 1266S-1 272S.
12. Guttormsen AB, Svarstad E, Ueland PM, Refsum H. Elimination of homocysteine from plasma in subjects with end-stage renal failure. *Irish J Med Sci* 1 995;1 64 (Suppl. 1 5):8-9
13. Bostom AG, Lathrop L. Hyperhomocysteinemia in end-stage renal disease: Prevalence, etiology, and potential relationship to arteriosclerotic outcomes. *Kidney Int* 1997;52:10-20.
14. Desouza C, Keebler M, McNamara DB, Fonseca V. Drugs affecting homocysteine metabolism. *Review. Drugs* 2002;62 (4):605-616.
15. Refsum H, Ueland PM. Clinical significance of pharmacological modulation of homocysteine metabolism. *Review. Trends Pharmacol Sci* 1990; 11:411-416
16. Frantzen F, Faaren AL, Alfheim I, Nordhei AK. An enzyme conversion immunoassay for determining total homocysteine in plasma or serum. *Clin Chem* 1998;44:31 1-316.
17. Sundrehagen E. Axis Biochemicals ASA. Enzymatic assay for homocysteine and a kit therefor. EP 623174/US5631 127.
18. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. HHS Publication No. (CDC) 93-8395. Washington, DC: US Government Printing Office, May 1999 (4th edition).
19. Ubbink JB, Vermaak Hayward WJ, van der Merwe A, Becker PJ. The effect of blood sample aging and food consumption on plasma total homocysteine levels. *Clinica Chimica Acta* 1992;207:1 19-1 28.

20. Guttormsen AB, Schneede J, Fiskerstrand T, Ueland PM, Refsum H. Plasma concentrations of homocysteine and other aminothiols compounds are related to food intake in healthy human subjects. *J Nutr* 1994;124:1934-1941.
21. Ueland PM, Refsum H, Stabler SP, Malinow R, Andersson A, Allen RH. Total homocysteine in plasma or serum: Methods and clinical applications. *Clin Chem* 1993;39:1764-1779.
22. Nehler MR, Taylor LM Jr, Porter JM. Homocysteine as a risk factor for atherosclerosis: A review. *Cardiovasc Pathol* 1997;6:1-9.
23. Lussier-Cacan S, Xhignesse M, Piolot A, Selhub J, Davignon J, Genest J Jr. Plasma total homocysteine in healthy subjects: sex-specific relation with biological traits. *Am J Clin Nutr* 1996; 64:587-593.
24. Clarke R, Woodhouse P, Ulvik A, Frost C, Sherliker P, Refsum H, Ueland PM, Khaw K- T. Variability and determinants of total homocysteine concentrations in plasma in an elderly population. *Clin Chem* 1998;44:102-107.
25. Jacques PF, Selhub J, Bostom AG, Wilson PWF, Rosenberg IH. The effect of folic acid fortification on plasma folate and total homocysteine concentrations. *N Engl J Med* 1 999;340: 1449-1454.
26. Lawrence JM, Pettitt DB, Watkins M, Umekubo MA. Trends in serum folate after food fortification. *Lancet* 1 999;354:91 5-916.
27. Fiskerstrand T, Refsum H, Kvalheim G, Ueland PM, Homocysteine and other thiols in plasma and urine: Automated determination and sample stability. *Clin Chem* 1 993;39:263-271.

SYMBOLS USED WITH DEMEDITEC ASSAYS

Symbol	English	Deutsch	Français	Español	Italiano
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	European Conformity	CE-Konformitätskennzeichnung	Conformité aux normes européennes	Conformidad europea	Conformità europea
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Numéro de catalogue	Número de catálogo	Numero di Catalogo
	Lot. No. / Batch code	Chargen-Nr.	Numéro de lot	Número de lote	Numero di lotto
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
Distributed by	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore
Content	Content	Inhalt	Conditionnement	Contenido	Contenuto
Volume/No.	Volume / No.	Volumen/Anzahl	Volume/Quantité	Volumen/Número	Volume/Quantità