
# Product information

Information about other products is available at: [www.demeditec.com](http://www.demeditec.com)

## User's Manual

# Aflatoxin B<sub>1</sub> Immuno- affinity Column

Sample Preparation and Concentration for the Quantitative Analysis  
of Aflatoxin B<sub>1</sub> with ELISA or HPLC

**REF**

DEAB1A01

96

Capacity	>200 ng
Recovery (spiked samples)	97 %

## 1. GENERAL INFORMATION

Aflatoxins belong to the class of mycotoxins. Chemically they are defined as difuranocyclopentanocumarines or difuranopentanolidocumarines, i.e. aflatoxins contain a dihydrofuran or a tetrahydrofuran ring, to which a substituted cumarin system is condensed. Out of about 20 known aflatoxins, the moulds *Aspergillus flavus* and *A. parasiticus* produce exclusively aflatoxin B<sub>1</sub>, B<sub>2</sub>, G<sub>1</sub> and G<sub>2</sub>, and all the other aflatoxins are derivatives of these four. The derivatives are developed either by metabolism in humans, animals and microorganisms or by environmental reactions.

Aflatoxins belong to the strongest mycotoxins, which act primarily in a hepatotoxic and carcinogenic way. The four main aflatoxins show a different toxicity. B<sub>1</sub> is without doubt the most toxic aflatoxin, followed by G<sub>1</sub>, B<sub>2</sub> and G<sub>2</sub>. Aflatoxin B<sub>1</sub>, however, does not show a direct toxic action. In the process of biotransformation in the liver, the lipophilic toxin is epoxidated and transformed into an active derivative, the so-called aflatoxin B<sub>1</sub>-2,3-epoxid. This highly reactive epoxid is able to react with nucleophilic regions of macromolecules. Amongst others this metabolite of aflatoxin B<sub>1</sub> binds covalently to the N-7 atom of the guanine bases of DNA. This covalent bond causes an inhibition of the DNA replication, the RNA synthesis and mutations.

Both chronic and acute intoxications are effected by aflatoxins. There are only few documented reports about acute intoxications, which are caused by uptake of mycotoxins. Of special importance for human beings are the chronic intoxications by aflatoxins. To the diseases, which develop after such chronic intoxications, belong primary liver carcinoma, hepatitis, Reye's syndrome and Kwashiorkor. Besides the generation of primary liver carcinoma, aflatoxins are presumably also responsible for other sorts of tumors, like intestinal cancer.

Contaminations with aflatoxins occur mostly with nuts and grain. In most cases aflatoxins penetrate the human body via the food. Aflatoxins are stable to heat and are only partly destroyed by boiling. In order to protect people against aflatoxin-induced diseases, there is a need for the quantitative and qualitative control of endangered foodstuff, besides appropriate hygienic precautions, which avoid the formation of aflatoxins.

The **Aflatoxin B<sub>1</sub> Immunoaffinity Column (IAC)** is a helpful tool in specific isolation and enrichment of aflatoxin B<sub>1</sub> from different matrices. After the IAC procedure the eluate can be quantitatively analyzed by various analytical techniques like the **Aflatoxin B<sub>1</sub> ELISA** or HPLC.

## 2. PRINCIPLE OF THE METHOD

The **Aflatoxin B<sub>1</sub> Immunoaffinity Column** is based on the principle of an antibody-antigen reaction. A polyclonal antibody raised against aflatoxin B<sub>1</sub> is bound on a sepharose gel, which is fixed between two frits in a polypropylene column. By passing the sample solution through the column, aflatoxin B<sub>1</sub> is bound selectively by the antibody fixed on the sepharose gel. Unbound matrix material is removed in a wash step. Bound aflatoxin B<sub>1</sub> is eluted by a small amount of methanol (100%). This eluate is used for further quantification by ELISA, HPLC or any other appropriate analytical method.

## 3. PRECAUTIONS

Full compliance of the following good laboratory practices (GLP) will determine the reliability of the results:

1. Prior to beginning of the immunoaffinity column cleanup procedure, bring all reagents to room temperature (20-25°C).
2. Once the cleanup procedure has been started, all subsequent steps should be completed without interruption.
3. Use a separate disposable tip for each specimen to prevent cross-contamination.
4. Do not use reagents after expiration date.
5. Check both precision and accuracy of the laboratory equipment used during the procedure (micropipets, etc.).

#### 4. HEALTH AND SAFETY INSTRUCTIONS

1. Do not smoke or eat or drink or pipet by mouth in the laboratory.
2. Wear disposable gloves whenever handling patient specimens or any other material of biological origin. Note that aflatoxins are very toxic substances, which may cause cancer and other serious diseases due to inhalation, ingestion or skin contact.
3. Handling and disposal of chemical products must be done according to good laboratory practices (GLP).

#### 5. CONTENT

10 immunoaffinity columns for the cleanup and concentration of aflatoxin B<sub>1</sub> containing samples.

#### 6. ADDITIONAL INSTRUMENTATION AND REAGENTS (NOT PROVIDED)

##### Instrumentation

- Syringe
- Pump unit
- Balance
- Homogenizer (mixer, magnetic stirrer)
- Clamp
- Centrifuge, filtering paper
- Funnel, beakers
- Micropipets

##### Reagents

- Double distilled/deionised water
- Methanol (100%)
- 0.01 M PBS
- Hydrochloric acid (1 N)
- Hydrochloric acid, concentrated (36%)
- Sodium hydroxide (1 N)

#### 7. SAMPLE PREPARATION

Basically all aqueous solutions can be used as samples. Make sure that the pH is at 7.0±0.5. If this is not the case, use hydrochloric acid (1 N) or sodium hydroxide (1 N) for pH adjustment. If the sample is not clear, it must be filtered or centrifuged for not risking a blockade of the column.

##### Solid samples

- An appropriate amount of solid sample (e.g. grain or nuts) is crushed in a mortar or mixer to produce a fine to medium-fine powder. If flour is used, the first step can be omitted.
- 2 g of this powder are extracted with 10 mL methanol / double distilled water (70/30 v/v). This mixture is agitated for 30 minutes on a horizontal shaker (120 / minute). Alternatively also a magnetic stirrer can be used.
- The mixture is filtrated or centrifuged afterwards for 5 minutes at 3000 g.
- The filtrate is diluted 1:10 in 0.01 M PBS (pH 7).

##### Liquid samples

- Samples with a high fat content (e.g. full cream milk) should be centrifuged prior to use (2-8°C, 3000g, 15 minutes) to separate the fat.
- After removal of the fat layer the sample can be applied directly onto the column. If the sample is highly viscous it can be diluted 1:2 to 1:10 in 0.01M PBS.
- Protein rich samples can be deproteinated prior to use. 20 µL concentrated hydrochloric acid (36%) is added to 5 mL sample and incubated for 45 minutes at 60°C. The sample is centrifuged for 10 minutes at 3000 g and the pH is adjusted with sodium hydroxide solution (1 N).

## 8. PROCEDURE

1. For each sample use a new column.
2. Let column and sample reach ambient temperature before use.
3. Take off the upper cap of the column and cut the bottom of the cap horizontally using a knife. Then place the opened cap again onto the column.
4. Connect the column to the syringe and fix the whole combination with a clamp. Place a beaker under the column.
5. Fill prepared sample into the syringe. The necessary sample volume depends on the sensitivity to be reached with the following quantification technique. Uncap the lower end of the column.
6. The sample is passed through the column using gravity. Application of slight pressure might be helpful until the first drop appears at the lower end of the column.
7. After the binding the column is washed by adding 10 mL distilled water to the syringe and passing it through the column. Remove any unbound water from the gel by passing quickly 1-2 times air through the column using the pump unit.
8. Place a clean vial under the column. Add 1 mL methanol (100%) into the dry syringe and let it pass through under gravity. In case the methanol should not pass through immediately, push the first drop through the column using the pump unit. Do not elute in less than 1 minute!
9. Push some air through the column to elute the last drop of methanol.

## 9. PERFORMANCE

### Capacity

The capacity of the **Aflatoxin B<sub>1</sub> Immunoaffinity Column** was determined to be >200 ng.






### Recovery

The recovery of spiked samples was determined to be 97 %.

## 10. REFERENCES

1. Reiß, J. (ed.): Mykotoxine in Lebensmitteln; Gustav Fischer Verlag, 1981.
2. Heathcote, J.G. und Hibbert, J.R.(eds): Aflatoxins: Chemical and biological aspects; Elsevier, Amsterdam, 1978.
3. Betina, V. (ed): Aflatoxins, sterigmatocystin and versicolorin; In: Mycotoxins, Band 9, 1. Auflage, Elsevier, Amsterdam, 1989.
4. Teuscher, E. und Lindequist, U. (eds): Biogene Gifte - Biologie, Chemie, Pharmakologie; 2. Auflage, Gustav Fischer Verlag, 1994.

## 11. WORKING SCHEME FOR ABSORPTION AND ELUTION

A	B	C	D	E
 <p>Syringe Extract Column</p>	 <p>Pump unit</p>		 <p>MeOH</p>	
Absorption of sample under gravity	Washing step 10 mL dist. Water	Removing of unbound water by passing air through IAC	Elution of Aflatoxin B <sub>1</sub> with 1 mL MeOH under gravity	ELISA or HPLC Analysis

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