Product information

Information about other products is available at: www.demeditec.com

User's Manual

Zearalenone Immuno- affinity Column

Sample Preparation and Concentration for the Quantitative Analysis
of Zearalenone with ELISA or HPLC

DEZEAA01

96

Capacity > 2000 ng
Recovery (spiked samples) 97 %

1. GENERAL INFORMATION

Zearalenone belongs to the group of mycotoxins. Similar to other mycotoxins, zearalenone is a degradation product of the secondary metabolism. It is produced by different moulds of the *Fusarium* genus. These moulds infect grain and other types of food like peanuts and beans already during their growth. When a considerable amount of zearalenone contaminated feed is taken up by cows, it can also be detected in their milk. Even in beer it could be found. Zearalenone shows a strong estrogen-like activity. Thus zearalenone can cause an enlargement of the uterus, diminution of the ovarian glands and even infertility. Zearalenone is one of the main contaminants of farm products, which can be taken up by humans and animals.

In order to protect humans of illnesses caused by zearalenone, a qualitative and quantitative control of threatened food is essential, together with appropriate hygienic measures, which help to prevent the development of zearalenone.

The **Zearalenone Immunoaffinity Column (IAC)** is a helpful tool in specific isolation and enrichment of zearalenone from different matrices. After the IAC procedure the eluate can be quantitatively analyzed by various analytical techniques like the **Zearalenone ELISA** or HPLC.

2. PRINCIPLE OF THE METHOD

The **Zearalenone Immunoaffinity Column** is based on the principle of an antibody-antigen reaction. A polyclonal antibody raised against zearalenone is bound on a sepharose gel, which is fixed between two frits in a polypropylene column. By passing the sample solution through the column, zearalenone is bound selectively by the antibody fixed on the sepharose gel. Unbound matrix material is removed in a wash step. Bound zearalenone is eluted by a small amount of methanol (100 %). This eluate is used for further quantification by ELISA, HPLC or any other appropriate analytical method.

3. PRECAUTIONS

Full compliance of the following good laboratory practices (GLP) will determine the reliability of the results:

1. Prior to beginning of the immunoaffinity column cleanup procedure, bring all reagents to room temperature (20-25 °C).
2. Once the cleanup procedure has been started, all subsequent steps should be completed without interruption.
3. Use a separate disposable tip for each specimen to prevent cross-contamination.
4. Do not use reagents after expiration date.
5. Check both precision and accuracy of the laboratory equipment used during the procedure (micropipets, etc.).

4. HEALTH AND SAFETY INSTRUCTIONS

1. Do not smoke or eat or drink or pipet by mouth in the laboratory.
2. Wear disposable gloves whenever handling patient specimens or any other material of biological origin. Note that zearalenone is a very toxic substance, which may cause cancer and other serious diseases due to inhalation, ingestion or skin contact.
3. Handling and disposal of chemical products must be done according to good laboratory practices (GLP).

5. CONTENT

10 immunoaffinity columns for the cleanup and concentration of zearalenone containing samples.

6. ADDITIONAL INSTRUMENTATION AND REAGENTS (NOT PROVIDED)

Instrumentation

- Syringe
- Pump unit
- Balance
- Homogenizer (mixer, magnetic stirrer)
- Clamp
- Centrifuge, filtering paper
- Funnel, beakers
- Micropipets

Reagents

- Double distilled/deionised water
- Methanol (100 %)
- 0.01 M PBS
- Hydrochloric acid (1 N)
- Hydrochloric acid, concentrated (36 %)
- Sodium hydroxide (1 N)

7. SAMPLE PREPARATION

Basically all aqueous solutions can be used as samples. Make sure that the pH is at 7.0 ± 0.5 . If this is not the case, use hydrochloric acid (1 N) or sodium hydroxide (1 N) for pH adjustment. If the sample is not clear, it must be filtered or centrifuged for not risking a blockade of the column.

Solid samples

- An appropriate amount of solid sample (e.g. grain or nuts) is crushed in a mortar or mixer to produce a fine to medium-fine powder. If flour is used, the first step can be omitted.
- 2 g of this powder are extracted with 10 mL methanol / double distilled water (60/40 v/v). This mixture is agitated for 30 minutes on a horizontal shaker (120 / minute). Alternatively also a magnetic stirrer can be used.
- The mixture is filtrated or centrifuged afterwards for 5 minutes at 3000 g.
- The filtrate is diluted 1:10 in 0.01 M PBS (pH 7).

Liquid samples

- Samples with a high fat content (e.g. full cream milk) should be centrifuged prior to use (2-8°C, 3000g, 15 minutes) to separate the fat.
- After removal of the fat layer the sample can be applied directly onto the column. If the sample is highly viscous it can be diluted 1:2 to 1:10 in 0.01 M PBS.
- Protein rich samples can be deproteinated prior to use. 20 µL concentrated hydrochloric acid (36%) is added to 5 mL sample and incubated for 45 minutes at 60°C. The sample is centrifuged for 10 minutes at 3000 g and the pH is adjusted with sodium hydroxide solution (1 N).

8. PROCEDURE

1. For each sample use a new column.
2. Let column and sample reach ambient temperature before use.
3. Take off the upper cap of the column and cut the bottom of the cap horizontally using a knife. Then place the opened cap again onto the column.
4. Connect the column to the syringe and fix the whole combination with a clamp. Place a beaker under the column.
5. Fill prepared sample into the syringe. The necessary sample volume depends on the sensitivity to be reached with the following quantification technique. Uncap the lower end of the column.
6. The sample is passed through the column using gravity. Application of slight pressure might be helpful until the first drop appears at the lower end of the column.
7. After the binding the column is washed by adding 10 mL distilled water to the syringe and passing it through the column. Remove any unbound water from the gel by passing quickly 1-2 times air through the column using the pump unit.
8. Place a clean vial under the column. Add 1 mL methanol (100%) into the dry syringe and let it pass through under gravity. In case the methanol should not pass through immediately, push the first drop through the column using the pump unit. Do not elute in less than 1 minute!
9. Push some air through the column to elute the last drop of methanol.

9. PERFORMANCE**Capacity**

The capacity of the **Zearalenone Immunoaffinity Column** was determined to be >2000 ng.



Recovery

The recovery of spiked samples was determined to be 97%.

10. REFERENCES

1. Reiß, J. (ed.): Mykotoxine in Lebensmitteln; Gustav Fischer Verlag, 1981.
2. Ueno, I. (ed): Trichothecenes: chemical, biological and toxicological aspects. Elsevier Science, ISBN 0-444-99661-3 (1983).
3. Scott, Peter, M. Mycotoxins. JAOAC, Vol. 70, No. 2, 276 - 281 (1987).
4. Teuscher, E. und Lindequist, U. (eds): Biogene Gifte - Biologie, Chemie, Pharmakologie; 2. Auflage, Gustav Fischer Verlag, 1994.

11. WORKING SCHEME FOR ABSORPTION AND ELUTION

A	B	C	D	E
 <p>Syringe Extract Column</p>	 <p>Pump unit</p>		 <p>MeOH</p>	
Absorption of sample under gravity	Washing step 10 mL dist. Water	Removing of unbound water by passing air through IAC	Elution of Zearalenone with 1 mL MeOH under gravity	ELISA or HPLC Analysis

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